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
Operation and Maintenance Manual

MOD. IG-65C

Mom-41
Revisión: 3
Fecha: 2023/03/21

Tecnología 2, Nave 7
08780 Pallejà
Barcelona - Spain
☎ (34) 93 334 01 12
FAX (34) 93 440 25 64
comercial@ordisi.com
www.ordisi.com

ÍNDICE

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
1.- INTRODUCTION

We would like to thank you for the trust you have placed in **ORDISI, S.A.** when you purchased a lamp from the **IGLUX** series. Your lamp model **IG-65C**, has specially designed for a recognition, for any medical specialty, but mainly in places where a total absence of calorific radiation is required, conserving very good quality of light.

Please read this manual carefully in order to ensure many years of trouble free use. If you have any requires about your **IG-65C**, please contact your nearest ORDISI distributor or contact us to our head office directly:

ORDISI, S.A.
Tecnología 2, Nave 7
08780 Pallemà (Barcelona) - SPAIN
Tel.: +34 93 334 01 12
Fax.: +34 93 440 25 64
e-mail: comercial@ordisi.com

2.- SYMBOLS AND ICONS USED IN THE PACKAGING

SYMBOL	MEANING
	SERIAL NUMBER
	PRODUCT NAME
	MANUFACTURER ADDRESS
	CE MARKING
	YEAR OF MANUFACTURE
	RECYCLE
	SEE THE USER'S GUIDE
1x2A/250 	FUSE CAPACITY
	CAUTION

	GROUNDING
	TEMPERATURE RANGE OF STORAGE AND TRANSPORT

3. - MANUFACTURER'S LIABILITIES

ORDISI, S.A. is responsible to whatever is related to the equipment safety, reliability and operation but only if:

- Installation is carried out in strict accordance with this manual.
- The readjustment or repair operations are carried out by persons and/or representative companies authorized by ORDISI, S.A.
- The electrical installation of the corresponding premises complies with the C.E.I. (U.N.E.) and with the Low Voltage Regulation.
- The equipment is operated in accordance with the instructions in this manual.
- The accessories and consumables that must be replaced periodically are original accessories and consumables and/or approved and authorized by ORDISI, S.A.

4. - PARTS DESCRIPTION

1. Handle
2. LED's module
3. Headlight
4. Fork
5. Switch
6. Moving arm
7. Horizontal arm
8. Vertical arm
9. Trim


Fig. 1 – Parts of the IG-65C lamp.

5. - INSTALLATION

5.1. - UNPACKING

If upon receiving the equipment you notice that the packaging is damaged, check that the device has not been affected. If so, notify the sender immediately.

5.2. - MOUNTING

5.2.1. - INSTALLATION OF THE MOUNTING PLATE

The assembly and installation of the IG-65C series lamp must take place in accordance with the Assembly Instructions in this manual and be carried out only by ORDISI personnel or by authorized assembly companies. Incorrectly mounting it can cause it to come loose during operation, potentially causing serious injuries to both the patient and healthcare personnel.

The type of anchorage described is designed for reinforced concrete ceilings executed according to DIN 1045. In other types of ceilings, the fastening must be projected on site.

The lamp has been designed for a ceiling height of about 275 cm. For different heights, adapt the vertical arm, see point 5.2.2.


Fig. 2. Dimensions of the holes in the ceiling.

To fix the anchor plate to the ceiling, mark the four holes using the supplied template. Four Ø10 holes and a depth of 80mm. Once the holes are made, introduce the plug No. 10 and attach the mounting plate No. 11 using washers No. 12 and screws No. 13. Level the plate.


Fig.3 Components of the mounting plate.


Fig.4 Template for making the holes.

5.2.2. - ARM INSTALLATION

The lamp is designed for a ceiling with a height of about 275 cm high. The height of the lamp will depend on the work area: the headlight should remain at a distance of approximately 100 cm from its workspace. If necessary, you can adjust the height cutting the vertical arm No. 8.

$$H_{\text{vertical arm}} = H_{\text{ceiling}} - (92 + H_{\text{head light to floor}})$$

Measurements in cm


Fig.5- Dimensions of the lamp.

Once mounted the ceiling plate, assemble the horizontal arm No.7 to vertical arm No. 8, join the electrical connectors No. 14 and No. 15 and fix by washers No. 16 and screws No. 17 (See Fig. 6).


Fig. 6 – Assembly of the arms.

After that, assemble the headlight No.3 with the arm No.6. Hold the fork through the screws No.13, and adjust it until you can move the headlight with gentle movements (See Fig. 7).


Fig. 7– Fork and headlight assembly

5.2.3. - ELECTRICAL INSTALLATION

The electrical installation of the room where the equipment is going to be used must comply with the UNE standards (CEI 60364-710), as well as with the Low Voltage Regulations. The electrical installation in the premises must include a fuse, as well as a network switch for the simultaneous separation of the lamp from the supply group.

It is essential that the lamp is properly grounded. Connect the network to the strip located on the anchor plate. Admitted voltage of 100-240V AC.

Once the connections have been made, check that the equipment is working properly and mount trim no. 9.

WARNING: To avoid the risk of electric shock, this equipment must only be connected to a mains supply with a protective earth connection.

5.3. - TURNING ON THE LAMP

1. Turn on the switch of the room. A red light should appear in the keypad when the lamp is connected.
2. Turn on the lamp, using the keypad No.5, and adjust the light intensity with the + and – buttons. (See Fig. 8).
3. To change the position of the lamp, do it by means of handle No. 1. Move the lamp gently, avoiding sudden movements.
4. The arms of the lamp are adjusted at the factory. If it is necessary to readjust, see point 7.2


Fig. 8 – Keypad

- | |
|---|
| <p>A – On / Off button
 B – On / Off Led indicator
 C – Luminous intensity adjustment buttons
 D – Luminous intensity Leds indicators</p> |
|---|

5.4. - RECOMENDATIONS FOR USE

Move the lamp gently, avoiding sudden movements.

If the lamp is not in use do not leave it turned on, since apart from wasting electricity you will reduce the useful life of the LEDs.

This equipment has been designed for its use in hospital environments. Therefore, the use of the lamp is reserved for specialized medical staff, though not special training is required.

6. - CLEANING

Before carrying out any operation, make sure that the lamp is off.

The lamp must be cleaned regularly with a cloth moistened with water, taking care that it does not penetrate the interior of the headlight. For persistent stains; e.g. blood, disinfectants, etc. you can use a little neutral soap. Do not use under any circumstances products that contain acetone.

Dry the lamp well before using it. If liquids have penetrated inside, you must wait the precise time until they evaporate.

The handle located in the headlight is sterilisable by means of steam at a temperature of 121°C with an exposure time of 20 minutes.

7. - ADJUSTMENTS AND MAINTENANCE

The electrical diagrams can be requested by the maintenance services to help them in their work. Please, if necessary, contact your distributor and provide him with the model and serial number so that he can provide you with these diagrams.

IMPORTANT: Before carrying out any repair/maintenance operation on the lamp, you must first make sure that the equipment is isolated from the electrical network. To do this, the plug of the lamp must be unplugged from the mains socket.

7.1. - PERIODIC CHECKS

Periodically, following the Hospital Protocol, the Maintenance Services must carry out the following verifications:

- 1) Check the general operation of the lamp.
- 2) Check mobility of the arm and dome. See point 7.2.
- 3) Check the continuity of the earth connection and its correct connection.
- 4) Check the general condition (cleanliness, etc.).

7.2. - ADJUSTMENT OF THE ARM

Check the adjustment of the ball joints, both of the swinging arm and of the fork and headlight. If necessary, adjust them: they should move effortlessly and stay in the desired position. Excessive tightening of the screws can block the movement. Adjust by means of screws no. 18, 19, 20 or 21.


Fig. 9a – Arm adjustment.


Fig. 9b – Headlight and fork adjustment.

7.3. - CHANGE THE FUSES

If the switch is in the "ON" position and the lamp connected to the mains, but does not work, check the condition of the fuse by following these instructions:

- 1) Remove trim No. 9.
- 2) Pull the fuse holder and remove the fuse.
- 3) If it is damaged, replace the fuse with one of the same characteristics (250VF2A).
- 4) Reassemble the fuse holder.
- 5) Mount trim No.9.

8.- ENVIRONMENTAL CONDITIONS

OPERATION		
	Min.	Max.
Temperature	10°C	30°C
Environmental humidity	30%	75%
TRANSPORT / STORAGE		
	Min.	Max.
Temperature	-10°C	50°C
Environmental humidity	20%	90%

9.- WASTE

The IG-65C lamp does not contain any dangerous product. All the components should be properly recycled at the end of their lifetime.

The electrical material must be recycled according to the proper procedure. The rest material must be removed according to their composition.

10. - TECHNICAL CHARACTERISTICS

Luminous intensity at 100 cm	65.000 lux ($\pm 10\%$)
Illumination intensity adjustment	20% - 100%
Light source	LED
Average LEDs Lifetime	60.000 h
Colour temperature	4.500 °K ($\pm 5\%$)
Electrical connection	100-240V / 47-63Hz
Power consumption	42 VA / 17 W
Headlight diameter	$\varnothing$ 300 mm
Field diameter d_{10}	$\varnothing$ 250 cm
Field diameter d_{50}	$\varnothing$ 130 cm
Colour	RAL 9002
Approximate lamp Weight	12 Kg


Fig. 10 – General dimensions

11. - OTHER MODELS AND WEIGHTS

MODEL	WEIGHT
IG-65M MOVABLE LAMP	13,5 kg
IG-65MB MOVABLE LAMP WITH EMERGENCY EQUIPMENT	19,5 kg
IG-65W WALL LAMP	8,5 kg
IG-65C CEILING LAMP	12 kg
IG-65DC CEILING LAMP WITH DOUBLE ARM	21 kg


Fig. 11 – Despiece de la cúpula / Exploded view of the headlight


Fig. 12 – Despiece de brazos y anclaje a techo / Exploded view of the arms and ceiling mount.