

AQUASOFT* GLOVES: A NEW LEVEL OF COMFORT, PROTECTION AND VALUE

AQUASOFT* COMFORT

AQUASOFT* Gloves are comprised of a unique nitrile formulation that provides users with both comfort and tactile sensitivity.

Like our entire nitrile glove family, AQUASOFT* Nitrile Exam Gloves are easy to don, with a silky, snug fit.

Key features include:

- Preferred tactile sensitivity
- Textured fingertips
- Beaded cuffs
- Ambidextrous.

They're ideal for use in applications needing superior tactile sensitivity and tasks which require the quality of nitrile, with added affordability.

AQUASOFT* PROTECTION

HALYARD* AQUASOFT* Exam Gloves provide you with the confidence of increased protection by providing resistance to the following chemotherapy drugs:



The following chemotherapy drugs and concentrations had NO breakthrough detected up to 240 minutes:

- | | |
|--|--|
| Arsenic Trioxide (Trisenox) (1.0 mg/ml) | Fluorouracil (50.0 mg/ml) |
| Bleomycin (15.0 mg/ml) | Gemcitabine (Gemzar) (38.0 mg/ml) |
| Bortezomib (Velcade) (1.0 mg/ml) | Idarubicin (1.0 mg/ml) |
| Busulfan (6.0 mg/ml) | Ifosfamide (IFEX) (50.0 mg/ml) |
| Carboplatin (10.0 mg/ml) | Irinotecan (20.0 mg/ml) |
| Chloroquine (50.0 mg/ml) | Mechlorethamine HCl (1.0 mg/ml) |
| Cisplatin (1.0 mg/ml) | Melphalan (5.0 mg/ml) |
| Cyclophosphamide (Cytoxan) (20.0 mg/ml) | Methotrexate (25.0 mg/ml) |
| Cyclosporine (100.0 mg/ml) | Mitomycin C (0.5 mg/ml) |
| Cytarabine (100.0 mg/ml) | Mitoxantrone (2.0 mg/ml) |
| Dacarbazine (DTIC) (10.0 mg/ml) | Oxaliplatin (2.0 mg/ml) |
| Daunorubicin (5.0 mg/ml) | Paclitaxel (Taxol) (6.0 mg/ml) |
| Docetaxel (10.0 mg/ml) | Paraplatin (10.0 mg/ml) |
| Doxorubicin HCl (2.0 mg/ml) | Retrovir (10.0 mg/ml) |
| Epirubicin (Ellence) (2.0 mg/ml) | Rituximab (Rituxan) (10.0 mg/ml) |
| Etoposide (Toposar) (20.0 mg/ml) | Topotecan HCl (1.0 mg/ml) |
| Fludarabine (25.0 mg/ml) | Vincristine Sulfate (Oncovin) (1.0 mg/ml) |

Carmustine (3.3 mg/ml): Not for use with Carmustine

ThioTEPA (10 mg/ml): Not for use with ThioTEPA

WARNING: Not for use with Carmustine, ThioTEPA

The AQUASOFT* Powder-Free Nitrile Exam Glove material was tested with the above chemotherapy drug concentrations in accordance with ASTM D 6978, Standard Practice for Assessment Of resistance of Medical Gloves to Permeation by Chemotherapy Drugs.¹

AQUASOFT* VALUE

AQUASOFT* Gloves are thinner than typical nitrile gloves. That means they take up less storage space – each inner box contains 300 gloves – and they generate less waste.



67% less storage space required with 300 count box†

† Compared to standard 100 count boxes, 300 count AQUASOFT* Exam Gloves saves space.

ASTM D6319 Standard Specification for Nitrile Examination Gloves for Medical Application.
ISO-10993-10 PRIMARY Skin Irritation & Sensitization.
ASTM F739 Standard Test Method for Permeation of Liquids and Gases through Protective Clothing Materials under Conditions of Continuous Contact.
ASTM F1671 Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System.



Infection Prevention Product Solutions
Halyard* Education Foundation
Tools & Best Practice
Commitment to Excellence and Sustainability
Experienced Customer Support
Knowledgeable Sales Force
Innovation & Technology

www.halyardhealth.com.au



HALYARD* AQUASOFT* NITRILE EXAM GLOVES

| Code | Description | Size | Packaging |
|-------|-------------------------------|---------|---------------------------------|
| 47634 | AQUASOFT* Nitrile Exam Gloves | X-Small | 300 eaches/box 10 boxes/case |
| 47635 | AQUASOFT* Nitrile Exam Gloves | Small | 300 eaches/box 10 boxes/case |
| 47636 | AQUASOFT* Nitrile Exam Gloves | Medium | 300 eaches/box 10 boxes/case |
| 47637 | AQUASOFT* Nitrile Exam Gloves | Large | 300 eaches/box 10 boxes/case |
| 47638 | AQUASOFT* Nitrile Exam Gloves | X-Large | 250 eaches/box 10 boxes/case |

Physical Properties²

Before Aging

Tensile Strength: **23MPa**
Ultimate Elongation: **500%**

After Aging

Tensile Strength: **28MPa**
Ultimate Elongation: **450%**

Target Specifications²

| Gauge Thickness Measurements | MM | MIL |
|------------------------------|------------|---------------|
| Middle Finger: | 0.08 | 3.1 |
| Palm: | 0.07 | 2.8 |
| Cuff: | 0.06 | 2.4 |
| Average Length: | 242 | (9.5") |

Quality Standards²

Exceeds current ASTM Standards for critical defects, tensile strength and puncture resistance.
Passes primary skin irritation testing in accordance with ISO 10993-10(2010).
Manufactured in accordance with Quality System ISO 13485.

AQL – ASTM D 5151²

Critical Defects Level **<1.0 AQL[†]**
† Acceptable Quality Level

CAUTION: The testing conditions used are intended to approximate the worst case conditions for clinical use. Testing was conducted on single layer glove material. It is the users' responsibility to determine the applicability of these gloves for their intended use with chemotherapy drugs.

Storage Recommendations: Store in a cool, dry place. Open box should be shielded from exposure to direct sunlight, intense artificial light, x-ray machines and other sources of ozone.



Halyard Australia Pty. Limited
Phone: 1800 664 227

www.halyardhealth.com.au

ALWAYS READ THE LABEL AND FOLLOW THE DIRECTIONS FOR USE.

Halyard Health makes no representation that its products described in this catalogue will provide complete protection from exposure to contaminants, fluids or micro-organisms or will eliminate the risk of contracting infectious diseases.

1. CLAIM-07590.
2. C11406.