

HydroTac®: Case Studies of Use

HARTMANN

Pam Spruce - Clinical Director TVRE Consulting.
Lindsey Bullough - Clinical Nurse Specialist – Tissue Viability, Wrightington, Wigan and Leigh Foundation Trust

Introduction

Foam dressings are now reported to be the most commonly used product in wound management^{1,2}. Non adhesive foam dressings are used less frequently than those with an adhesive border³, however there is still an indication for use where the patient may have a preference or a clinical indication (e.g vulnerable peri-wound skin) for this type of product.

Shaped foam dressings are also frequently used to improve conformability on areas of the body, which may be difficult to dress.

HydroTac®

HydroTac® is a NEW, unique foam dressing with AquaClear Technology that provides a combination of absorption and moisture donation. It provides a moist wound environment, by absorbing exudate but can release moisture when applied to a dry wound.

The interface of the dressing is impregnated with a hydrogel (AquaClear Technology) which prevents it from adhering to the wound bed, and facilitates painfree dressing removal.

In a recent evaluation undertaken on 20 patients by both nurses and podiatrists, HydroTac® was used on a range of acute and chronic wounds where a non-adhesive or shaped dressing was required. In 85% (n=17) of the patients, the wounds progressed with 20% (n=4) healing within the four-week evaluation period.

The results also indicated that:-

- In 100% of dressing changes (n=93) the dressing was easy to apply and remove.
- In 100% of applications (n=93) the dressing conformed well to the wound.
- In 100% of responses (n=93) the patients reported that the dressing was comfortable during wear and painless on removal.

HydroTac® was reported to manage exudate effectively, with dressing changes being undertaken every 3 days in 57% (n=47) of procedures, alternate days in 29% (n=27) and 5-7days in 20% (n=19). The peri-wound skin condition also improved in 55% of patients (n=11) where the tissue was damaged at baseline.

Case Study 1

The patient was a 46 year old female, who had a medical history of ischaemic heart disease and heart failure with uncontrolled oedema in her legs. A blistered area appeared on the gaiter area of her left leg which measured 9cm², from which there was a small amount of exudate.

Because of the friable skin condition, a non adhesive foam dressing was indicated to provide a moist wound environment, absorb the excess exudate and protect from further contamination and the patient also complained that the wound was painful

The wound had been present for 1 week with no treatment before HydroTac® was applied, with a tubular cotton bandage to secure in place. After 3 dressing changes (7 days) the blistered area was fully epithelialised with no trauma to the surrounding skin.



Case Study 2

The patient was a 72 year old female who presented with a grade 2 pressure ulcer on her left heel. Although the wound area was small measuring 1.5cm², the wound bed contained both slough and granulation tissue and the peri-wound skin was macerated.

After a full assessment of the patient, a programme of care was implemented which included pressure relief and local wound management using HydroTac®. The concave shape dressing was used to ensure that the dressing interface was in contact with the wound bed, and was held in situ with a retention bandage.


The wound was re-assessed every 3 days, and HydroTac® was reapplied. After 18 days the wound was fully healed.



Case Study 3

The patient was a 52 year old female who was admitted to hospital with a fractured neck of femur on the left leg. On examination, it was observed that she had a venous ulcer present, which measured 40cm². Although the wound bed was clean and granulating, there was a moderate amount of exudate and the peri-wound skin was macerated. The patient also complained of pain in the wound.

Because of the peri-wound skin condition, HydroTac was used on the ulcer, with a wool/retention bandage to secure in place. The wound was reassessed every 3-4 days, and HydroTac was re-applied. The patient became pain free and the skin condition returned to normal as the wound progressed. After 22 days of treatment the wound was reported to be healed.



Conclusion

In a small evaluation, the HydroTac® dressing was observed to be effective in managing patients with a range of acute and chronic wounds, which required either a non adhesive, or shaped dressing. The outcome of the evaluation demonstrated that the dressing facilitated wound progression in a high number of patients, was easy to use and highly acceptable to both patients and clinicians. It was comfortable during wear and removal, and an improvement in peri-wound skin condition was observed.

1. Bianchi J, Gray D, Timmons J, Meaurio S. Do all foam dressings have the same efficacy in the treatment of chronic wounds? Wounds UK (2011) 7(1): 62-67
2. Carter K. Hydropolymer dressings in the management of wound exudate. Br J Community Nurs (2003) 8: suppl 10-6
3. Peach V. Evaluating adhesive foam wound care dressings in clinical practice. Wounds UK 2012, Vol 8, No 3 53-54.