

INADINE™

(PVP-I) NON ADHERENT DRESSING

Patient Guide

What is it?

INADINE™ PVP-I Non-Adherent Dressing is a topical wound Dressing impregnated with an ointment containing 10% povidone iodine. The Dressing also contains polyethylene glycol and purified water.


What is it for?

INADINE™ Dressing is designed to protect the wound, even if infected. INADINE™ Dressings can be used for the management of ulcerative wounds and may also be used for the prevention of infection in minor burns and minor traumatic skin loss injuries. In heavily infected wounds, it may also be used in conjunction with systemic antibiotics.

How to Apply:


Refer to the instructions for use. Select a size of INADINE™ Dressing that is slightly larger than the wound. Always wash your hands before and after Dressing change. Prepare the wound according to appropriate wound management protocol. Ensure skin surrounding the wound is dry. Remove the Dressing from the pouch.

Dressing preparation




1. If necessary, **INADINE™** Dressing can be cut using clean scissors or folded to fit the wound bed.
2. It is best to cut **INADINE™** Dressing whilst one or two backing papers are in place.

Secondary Dressing application



1. Cover **INADINE™** Dressing with a secondary Dressing such as TIELLE™. Choose the secondary Dressing based on the level of wound fluid.
- Apply any secondary cover or retention/ compression bandage.

Dressing application



1. Remove one backing paper.
2. Peel off the Dressing from the remaining backing paper and apply **INADINE™** Dressing directly to the wound.

Dressing change and removal



1. When the **INADINE™** Dressing colour fades this indicates loss of antiseptic efficacy and the Dressing should be changed. It may be changed up to two times daily in the initial phase or with highly infected wounds or wounds producing large amounts of wound fluid.
2. Gently remove the secondary Dressing.
3. Removal of **INADINE™** Dressing may be helped by gentle stretching of the Dressing at diagonally opposite corners.

How often to Apply:

When INADINE™ Dressing turns white this indicates loss of antiseptic efficacy and therefore the Dressing should be changed. It may be changed up to two times daily in the initial phase or with highly infected wounds or wounds producing large amounts of exudate.

NOTE: This document is intended for patient self-care as guided by their HCP and is not intended to replace the Instructions for Use.

NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. The instructions for use are intended for healthcare professionals.