

**Coated
VICRYL®**














Coated VICRYL* (polyglactin 910) Suture














Product Characteristics	Coated VICRYL* Suture
Physical Material	Polyglactin 910 (copolymer of glycolide and lactide)
Construction	Braided (Monofilament in sizes 9/0 and 10/0)
Absorbable / Non-absorbable	Absorbable
Coating (if applicable)	50% Polyglactin 370 and 50% Calcium Stearate
Colours	Undyed – beige and Dyed – violet
Available Size Range	3 through 8/0
Tensile Strength (Wound Support)	Medium Term (4 weeks)
Tensile Strength Retention in Tissue	6/0 and larger: 75% of original strength remains at 2 weeks 50% of original strength remains at 3 weeks 25% of original strength remains at 4 weeks 7/0 and smaller: 75% of original strength remains at 2 weeks 40% of original strength remains at 3 weeks
Absorption Time	Absorbed in 56-70 days, 63 days average
Frequent Uses	General soft tissue approximation and/or ligation, including use in Ophthalmic procedures

*Trademark











Coated VICRYL® Suture (Polyglactin 910) Braided Suture

Needle length and description		Length	Colour	Gauge						
○ 1/2 CIRCLE ETHIGUARD BLUNT POINT				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
SHB 	26mm	8x45cm control release	undyed	JB864						
CTB-2 	26mm	70cm	violet				D9754			
MHD-1 	31mm	75cm	violet		W9984					
CTB-1 	36mm	8x45cm control release	undyed			JB840	JB841†			
CTB-1 	36mm	70cm 70cm 90cm	violet undyed undyed	JB944	JB259	JB340†		JB947		
CTD-1 	36mm	75cm 90cm	violet violet		W9993		W9994	W9995		
CTD 	40mm	75cm 90cm	violet violet			W9988		W9997		
TD-3 	45mm	90cm	violet				W9998	W9999		
CTX-B 	48mm	8x45cm control release	undyed				JB724D†	JB725†		
CTX-B 	48mm	90cm 90cm	violet undyed				JB978†	JB977†		
BP-1 	65mm	4x90cm	violet					D8659		
















Coated VICRYL® Suture (Polyglactin 910) Braided Suture

Needle length and description		Length	Colour	Gauge						
○ 3/8 CIRCLE ETHIGUARD BLUNT POINT				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
	65mm	100cm	violet				W9391			
▼ COMPOUND CURVE PRECISION POINT REVERSE CUTTING				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
	16mm	45cm	undyed	W9436	J656G†					
▼ 3/8 CIRCLE PRECISION POINT REVERSE CUTTING				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
	10.5mm	45cm 45cm	undyed undyed	J488G	J489G W9500T	W9501T				
	11mm 11mm	45cm	violet		W9537					
	13mm	45cm	undyed		J492G	J493G†	J494G†			
	16mm	45cm	undyed			J500G†				
	19mm	45cm 45cm	undyed undyed				J496H† W9515T	J497H†		
	24mm	45cm	undyed					J683H†		
▼ 3/8 CIRCLE REVERSE CUTTING				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
	19mm	45cm 70cm	violet undyed			J421H†	J422H†	J393H† J423H†		
	24mm	70cm 70cm	violet undyed				J451H† J441H†	J442H†	J443H†	
	30mm	70cm	undyed					J588H†	J589H†	
▼ 1/2 CIRCLE PRECISION POINT REVERSE CUTTING				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
	8mm	45cm	undyed			J503G†	J504G†			
▼ 1/2 CIRCLE REVERSE CUTTING				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
	22mm	45cm control release	violet							












Coated VICRYL® Suture (Polyglactin 910) Braided Suture

Needle length and description		Length	Colour	Gauge						
▼ 1/2 CIRCLE REVERSE CUTTING				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
OS-4 	22mm	70cm 75cm	undyed violet		J693H	J694H	J695H†	W9287		
X-1 	22mm	70cm 70cm	violet undyed	J460H† J458H†						
CP-2 	26mm	70cm	undyed		J869H†	J870H†	J871H†			
▼ 1/2 CIRCLE REVERSE CUTTING				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
CP-1 	36mm	70cm 70cm 90cm	violet undyed violet	J471H	J466H† J472H†	J467H† J267H†	J468H† J268H† J474H			
OS-6 	36mm	70cm	undyed				J535H†			
CP 	40mm	70cm 90cm	violet violet			J485H				
OS-8 	40mm	70cm	undyed			J699H				
CPX 	48mm	70cm	undyed				J569H†			
▼ STRAIGHT CUTTING				4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	3 6
KS 	60mm	70cm 70cm 75cm	violet undyed undyed	J662H†	J663H† W9717	J664H				
▲ 3/8 CIRCLE CONVENTIONAL CUTTING P NEEDLE				4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	3 6
PC-26 	26mm	45cm 75cm	undyed undyed	W9520T W9580T	W9521T W9581T	W9522T W9582T				









Coated VICRYL® Suture (Polyglactin 910) Braided Suture

Needle length and description		Length	Colour	Gauge						
▲ 3/8 CIRCLE CONVENTIONAL CUTTING PRIME				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 PC-1	13mm	45cm	undyed	J833G	J834G†	J835G†				
 PC-3	16mm	45cm 45cm	undyed undyed		J844G† W9505T	W9506T				
 PC-5	19mm	45cm 45cm 75cm	undyed undyed undyed			J823G† W9510T W9570T	J824G† W9511T W9571T			
▲ 1/2 CIRCLE CONVENTIONAL CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 CX-1	22mm	45cm	violet				W9472^			
▲ 1/2 CIRCLE CONVENTIONAL CUTTING STERNUM				0 3.5	1 4	2 5	3 6	4 6	5 7	7 9
 CCS-1	60mm	70cm	undyed				J808T			
⊙ 3/8 CIRCLE TAPER POINT				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
 BV100-3	3.9mm	15cm	violet	V402G						
 BV130-5	6.5mm	15cm	violet		J401G					
 BV-1	9.3mm	30cm	violet		W9577^					
⊙ 1/2 CIRCLE TAPER POINT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 RB-2	13mm	45cm	undyed	W9981	W9982					
 TF	13mm	70cm	undyed		J433H†					
 RB-1	17mm	70cm 70cm	violet undyed		J303H†	J304H† J214H†	J305H†	J306H†		
 CT-3	22mm	70cm	violet				J327H		J329H	
 SH-1	22mm	70cm	violet			J310H†	J311H†			
 SH	26mm	8x45cm control release	violet undyed				J774D† J864D†	J775D†		
 SH	26mm	70cm control release	violet				J784G			










Coated VICRYL® Suture (Polyglactin 910) Braided Suture

Needle length and description		Length	Colour	Gauge						
☉ 1/2 CIRCLE TAPER POINT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 SH	26mm	70cm 70cm 75cm	violet undyed violet			J315H†	J316H† J416H† W9120	J317H† J417H† W9121	J318H	
 SH	26mm 26mm	90cm	violet				J527H†			
 CT-2	27mm	8x45cm control release	violet					J726D	J727D	
☉ 1/2 CIRCLE TAPER POINT				4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	3 6
 CT-2	27mm	70cm 70cm	violet undyed		J332H† J232H†	J333H† J269H†	J334H† J270H†	J335H†		
 MH-1	31mm	75cm	violet					W9213		
 MO-5	31mm	75cm	violet					W9216		
 CT-1	36mm	8x45cm control release	undyed				J840D			
 CT-1	36mm	70cm 70cm 90cm 90cm	violet undyed violet undyed		J338H†	J339H† J259H† J345H† J945H†	J340H† J260H† J346H† J946H†	J341H† J261H† J347H†		
 MH	36mm	70cm 75cm	violet violet		J322H†	J323H†	J324H W9141			
 MO-4	36mm	70cm 75cm	violet undyed				J436H†	J437H†		
 CT	40mm	8x45cm control release	violet			J751D	J752D†	J753D†		




Coated VICRYL® Suture (Polyglactin 910) Braided Suture

Needle length and description		Length	Colour	Gauge						
1/2 CIRCLE TAPER POINT				4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	3 6
	CT 40mm	70cm 70cm 90cm 90cm 90cm	violet undyed violet undyed violet			J351H† J275H† J357H†	J352H† J358H† J958H†	J353H† J359H† J959H† W9431		
	MO-45 45mm	75cm	violet						W9246	
	CTX 48mm	70cm 90cm	violet violet			J363H† J369H†	J364H† J370H	J365H† J371H†		
	TP-1 65mm	4x70cm 70cm 135cm 135cm	violet violet undyed violet					J650G D5792	J649G J880T†	
	XLH 70mm	70cm	violet			J581G†		J583G		
1/2 CIRCLE TAPER POINT, 1/2 CIRCLE REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	CT-1 / CP-1 36mm 36mm	135cm	violet					J532H		
5/8 CIRCLE TAPER POINT				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
	UR-6 27mm	70cm	violet	D7973	J602H†	J603H†				
	UR-6 27mm 27mm	90cm	undyed							










Coated VICRYL® Suture (Polyglactin 910) Braided Suture

Needle length and description		Length	Colour	Gauge						
5/8 CIRCLE TAPER POINT				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
	36mm	70cm	violet		J375H†					
	40mm	70cm 70cm	violet violet		J381H W9165					
1/4 CIRCLE TAPERCUT				8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3
	5mm	25cm	violet	J408G						
1/2 CIRCLE TAPERCUT				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
	36mm	75cm 90cm 90cm	violet violet undyed	J516H J916H†	J517H†	J518H†	W9365 J519H†	W9357^		
1/2 CIRCLE TAPERCUT				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
	45mm	75cm	violet			W9376	W9377	W9378		
5/8 CIRCLE TAPERCUT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	10mm	30cm	violet	W9172				J532H		
J NEEDLE TAPERCUT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	31mm	75cm	violet							W9335
	22mm	75cm	violet				W9223			
	36mm	75cm	violet					W9220	W9221	

Coated VICRYL® Suture (Polyglactin 910) Braided Suture

Needle length and description		Length	Colour	Gauge						
▼ 3/8 CIRCLE MICROPOINT REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	G-6 8mm 8mm	45cm	violet	J551G						
◁ 1/4 CIRCLE MICROPOINT SPATULA				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	TG100-8 6.5mm 6.5mm	45cm	violet	J544G						
◁ 3/8 CIRCLE MICROPOINT SPATULA				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
	TG140-6 6.5mm 6.5mm	30cm	violet		J974G					
	TG140-8 6.5mm 6.5mm	30cm	violet			W9561				
	TG140-8 6.5mm 6.5mm	30cm 45cm	violet violet	V549G	J548G	J566G J546G				
◁ 1/2 CIRCLE MICROPOINT SPATULA				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
	TG160-6 5.5mm 5.5mm	30cm	violet		J975G					
◁ 1/2 CIRCLE MICROPOINT SPATULA				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
	TG160-8 5.5mm 5.5mm	45cm	violet			J576G				
◁ 1/4 CIRCLE SABRELOC CONVENTIONAL SPATULA				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
	S-29 8mm 8mm	30cm 45cm	violet violet				J556G J555G			
◁ 1/2 CIRCLE SABRELOC CONVENTIONAL SPATULA				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
	S-28 8mm 8mm	45cm	violet				J562G			
◁ 1/4 CIRCLE SABRELOC SPATULA				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
	S-14 8mm 8mm	20cm 45cm 45cm	violet violet undyed				W9562 J570G J670G	J571G J671G		
◁ 1/2 CIRCLE SABRELOC SPATULA				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
	S-22 8mm 8mm	45cm 45cm	violet undyed				W9566 W9756	W9567		
◁ 1/4 CIRCLE SABRELOC CENTRE POINT SPATULA				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
	RD-1 8mm 8mm	30cm	violet				J446G			

Coated VICRYL® Suture (Polyglactin 910) Braided Suture

Needle length and description		Length	Colour	Gauge						
 1/4 CIRCLE ULTIMA SPATULA				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
CS90-6 	7mm	10cm	violet	V980G						
 3/8 CIRCLE ULTIMA SPATULA				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
CS140-6 	6.5mm	10cm	violet	V960G						
 3/8 CIRCLE ULTIMA SPATULA				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
CS140-6 	6.5mm	10cm	violet	V950G						
CS140-8 	6.75mm 6.75mm	20cm	violet		J547G					
 1/2 CIRCLE ULTIMA SPATULA				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
CS160-6 	5.5mm 5.5mm	30cm	violet	V448G						
NON-NEEDED				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
LIGIPAK*	135cm 135cm 360cm	violet undyed undyed	J283G		J205G† J285G†	J206G† J286G†	J207G†	J208G†		
NON-NEEDED	6x45cm	violet					J106T†			
NON-NEEDED	12x45cm	violet		J103T†	J104T†	J105T†				
NON-NEEDED	135cm 135cm 150cm 150cm	violet undyed violet undyed		W9023 J643H*	J614H J606H	J615H† J607H† W9025	J616H†	J617H†	J618H†	W9027

