

**ETHILON®**
















## ETHILON\* Nylon Suture











<b>Product Characteristics</b>	<b>ETHILON* Suture</b>
Physical Material	Polyamide polymer and Nylon 6
Construction	Monofilament
Absorbable / Non-absorbable	Non-absorbable
Coating (if applicable)	N/A
Colours	Clear, black and green
Available Size Range	2 through 11/0
Tensile Strength (Wound Support)	Long term wound support
Tensile Strength Retention in Tissue	Non-absorbable, but loses about 10-15% strength per year
Absorption Time	Degrades at a rate of about 10-15% per year
Frequent Uses	Skin closure; retention; Plastic surgery; Ophthalmology and Microsurgery

\*Trademark


















# ETHILON Nylon

Needle length and description		Length	Colour	Gauge						
○ 1/2 CIRCLE ETHIGUARD BLUNT POINT				0 3.5	1 4	2 5	3 6	4 6	5 7	7 9
	36mm	100cm	black		<b>B020G</b>					
	48mm	150cm loop	black		<b>W729</b>					
	65mm	240cm loop	black							
▼ 3/8 CIRCLE PRECISION POINT REVERSE CUTTING				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
	11mm	45cm 45cm	black black undyed	<b>1696G</b>	<b>697G</b> <b>W1610T</b> <b>689G</b>					
	13mm	45cm 45cm 45cm	black black undyed			<b>698G</b> <b>W1611T</b> <b>690G</b>	<b>699G</b> <b>691G</b>			
 PLIABILISED	13mm	45cm	black		<b>1698G</b>					
	16mm	45cm	black		<b>1665G</b>	<b>1668G</b>		<b>W1625BH</b>		
	19mm	45cm 75cm 75cm	undyed black undyed				<b>1611G</b>	<b>W1685BH</b>	<b>D9635</b>	
 PLIABILISED	19mm	45cm	black			<b>1667G</b> <b>1666G</b>				
	24mm	45cm 75cm	black undyed				<b>1663H</b> <b>D9782</b>			
▼ 3/8 CIRCLE REVERSE CUTTING				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
	13mm	45cm	black		<b>667G</b>	<b>668G</b>				
	19mm	45cm	black				<b>662G</b>	<b>653G</b>		
 PLIABILISED	19mm	45cm	black			<b>661G</b>				











# ETHILON Nylon

Needle length and description		Length	Colour	Gauge						
<b>▼ 3/8 CIRCLE REVERSE CUTTING</b>				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
 FS-1	24mm	45cm 75cm	black black				<b>1629H</b>	<b>663G</b> <b>669H</b>		
 FS	26mm	45cm	black						<b>664G</b>	
 FSL	30mm	75cm	black				<b>1670BH*</b>	<b>1671H</b>		
 FSLX	40mm	75cm	black					<b>1673BH</b>	<b>1674H</b>	
<b>▼ 3/8 CIRCLE REVERSE CUTTING</b>				0 3.5	1 4	2 5	3 6	4 6	5 7	7 9
 LR	75mm	50cm	black			<b>460T</b>				
 LR	75mm 75mm	75cm	black		<b>489T</b>					
<b>▼ 1/2 CIRCLE PRECISION POINT REVERSE CUTTING</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 PS-4	16mm	45cm 45cm	black undyed			<b>1662G</b> <b>1603G</b>				
<b>▼ 1/2 CIRCLE REVERSE CUTTING</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 X-1 PLIABILISED	22mm	45cm	black				<b>642G</b>			
<b>▼ 1/2 CIRCLE REVERSE CUTTING</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 CP	40mm	100cm 150cm loop	black black						<b>W737</b> <b>D10109</b>	
 CPX	48mm	100cm	black							















# ETHILON Nylon

Needle length and description		Length	Colour	Gauge						
<b>▼ STRAIGHT CUTTING</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 KS	60mm	75cm	black			<b>626H*</b>	<b>627H</b>	<b>628H</b>		
<b>▲ 3/8 CIRCLE CONVENTIONAL CUTTING PRIME</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 PC-1	13mm	45cm	black	<b>1856G</b>						
<b>PLIABILISED</b>  PC-1	13mm	45cm	black							
 PC-3	16mm	45cm	black	<b>W1615T</b>	<b>1865H</b>					
<b>PLIABILISED</b>  PC-3	16mm	45cm	black							
 PC-5	19mm	45cm	black			<b>1894G</b>	<b>1893G</b>			
<b>PLIABILISED</b>  PC-5	19mm	45cm	black							
<b>▲ 3/8 CIRCLE CONVENTIONAL CUTTING</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 CFS-3	16mm	35cm	black							
<b>▲ 3/8 CIRCLE CONVENTIONAL CUTTING</b>				3	0 3.5	1 4	2 5	3 6	4 6	5 7
 CKP-3	90mm	2x50cm	black			<b>W2797</b>				
<b>WITH 40MM RETENTION TUBING</b>										
<b>● 3/8 CIRCLE TAPER POINT</b>				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
 BV130-3	3.8mm	12cm	black		<b>2820G</b>	<b>2819G</b>				
 BV75-3	3.8mm	12cm	black		<b>2870G</b>					
 BV50-3	4.1mm	13cm 13cm	black black	<b>2881G</b> <b>W2881</b>						
 BV130-4	5mm	12cm 15cm	black black		<b>2814G</b>	<b>2813G</b>	<b>2815G</b>			
 BV100-4	5.1mm	15cm	black		<b>2830G*</b>	<b>2829G</b>				
 BV75-4	5.1mm	12cm	black		<b>2850G*</b>					
 BV130-5	6.5mm	12cm 15cm	black black		<b>2810G</b>	<b>2809G</b>	<b>2808G</b>			
 BV-1	9.3mm 9.3mm	45cm	black				<b>W8170</b>			















# ETHILON Nylon

Needle length and description		Length	Colour	Gauge						
<b>1/2 CIRCLE TAPER POINT</b>				0 3.5	1 4	2 5	3 6	4 6	5 7	7 9
 MH-1	31mm 31mm	75cm 100cm	black black	<b>4784G</b>	<b>W745</b>					
 MH-1	31mm	150cm loop	black	<b>W739</b>						
 CT-1	36mm	150cm loop	black		<b>W743</b>					
 CT	40mm	100cm 100cm	black black		<b>4999G</b> <b>W747</b>					
<b>1/2 CIRCLE TAPER POINT</b>				0 3.5	1 4	2 5	3 6	4 6	5 7	7 9
 CT	40mm	150cm loop	black	<b>W740G</b>						
 CT	40mm	200cm loop	black	<b>W740G</b>						
 CTX	48mm	100cm	black		<b>4985G</b>	<b>4986G</b>				
 CTX	48mm	150cm loop	black	<b>W741</b>	<b>W748</b>					
 CTX	48mm	200cm loop	black		<b>W768</b>					
 TP-1	65mm	120cm loop	black	<b>L880G</b>						

# ETHILON Nylon

Needle length and description		Length	Colour	Gauge						
<b>1/2 CIRCLE TAPER POINT</b>				0 3.5	1 4	2 5	3 6	4 6	5 7	7 9
	65mm	150cm	black			<b>825G<sup>^</sup></b>				
<b>1/4 CIRCLE TAPERCUT</b>				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
V100-3 	3.5mm	15cm	black		<b>2894G</b>					
V75-3 	3.5mm	15cm	black	<b>2892G</b>						
V130-3 	3.8mm	15cm	black			<b>2897G</b>				
<b>1/4 CIRCLE MICROPOINT SPATULA</b>				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
TG100-8 	6.5mm 6.5mm	30cm	black			<b>1714G</b>				
<b>3/8 CIRCLE MICROPOINT SPATULA</b>				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
VAS100-4 	5.1mm	15cm	black			<b>2890G</b>				
TG140-6 	6.5mm 6.5mm	30cm	black	<b>W1780</b>	<b>W1770</b>	<b>W1769</b>				
TG140-8 	6.5mm 6.5mm	30cm	black		<b>7718G</b>					
<b>1/2 CIRCLE MICROPOINT SPATULA</b>				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
TG160-6 	5.5mm 5.5mm	15cm 30cm	black black		<b>7757G</b> <b>W1756</b>					
TG175-8 	7mm 7mm	30cm	black			<b>1716G</b>				
<b>1/4 CIRCLE SABRELOC SPATULA</b>				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
S-24 	8mm 8mm	45cm	black							<b>7731G</b>
<b>1/4 CIRCLE SABRELOC CENTRE POINT SPATULA</b>				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
RD-1 	8mm 8mm	45cm	black							<b>749G</b>
<b>1/4 CIRCLE ULTIMA SPATULA</b>				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
CS100-8 	6.7mm 6.7mm	30cm	black				<b>9039G</b>			
CS90-6 	7mm	10cm	black		<b>9800G</b>					

# ETHILON Nylon

Needle length and description		Length	Colour	Gauge						
 <b>3/8 CIRCLE ULTIMA SPATULA</b>				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
CS140-6 	6.5mm	10cm	black		<b>9060G</b>					
CS140-6 	6.5mm 6.5mm	30cm	black		<b>9003G</b>					
CS140-8 	6.75mm 6.75mm	30cm	black		<b>9033G</b>					
 <b>1/2 CIRCLE ULTIMA SPATULA</b>				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
CSM-6 	5mm	30cm	black		<b>9005G</b>					
 <b>1/2 CIRCLE ULTIMA SPATULA</b>				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
CS160-6 	5.5mm	10cm	black		<b>9061G</b>					
CS160-6 	5.5mm 5.5mm	15cm	black		<b>9001G</b>					
CS160-6 	5.5mm	20cm	black	<b>9016G</b>						
CS160-6 	5.5mm 5.5mm	30cm	black		<b>9000G</b>					
CS160-8 	5.5mm 5.5mm	30cm	black		<b>9030G</b>					
CS175-6 	5.5mm 5.5mm	30cm	black		<b>9002G</b>					
CS175-8 	7mm 7mm	30cm	black		<b>9032G</b>	<b>9036G</b>	<b>9038G</b>			
<b>NON-NEEDED</b>				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
NON-NEEDED		2x100cm	black				<b>F6365G</b>			