

MONOCRYL®
















MONOCRYL* (poliglecaprone 25) Suture

Product Characteristics	MONOCRYL* Suture
Physical Material	Poliglecaprone 25 (copolymer of glycolide and epsilon-caprolactone)
Construction	Monofilament
Absorbable / Non-absorbable	Absorbable
Coating (if applicable)	N/A
Colours	Undyed – beige and Dyed – violet
Available Size Range	2 through 6/0
Tensile Strength (Wound Support)	Medium Term (3-4 weeks)
Tensile Strength Retention in Tissue	Undyed - 50-60% of original strength remains at 1 week - 20-30% of original strength remains at 2 weeks Dyed - 60-70% of original strength remains at 1 week - 30-40% of original strength remains at 2 weeks
Absorption Time	Completely absorbed in 91-119 days














*Trademark

MONOCRYL®






(Poliglecaprone 25) Monofilament Suture

Needle length and description		Length	Colour	Gauge						
○ 1/2 CIRCLE ETHIGUARD BLUNT POINT				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
	26mm	70cm	undyed	YB416						
	36mm	70cm 90cm	violet undyed	YB944		YB946	YB947			
	48mm	90cm	undyed			YB978H				
	70mm	90cm	violet				w3709			
▼ 3/8 CIRCLE PRECISION POINT REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	10.5mm	45cm	undyed	Y489G†	Y490G†					
	13mm	45cm	undyed	Y492G†	Y493G†	Y494G†				
	19mm	45cm 70cm	undyed undyed			Y496G† Y426H†	Y427H†			
	24mm	70cm 90cm	undyed undyed			Y935H† D8797	Y936H†			
	26mm	70cm	undyed				W3213			
▼ 3/8 CIRCLE REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	26mm	70cm	undyed					W3327		
▼ 1/2 CIRCLE PRECISION POINT REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	8mm	45cm	undyed			D8539				
	10.5mm	45cm	undyed	Y510G						
▼ 1/2 CIRCLE REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	36mm	70cm	undyed						Y267H	

Coated MONOCRYL® (Poliglecaprone 25) Monofilament Suture

Needle length and description		Length	Colour	Gauge						
▼ STRAIGHT CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 KS	60mm	70cm	undyed				Y523H†			
▲ 3/8 CIRCLE CONVENTIONAL CUTTING PRIME				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 PC-1	13mm	45cm	undyed	Y833G	Y834G†	Y835G†				
 PC-3	16mm	45cm	undyed		Y844G†	Y845G†				
 PC-5	19mm	45cm	undyed		Y822G†	Y823G†				
◎ 1/2 CIRCLE TAPER POINT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 TF	13mm	70cm	undyed	Y432H						
 RB-1	17mm	70cm 70cm	violet undyed		Y303H	Y304H†	Y305H Y215H			
 JB	26mm	70cm	violet					W3665		
 SH	26mm	70cm 70cm	violet undyed			Y315H†	Y316H†	Y317H† Y417H†		
 SH	26mm 26mm	90cm	violet				Y527H			
 CT-2	27mm	70cm	violet					Y333H†	Y334H	Y335H
 CT-1	36mm	90cm	violet					Y345H†	Y346H†	Y347H†
 MH	36mm	70cm 70cm	undyed violet					Y227H W3463		
 CT	40mm	90cm	violet					Y357H	Y358H†	Y359H

Coated MONOCRYL® (Poliglecaprone 25) Monofilament Suture

Needle length and description	Length	Colour	Gauge						
1/2 CIRCLE TAPER POINT			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 <p>CTX</p> <p>48mm</p>	90cm	violet					Y397H		Y399H
5/8 CIRCLE TAPER POINT			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 <p>UR-6</p> <p>27mm</p>	70cm	violet				Y604H†	Y605H†	Y606H†	
5/8 CIRCLE TAPER POINT			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 <p>UR-6</p> <p>27mm 27mm</p>	70cm 90cm	undyed violet				D9645 D9503	D8592		
1/2 CIRCLE TAPERCUT			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 <p>V-34</p> <p>36mm</p>	90cm	undyed					Y917H		
1/2 CIRCLE TAPER POINT, 1/2 CIRCLE REVERSE CUTTING			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 <p>CT-1 / CP-1</p> <p>36mm 36mm</p>	135cm	undyed					Y932H		
Needle length and description	Length	Colour	Gauge						
NON-NEEDEDLED			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
NON-NEEDEDLED	6X45cm	violet					Y909G	Y9110G	Y911G

