

**PDS™**
















## PDS™ II (polydioxanone) Suture

Product Characteristics	PDS™ II Suture
Physical Material	Polydioxanone
Construction	Monofilament
Absorbable / Non-absorbable	Absorbable
Coating (if applicable)	N/A
Colours	Undyed – clear and Dyed – violet
Available Size Range	2 through 9/0
Tensile Strength (Wound Support)	Long Term (8 weeks)
Tensile Strength Retention in Tissue	<p>4/0 and smaller:  60% of original strength remains at 2 weeks  40% of original strength remains at 4 weeks  35% of original strength remains at 6 weeks</p> <p>3/0 and larger:  80% of original strength remains at 2 weeks  70% of original strength remains at 4 weeks  60% of original strength remains at 6 weeks</p>
Absorption Time	Completely absorbed in 182-238 days
Frequent Uses	Abdominal and thoracic closure; subcutaneous tissue; Colon and Rectal surgery; Orthopaedic and Plastic Surgery















\*Trademark

# PDS™ II















## (Polydioxanone) Monofilament Suture

Needle length and description		Length	Colour	Gauge						
○ 1/2 CIRCLE ETHIGUARD BLUNT POINT				0 3.5	1 4	2 5	3 6	4 6	5 7	7 9
	36mm	70cm 90cm	violet violet	<b>ZB340</b>	<b>ZB347†</b>					
	45mm	90cm 90cm	violet violet		<b>W9965T</b>					
	48mm	90cm	violet	<b>ZB370†</b>	<b>ZB371†</b>					
	48mm	150cm loop	violet		<b>Z9967T</b>					
	65mm	244cm loop	violet		<b>D8892†</b>					
▼ 3/8 CIRCLE PRECISION POINT REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	11mm	45cm	undyed	<b>Z489G†</b>						
	13mm	45cm 45cm	violet undyed		<b>Z463G†</b> <b>Z493G†</b>	<b>Z494G†</b>				
	16mm	45cm	undyed			<b>W9873T</b>				
	19mm	45cm 45cm 45cm	violet undyed undyed		<b>Z513G†</b> <b>Z496G†</b> <b>W9874T</b>	<b>Z497G†</b>				
	24mm	45cm	undyed				<b>Z683G†</b>			
▼ 3/8 CIRCLE REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	19mm	45cm 70cm	undyed undyed		<b>Z421H†</b>	<b>W9950T</b> <b>Z422H†</b>	<b>Z423H†</b>			
	24mm	70cm	undyed			<b>Z441H†</b>	<b>Z442H†</b>	<b>Z443H†</b>		
	26mm	45cm	undyed				<b>W9957T</b>			

# PDS™ II (Polydioxanone) Monofilament Suture









Needle length and description		Length	Colour	Gauge						
<b>▼ 1/2 CIRCLE PRECISION POINT REVERSE CUTTING</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 P-2	8mm	45cm	undyed		Z503G <sup>†</sup>	Z504G				
 PS-4	16mm	45cm	undyed			Z507G <sup>†</sup>				
<b>▼ 1/2 CIRCLE REVERSE CUTTING</b>				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
 OS-4	22mm	70cm 70cm	violet violet				W9286T			
 CP-2	26mm	70cm	violet	Z968H	Z969H <sup>†</sup>	Z970H <sup>†</sup>				
 CP-1	36mm	70cm	violet		Z466H <sup>†</sup>	Z467H <sup>†</sup>	Z468H <sup>†</sup>			
 OS-6	36mm	70cm	undyed				Z535H <sup>†</sup>			
 CP	40mm	70cm	violet					Z195T		
<b>▼ STRAIGHT CUTTING</b>				4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	3 6
 KS	60mm	70cm	undyed	Z662H <sup>†</sup>	Z663H <sup>†</sup>					
<b>▲ 3/8 CIRCLE CONVENTIONAL CUTTING PRIME</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 PC-1	13mm	45cm	undyed	Z833G <sup>†</sup>						
 PC-3	16mm	45cm	undyed		Z844G <sup>†</sup>					
 PC-5	19mm	45cm 45cm	violet undyed			W9615T Z823G <sup>†</sup>	Z824G <sup>†</sup>			
 PC-25	26mm	45cm 70cm	undyed undyed			W9870T	W9625T			
<b>▲ 3/8 CIRCLE TAPER POINT</b>				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
 BV-1	9.3mm 9.3mm	60cm 75cm	violet violet	Z155H Z135H	Z117H <sup>†</sup>					
 BV	10.5mm 10.5mm	75cm	violet		Z123H <sup>†</sup>					

# PDS™ II (Polydioxanone) Monofilament Suture





Needle length and description		Length	Colour	Gauge						
<b>▲ 3/8 CIRCLE TAPER POINT</b>				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
	13mm 13mm	75cm	violet		<b>Z127H†</b>					
<b>● 1/2 CIRCLE TAPER POINT</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	13mm 13mm	75cm	violet	<b>Z149H†</b>						
	13mm	70cm	violet	<b>Z432H†</b>						
	14mm	70cm	violet					<b>Z114H†</b>		
	17mm	70cm	violet		<b>Z303H†</b>	<b>Z304H†</b>	<b>Z305H†</b>			
	17mm	75cm	violet		<b>Z320H†</b>					
	22mm	70cm	violet			<b>Z310H</b>				
	26mm	8x45cm control release	violet			<b>Z773D</b>	<b>Z774D</b>	<b>Z775D</b>		
	26mm	70cm	violet			<b>Z315H†</b>	<b>Z316H†</b>	<b>Z317H†</b>		
	26mm 26mm	90cm	violet				<b>D6264</b>			
	27mm	70cm	violet					<b>Z333H†</b>	<b>Z334H†</b>	
	22mm	70cm	violet							
<b>● 1/2 CIRCLE TAPER POINT</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	36mm	70cm 90cm	violet violet				<b>Z338H†</b>	<b>Z339H†</b>	<b>Z340H†</b> <b>Z346H†</b>	<b>Z341H†</b> <b>Z347H†</b>
	40mm	70cm 90cm	violet violet					<b>Z357H†</b>	<b>Z352H†</b> <b>Z358T†</b>	<b>Z353H†</b> <b>Z359T†</b>

See IFU for additional information, † Plus Code Available

# PDS™ II (Polydioxanone) Monofilament Suture

Needle length and description		Length	Colour	Gauge						
● 1/2 CIRCLE TAPER POINT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	MO-2 40mm	150cm loop	violet						<b>W9236T</b>	
	CTX 48mm	90cm	violet						<b>Z370T†</b>	<b>Z371T†</b>
	CTX 48mm	150cm loop	violet							<b>W9262T</b>
	TP-1 65mm	75cm 150cm loop 244cm loop	violet violet violet						<b>Z991G</b>	<b>D7819</b> <b>Z880G†</b>
	XLH 70mm	70cm 135cm	violet violet							<b>Z583G</b> <b>Z877G</b>
● 5/8 CIRCLE TAPER POINT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	UR-6 27mm	70cm	violet						<b>D7185</b>	
	MF-110 110mm	150cm loop	violet							
● STRAIGHT TAPER POINT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	ST-4 19mm 19mm	135cm	violet			<b>Z420G†</b>				

# PDS™ II (Polydioxanone) Monofilament Suture

Needle length and description		Length	Colour	Gauge						
<b>1/2 CIRCLE TAPERCUT</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	V-38 31mm	70cm	violet							<b>W9333T</b>
	V-34 35mm	90cm	violet							<b>W9352T</b>
	V-37 40mm	240cm loop	violet							<b>W9394</b>
<b>STRAIGHT TROCAR POINT</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	STP-10 254mm 254mm	70cm	violet					<b>Z997G</b>		

