

SILK



















Silk Suture













Product Characteristics	Silk Suture
Physical Material	Natural protein fibre of raw silk
Construction	Braided
Absorbable / Non-absorbable	Non-absorbable
Coating (if applicable)	Bees wax
Colours	Black and white
Available Size Range	5 through 9/0
Tensile Strength (Wound Support)	Medium to Long Term
Tensile Strength Retention in Tissue	Indefinite
Absorption Time	Non-absorbable
Frequent Uses	Most body tissues for ligating and suturing; General surgery; Ophthalmology and Plastic surgery

SILK









Braided Silk Suture

Needle length and description		Length	Colour	Gauge						
▼ 3/8 CIRCLE PRECISION POINT REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
P-1 	10.5mm	45cm	black	639G						
P-3 	13mm	45cm	black		640G					
PS-3 	16mm	45cm	black	1675H	1678H*					
PS-2 	19mm	45cm 45cm	black black		1676G	1677G	1679H			
PS-1 	24mm	45cm	black				1684G			
PS 	26mm	45cm	black			W2510T	W2511T	1588H		
PSL 	30mm	75cm	black					486H		
PSLX 	36mm	75cm	black					W2524T		
▼ 3/8 CIRCLE REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
C-2 	13mm	75cm	black							
C-2 	13mm 13mm	45cm	black	714G						
C-3 	13mm	30cm 45cm	black black	711G		W811				
C-3 	13mm 13mm	30cm 45cm	black undyed			735G W816				
FS-2 	19mm	45cm	black			683G				
FS-1 	24mm	45cm	black			629G	684G			
FS 	26mm	45cm	black					685G		
FSL 	30mm	45cm 75cm	black black				675H	677G 679H	680H	

















SILK Braided Silk Suture

Needle length and description		Length	Colour	Gauge						
▼ 3/8 CIRCLE REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 FSLX	36mm	75cm	black					1689H		
 FSLY	45mm	75cm	black						W782H	
 LR-60	60mm	75cm	black							W562H
▼ 1/2 CIRCLE PRECISION POINT REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 PS-4	16mm	45cm	black			1687G				
▼ 1/2 CIRCLE REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 X-1	22mm	8x45cm control release	black							
 X-1	22mm	45cm	black				632G	737G		
 CP-1	36mm	75cm	black					443H		
 CP	40mm	75cm	black							4796G
▼ STRAIGHT CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 KS	60mm	75cm	black				622H	623H	624H*	
 KS-75	75mm	100cm	black						W767	
▲ 3/8 CIRCLE CONVENTIONAL CUTTING PRIME				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 PC-3	16mm	45cm	black		W2500T	W2501T	W2502T			
▲ 3/8 CIRCLE CONVENTIONAL CUTTING PRIME				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 PC-3	16mm	4x45cm	3 black 1 white							














SILK Braided Silk Suture

Needle length and description		Length	Colour	Gauge						
▲ 3/8 CIRCLE CONVENTIONAL CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 SB-15	16mm	45cm	black	W529H	1925G					
 CFS-3	16mm	75cm 75cm	black undyed		W500H	W501H W508	W502H W509	W503H		
 CFS-3	16mm 16mm	38cm	black				W642			
 SB-25	26mm	45cm	black		1925G	W536H	W537H			
 CP-80	80mm	75cm	black				W260	W261		
▲ 3/8 CIRCLE CONVENTIONAL CUTTING				2/0 3	0 3.5	1 4	2 5	3 6	4 6	5 7
 CKP-3	90mm	2x50cm	black			W2793				
 CKP-3	90mm	100cm	black	W791	W792	W793G	W794			
▲ 1/2 CIRCLE CONVENTIONAL CUTTING				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
 CX-1	22mm	45cm	black	W577H						

SILK Braided Silk Suture

Needle length and description		Length	Colour	Gauge						
3/8 CIRCLE TAPER POINT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	BV-1 9.3mm 9.3mm	45cm	black							
	BV 10.5mm	75cm	black		W595					
	C-1 13mm 13mm	45cm	black	706G						
	BB 17mm	75cm	black		K880H	K881H				
	BB-2 22mm	75cm	black				W598H			
1/2 CIRCLE TAPER POINT				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
	TF-4 8mm	45cm	black		N272H					
	TF 13mm	8x45cm control release	black		C084D					
	TF 13mm	45cm	black	N266H						
	RB-1 17mm	75cm	black	K870H	K871H	K872H	K873H			
	RB-1 17mm 17mm	60cm	black		K571H	K572H				
	SH-2 20mm	75cm	black		W586H	W587H				
	SH-1 22mm	8x45cm control release	black			C003D				
	CT-3 22mm	50cm	black				413H			
	SH 26mm	8x45cm control release	black			C013D	C012D			
	SH 26mm	75cm 120cm	black black		K831H	K832H	K833H D8342	K834H	K835H	
1/2 CIRCLE TAPER POINT				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
	CT-2 27mm	45cm control release	black							

SILK Braided Silk Suture

Needle length and description		Length	Colour	Gauge						
☉ 1/2 CIRCLE TAPER POINT				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
 MO-5	31mm	75cm	black						W723H	
 CT-1	36mm	8x45cm control release	black					C021D		
 CT-1	36mm	75cm	black				423H	424H		
 MH	36mm	75cm 75cm	black black			W583	K843H	K844H	K845H	K846H
☽ 1/2 CIRCLE TAPERCUT				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
 V-5	17mm	75cm	black		K952H					
 V-7	26mm	75cm	black			K963G	K964H			
▼ 3/8 CIRCLE MICROPPOINT REVERSE CUTTING				8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3
 G-6	8mm 8mm	45cm	black		768G	769G				
 G-7	8mm 8mm	45cm	blue	W818						
 G-1	10.5mm 10.5mm	45cm	black			780G				
 G-3	13mm	45cm	black					789G		
▼ 3/8 CIRCLE MICROPPOINT REVERSE CUTTING				8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3
 G-3	13mm 13mm	45cm	black					783G		
☽ 3/8 CIRCLE MICROPPOINT SPATULA				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
 TG140-6	6.5mm 6.5mm	30cm	black		W1782					
 TG140-8	6.5mm 6.5mm	30cm 45cm	black black			7733G	1732G			

SILK Braided Silk Suture

Needle length and description		Length	Colour	Gauge						
NON-NEEDED				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
LABYRINTH'	6x45cm	black						A186H	A187H	
LABYRINTH'	12x45cm	black	A182H	A183H	A184H	A185H'				
LABYRINTH'	6x75cm	black						A306H	A307H	
LABYRINTH'	12x75cm	black		A303H'	A304H	A305H				
LIGAPAK'	335cm	black								
NON-NEEDED	13x60cm 13x60cm	black black		W211	W212	W213 SA75H	W214	W215		
NON-NEEDED	10x75cm	black		SA83H	SA84H	SA85H	SA86G	SA87G		
NON-NEEDED	2x150cm	black								SA8H
NON-NEEDED	180cm	black			W192	W193	W194G	W195		
SUTUPAK'	2x150cm	black				SA5H				
NON-NEEDED				0 3.5	1 4	2 5	3 6	4 6	5 7	7 9
NON-NEEDED	2x150mm	2x150cm	black					SA10G		