

**Surgical
Gut – Chronic**














Chromic Gut Suture













Product Characteristics	Chromic Gut Suture
Physical Material	Collagen derived from sheep mucosa and/or beef serosa, treated with chromic salts to slow absorption
Construction	Virtual monofilament
Absorbable / Non-absorbable	Absorbable
Coating (if applicable)	N/A
Colours	Natural – tan (dyed blue for Ophthalmic surgery)
Available Size Range	3 through 7/0
Tensile Strength (Wound Support)	Short to Medium Term
Tensile Strength Retention in Tissue	21-28 days
Absorption Time	Absorbed in 90 days
Frequent Uses	Used in tissues that do not need extended wound support and Ophthalmology

CHROMIC GUT









Surgical Gut Suture - Chromic

Needle length and description		Length	Colour	Gauge						
▼ 3/8 CIRCLE PRECISION POINT REVERSE CUTTING				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
 P-3	13mm	45cm	chromic	687G						
▼ 3/8 CIRCLE REVERSE CUTTING				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
 M-2	13mm	30cm	chromic		725G					
 FS-2	19mm	70cm	chromic		635H	636H				
 FS-1	24mm	70cm	chromic				637H			
▼ 1/2 CIRCLE PRECISION POINT REVERSE CUTTING				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
 PS-5	13mm	45cm	chromic		1653G					
▼ 1/2 CIRCLE REVERSE CUTTING				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
 J-1	17mm	70cm	chromic		U213H					
 CP-2	26mm	70cm	chromic				887H			
 CP	40mm	70cm	chromic						805H	
▼ STRAIGHT CUTTING				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
 KS	60mm	75cm	chromic			654H				
▼ 1/2 CIRCLE CONVENTIONAL CUTTING STERNUM				3/0 3	2/0 3.5	0 4	1 5	2 6	3 7	4 8
 CCS-1	60mm	75cm	chromic							
● 3/8 CIRCLE TAPER POINT				4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6	3 7
 BB	17mm	70cm	chromic	4887G						






CHROMIC GUT Surgical Gut Suture - Chromic

Needle length and description		Length	Colour	Gauge						
☉ 1/2 CIRCLE TAPER POINT				4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6	3 7
RB-1 	17mm	70cm	chromic	U203H	U204H					
CT-3 	22mm	70cm	chromic		892H	893H				
SH-1 	22mm	70cm	chromic	G181H	G182H					
SH 	26mm	70cm	chromic	G121H	G122H	G123H	G124H			
SH 	26mm 26mm	90cm	chromic			G173H				
CT-2 	27mm	70cm	chromic			883H	884H	885H		
MO-5 	31mm	70cm	chromic					4995H		
CT-1 	36mm	70cm 90cm	chromic chromic		810H	811H	812H	813H 925H	814H	
MH 	36mm	70cm	chromic		G126H	G127H	G128H			
MO-4 	36mm	70cm	chromic				4924H	4925H		
CT 	40mm	70cm 90cm	chromic chromic			801H	802H 914H	803H 915H		
LH 	40mm	70cm	chromic			G133H				

CHROMIC GUT Surgical Gut Suture - Chromic

Needle length and description		Length	Colour	Gauge						
☉ 1/2 CIRCLE TAPER POINT				4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6	3 7
 MO-2	40mm	70cm	chromic							
 CTX	48mm	70cm 90cm	chromic chromic			863H	904H	865H 905H		
☽ 1/2 CIRCLE TAPER POINT, 1/2 CIRCLE REVERSE CUTTING				4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6	3 7
 CT-1 / CP-1	36mm 36mm	150cm	chromic			983H*				
☉ 5/8 CIRCLE TAPER POINT				4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6	3 7
 UR-6	27mm	70cm	chromic			N878H				
☾ 1/2 CIRCLE TAPERCUT				4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6	3 7
 V-7	26mm	75cm	chromic							
 V-34	36mm	90cm	chromic			943H	944H	945H		
▲ TROCAR POINT				4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6	3 7
 V-38	31mm	75cm	chromic							
 V-37	40mm	75cm	chromic							

CHROMIC GUT Surgical Gut Suture - Chromic

Needle length and description	Length	Colour	Gauge						
▲ TROCER POINT			4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6	3 7
 V-40 48mm	75cm	chromic							
▼ 3/8 CIRCLE MICROPOINT REVERSE CUTTING			6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4	1 5
 G-1 10.5mm	45cm	chromic	795G						
 G-1 10.5mm 10.5mm	45cm	chromic	790G						
◀ 1/4 CIRCLE SABRELOC CONVENTIONAL SPATULA			7/0 0.5	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4
 S-4 8mm 8mm	45cm	chromic				1793G			
◀ 1/4 CIRCLE SABRELOC SPATULA			7/0 0.5	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4
 S-14 8mm 8mm	45cm 45cm	chromic blue		1791G*					
NON-NEEDED			4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6	3 7
NON-NEEDED	150cm	chromic	S111H	S112H	S113H	S114H	S115H	S116H	