























PLAIN GUT

Surgical Gut Suture - Plain

Needle length and description		Length	Colour	Gauge						
▼ 3/8 CIRCLE PRECISION POINT REVERSE CUTTING				7/0 0.5	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4
 P-3	13mm	45cm	undyed			686G	1644G			
 PS-2	19mm	45cm	undyed			1626G	1627H			
▼ 3/8 CIRCLE REVERSE CUTTING				7/0 0.5	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4
 FS-2	19mm	70cm	undyed			H820G	H821H	H822H		
 FS-1	24mm	70cm	undyed					H810H	H823H	
▼ 1/2 CIRCLE REVERSE CUTTING				7/0 0.5	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4
 M-1	13mm	30cm	undyed				742G			
 CP-1	36mm	70cm	undyed						898H	
▼ STRAIGHT CUTTING				7/0 0.5	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4
 SC-1	13mm	45cm	undyed				1824H			
▲ 3/8 CIRCLE CONVENTIONAL CUTTING PRIME				7/0 0.5	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4
 PC-1	13mm	45cm	undyed			1915G				
FAST ABSORBING PLAIN GUT  PC-3	16mm	38cm	undyed							
 PC-25	26mm	75cm	undyed							
◎ 1/2 CIRCLE TAPER POINT				7/0 0.5	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4
 RB-1	17mm	70cm	undyed				U207H			
 CT-3	22mm	70cm	undyed					N862H	N863H	
◎ 1/2 CIRCLE TAPER POINT				7/0 0.5	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4
 SH	26mm	70cm	undyed				G321H	Z339H†		
 CT-1	36mm	70cm	undyed							844H

PLAIN GUT Surgical Gut Suture - Plain

Needle length and description		Length	Colour	Gauge						
☉ 1/2 CIRCLE TAPER POINT				7/0 0.5	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4
 CT-1	36mm	70cm	undyed							844H
 CT	40mm	70cm	undyed							
 CTX	48mm	70cm	undyed							872H
 XLH	70mm	70cm	undyed							53T
☽ 1/2 CIRCLE TAPERCUT				7/0 0.5	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4
 V-7	26mm	45cm	undyed					849G		
▼ 3/8 CIRCLE MICROPOINT REVERSE CUTTING				7/0 0.5	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4
 G-1	10.5mm	45cm	undyed		774G					
 G-1	10.5mm 10.5mm	45cm	undyed		770G					
☞ 1/4 CIRCLE SABRELOC SPATULA				7/0 0.5	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4
 S-14	8mm 8mm	45cm	undyed			1772G				
NON-NEEDED				7/0 0.5	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4
NON-NEEDED		150cm	undyed	S102H	S103H		S105H			