

Surgical Stainless Steel












Surgical Steel Suture

Product Characteristics	Surgical Steel Suture
Physical Material	316L Stainless Steel
Construction	Braided and monofilament
Absorbable / Non-absorbable	Non-absorbable
Coating (if applicable)	N/A
Colours	Steel
Available Size Range	7 through 10/0 (No 4/0 Variant Available)
Tensile Strength (Wound Support)	Long term wound support
Tensile Strength Retention in Tissue	Indefinite
Absorption Time	Non-absorbable
Frequent Uses	Abdominal wall and skin closure; retention; tendon repair; Orthopaedic and Neurosurgery; Sternal closure in Cardiovascular surgery










STEEL

Stainless Steel Wire

Needle length and description	Length	Colour	Gauge						
○ 1/2 CIRCLE ETHIGUARD BLUNT POINT			1 4	2 5	3 5	4 6	5 7	6 8	7 9
 BP-1 65mm	2x75cm	undyed					D8846		
▼ 3/8 CIRCLE REVERSE CUTTING			4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	3 5
 FS-2 19mm	45cm	undyed	603G						
 FS-1 24mm	45cm	undyed		631G					
 FSL 30mm	45cm	undyed							
 FSLX 40mm	45cm	undyed			618G				
 LS-1 45mm	50cm	undyed						D3000	
 LS-1 45mm	4x50cm	undyed						M400G	
▼ 1/2 CIRCLE REVERSE CUTTING			4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	3 5
 CP-1 36mm	50cm	undyed		497G					
▼ 1/2 CIRCLE REVERSE CUTTING			1 4	2 5	3 5	4 6	5 7	6 8	7 9
 CPXX 55mm	75cm	undyed							







STEEL

Stainless Steel Wire

Needle length and description	Length	Colour	Gauge							
▲ 1/2 CIRCLE CONVENTIONAL CUTTING STERNUM			1 4	2 5	3 5	4 6	5 7	6 8	7 9	
CCS 	48mm	2x45cm	undyed					M657G		
CCS 	48mm	4x45cm	undyed				M652G	M653G	M654G	M655G
CCS-1 	60mm	2x75cm	undyed							
● 1/2 CIRCLE TAPER POINT			4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	3 5	
SH-2 	20mm	75cm	undyed							
CT-2 	27mm	45cm	undyed					D9090		
MH-1 	31mm	75cm	undyed		W311					
CT 	40mm	2x45cm	undyed							D9598
▼ 1/2 CIRCLE TAPERCUT			1 4	2 5	3 5	4 6	5 7	6 8	7 9	
V-37 	40mm	4x45cm	undyed	M660G						
V-40 	48mm 48mm	45cm	undyed					D8405		

STEEL

Stainless Steel Wire

Needle length and description	Length	Colour	Gauge						
1/2 CIRCLE TAPERCUT			1 4	2 5	3 5	4 6	5 7	6 8	7 9
 V-40 48mm	2x45cm	undyed					M656G		
 V-40 48mm	4x45cm	undyed				M651G	M650G	M649G	
 V-56 60mm	45cm	undyed							
 V-56 60mm	4x75cm	undyed					M435G		D9598
TROCAR POINT			1 4	2 5	3 5	4 6	5 7	6 8	7 9
 TR-55 55mm	75cm	undyed							
3/8 CIRCLE TAPERCUT, STRAIGHT CUTTING			4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	3 5
 V-26 / TS 26mm 48mm	35cm	undyed							

