

MEDISAFE PLUS PRODUCTS

VOLUME 1 2025



IMAGE	SKU#	DESCRIPTION	POWER CONSUMP.	REFRIGERANT	FLYER
	MEDISAFE PLUS 80 Dimensions (mm) 500w x 560d x 830h Weight (kg) 43	Gross Capacity 80 L +2C to +8C holding temperature 1 Door / 3 shelves Integrated AQ Box Monitoring System	1.05 kWh / 24H	R600a	DOW NLOAD
	MEDISAFE PLUS 420 Dimensions (mm) 600w x 700d x 2200h Weight (kg) 92	Gross Capacity 420 L +2C to +8C holding temperature 1 Door / 6.5 shelves Integrated AQ Box Monitoring System	2.46 kWh / 24H	R290	DOW NLOAD
	MEDISAFE PLUS 650 Dimensions (mm) 745w x 745d x 2200h Weight (kg) 175	Gross Capacity 610 L +2C to +8C holding temperature 1 Door / 5 shelves Integrated AQ Box Monitoring System	2.36 kWh / 24H	R600a	DOW NLOAD
	MEDISAFE PLUS 1000 Dimensions (mm) 1130w x 745d x 2200h Weight (kg) 230	Gross Capacity 980 L +2C to +8C holding temperature 2 Doors / 10 shelves Integrated AQ Box Monitoring System	3.79 kWh / 24H	R290	DOW NLOAD
	UPS & Battery Back Up UPS Dimensions (mm) 192w x 428d x 318h Weight (kg) 23 Battery Back Up (mm) 192w x 428d x 318h Weight (kg) 41	Approx Holding times: Medisafe Plus 420 ~ 13 hours AQ Medical 46 ~ 22 hours	Battery Type 12V7Ah Batter Voltage 72 VDC	-	DOW NLOAD



Integrated AQ Box
Monitoring System –
scan QR code to
claim first year
subscription



QCPD
Accredited



2-year warranty
on Medisafe Plus



1800 658 120



Medisafeplus.com.au



sales@avemquirks.com.au



IMAGE	SKU#	DESCRIPTION	POWER CONSUMP.	REFRIGERANT	FLYER
	AQ MEDICAL 46 Dimensions (mm) 430w x 475d x 543h Weight (kg) 22	Gross Capacity 41 L +2C to +8C holding temperature 1 Door / 1.5 shelves Optional AQ Box Monitoring System	0.12 kWh / 24H	R600a	DOW NLOAD
	AQ MEDICAL 110 Dimensions (mm) 595w x 595d x 820h Weight (kg) 41.5	Gross Capacity 145 L +2C to +8C holding temperature 1 Door / 4 shelves Optional AQ Box Monitoring System	1.26 kWh / 24H	R600a	DOW NLOAD
	AQ MEDICAL 280 Dimensions (mm) 584w x 584d x 1575h Weight (kg) 79	Gross Capacity 280 L +2C to +8C holding temperature 1 Door / 6 shelves Optional AQ Box Monitoring System	2.6 kWh / 24H	R600a	DOW NLOAD
	AQ MEDICAL 350 Dimensions (mm) 580w x 580d x 1881h Weight (kg) 89	Gross Capacity 350 L +2C to +8C holding temperature 1 Doors / 6 shelves Optional AQ Box Monitoring System	3.3 kWh / 24H	R600a	DOW NLOAD
	AQ MEDICAL 1000 Dimensions (mm) 1200w x 697d x 1998h Weight (kg) 190	Gross Capacity 960 L +2C to +8C holding temperature 2 Doors / 8 shelves Optional AQ Box Monitoring System	22.47 kWh / 24H	R290	DOW NLOAD



QCPP
Accredited



24-month
warranty

**GREENHOUSE & ENERGY
MINIMUM
STANDARDS
REGULATOR**

Registered as per the
Greenhouse and Energy
Minimum Standards Act 2012



1800 658 120



Medisafeplus.com.au



sales@avemquirks.com.au



MEDISAFE BIOMEDICAL PRODUCTS

VOLUME 1 2025



IMAGE	SKU#	DESCRIPTION	POWER CONSUMP.	REFRIGERANT	FLYER
	MEDISAFE ULTF-C74i Dimensions (mm) 550w x 550d x 921h Weight (kg) 73	Gross Capacity 74 L -20C to -86C holding temperature Max Ambient/Humidity 30C / 65%	6.56 kWh / 24H	Nature R	DOW NLOAD
	MEDISAFE BIOMEDICAL 92 Dimensions (mm) 595w x 585d x 808h Weight (kg) 61	Gross Capacity 92 L -20C to -86C holding temperature Max Ambient/Humidity 25C / 65%	6.7 kWh / 24H	Nature R	DOW NLOAD



24-month
warranty



1800 658 120



Medisafeplus.com.au



sales@avemquirks.com.au

