



Hamilton Latex Surgical Glove



- Natural Rubber Latex
- Fingertip Textured
- Curved Fingers
- Sterilised & Powder-Free
- Tested for Resistance to Chemical Permeation



A brand by
D.MUN



Hamilton Latex Surgical Glove

The GloveOn® Hamilton Surgical Glove is defined by its high level of comfortability and flexibility. Made with natural rubber latex, its feel and fit allows surgeons to perform to the best of their ability.

- ✔ Provides exceptional elasticity, dexterity and grip
- ✔ Naturally curved shape enables a better fit to the wearer's hands and reduce muscle fatigue
- ✔ Sterilised using gamma radiation which provides protection with chemotherapy drugs.
- ✔ Protection against particulate radioactive contamination (EN421:2010)
- ✔ Tested in accordance with AS/NZ 4179
- ✔ Hand Specific and comes in sizes 5.5-9.0



EN ISO 374-1:2016/Type B



EN ISO 374-5:2016



Tested in accordance with EN ISO 374-5:2016
Resistance to Bacteria and Fungi - pass
Resistance to Viruses - pass

EN421:2010



Protection against particulate radioactive contamination (Excluding Clause 4.3)



A brand by
D.MUN