


PE-1204

Bluetooth ECG




The tablet not included

• Product Features


- Simultaneous 12 leads ECG acquisition and display
- Lightweight and ultra-compact design
- Sampling box supports wired and wireless ECG acquisition
- Compatible with Windows and Android operating system
- Automatic ECG measurements and interpretation
- Data access with HIS/ EMR/ PACS based on DICOM/ FTP standard protocol
- Multi- format reports export, such as PDF, JPEG, BMP, EM-XML, SCP*, DICOM, FDA-XML*
- The comprehensive patient case management supports quick retrieve, query or combined query, delete, edit, upload and export
- Patient information can be retrieved from the worklist*

Items with * are optional

• Software Interface




12 Leads Mode



Lead Off Alarm

• Flexible and Convenient Data Transmission Solution



• Multiple Report Formats with Automatic ECG Measurements and Interpretation

