


OtoClear® Ear Irrigation Tip

Bionix OtoClear® Ear Irrigation allows healthcare professionals to utilize a soft, flared, 3 streamered tip that eliminates the mess, improves patient safety, and is compatible with a variety of delivery devices to meet the personal preference of every facility.



PROBLEM

Traditional ear irrigation shoots water directly into the ear canal, which can be dangerous and run the risk of causing perforation to the eardrum.

SOLUTION

OtoClear® Ear Irrigation Tips direct water toward the ear canal wall, with three divergent streams, instead of directly at the Tympanic Membrane. The tip is also designed with a soft, flexible material for patient comfort and flared design to prevent backslash and over insertion.

RESULTS

A clean and safe procedure that improves efficiency and patient experience.

“Our office is very happy with the soft (OtoClear) tips. They are easier to maneuver in the ear canal and more comfortable for the patient.” - Laura Griffith, MA

POTENTIAL REIMBURSEMENT (U.S. ONLY)

CPT Code/Description

69209 - Removal impacted cerumen irrigation/lavage unilateral



OTOCLEAR TIPS ARE
SINGLE USE ONLY

R_x ONLY

OTOCLEAR TIPS ARE
MADE IN THE U.S.A.



INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS

800.551.7096 • www.bionix.com

©2019 Bionix Development Corp. | US Patent # 6,706,023 | MKT-0041 | Rev. B

OtoClear® Ear Irrigation Tip

PRODUCT OFFERING



OtoClear WaterPik® Kit

Recommended for high volume facilities

- A completely portable irrigation delivery device that can be used for procedures anywhere in or out of the office.
- Provides effective, pulsed irrigation.
- Higher volume than manual devices.
- Electronic pressure control system with 3 settings.
- Proven safety & effectiveness. (Children's Hospital of Pittsburgh Study)
- Rechargeable battery.
- 210 ml self-contained water reservoir.

#7280 - 1 Portable WaterPik®, 40 OtoClear tips, 1 Ear Basin

#7245 - 1 Portable WaterPik®, 5 OtoClear tips, 1 Ear Basin



OtoClear SprayWash Kit

- Manual, intermittent lavage.
- Delivers greater volume compared to syringes or other manual devices.
- Efficient design allows for two-handed procedure.
- 500ml bottle volume.

#7290 - 1 SprayWash bottle, 20 OtoClear tips, 1 Ear Basin

#7295 - 1 SprayWash bottle, 5 OtoClear tips, 1 Ear Basin



OtoClear Tabletop WaterPik® Kit

- Provides effective, pulsed irrigation.
- Higher volume than manual devices.
- Manual pressure control with 6 settings.
- 1035ml water reservoir.

#7240 - 1 Tabletop WaterPik®, 5 OtoClear tips, 1 Ear Basin



OtoClear Aquabot® Kit

- Delivers greater volume compared to syringes or other manual devices.
- Efficient design allows for two-handed procedure.
- Provides continuous stream without manually pumping.
- 650ml bottle volume.

#7270 - 1 Aquabot® bottle, 20 OtoClear tips, 1 Ear Basin

#7275 - 1 Aquabot® bottle, 5 OtoClear tips, 1 Ear Basin



OtoClear Ear Irrigation Tips

- Directs water to the ear canal wall, with three divergent streams, instead of directly at the Tympanic Membrane.
- Flared design prevents backslash and over insertion.

#7200 - OtoClear Ear Irrigation Tips (40/box)

#7202 - OtoClear Ear Irrigation Tips (20/box)

OtoClear Tip US Patent # 6,706,023.

Waterpik® is a registered trademark of Water Pik, Inc. The OtoClear® Ear Irrigation Tip is not affiliated with Water Pik Technologies, Inc.

Aquabot® is registered under Lunatech®, used under permission of Lunatech®, Inc. All rights reserved.



INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS

800.551.7096 • www.bionix.com