

Wide Bore Suction Tubing



Fairmont Medical offers an extensive range of Wide Bore Suction Tubing solutions.

As well as offering Australian made standard suction tubing, we also offer speciality tubing including ENT suction tubing, soft tubing, wide bore tubing, irrigation tubing, Frazier tubing and Caws tubing. All tubing is available in multiple lengths, and can be tailored to suit specific customer requirements.

- Phthalate free PVC complies to European standards
- Range of connectors available for different applications
- Ribbed to allow easy release from packaging
- Single patient use




Wide Bore Liposuction Tubing

Wide Bore Suction Tubing


Wide Bore Tubing - Sterile, Disposable

- 2 White Hexagonal Connectors

Product	Product Code	Length	Tube ID	Tube OD	Box Qty
	DST3402	2m	6.9mm	10.5mm	60
	DST3403	3m	6.9mm	10.5mm	60
	DST3445	4.5m	6.9mm	10.5mm	40
	DST3406	6m	6.9mm	10.5mm	30
	DST3408	8m	6.9mm	10.5mm	25

Wide Bore Liposuction Tubing - Sterile, Disposable

- Blue Connectors

Product	Product Code	Length	Tube ID	Tube OD	Box Qty
	DST3503	3m	9.8mm	15mm	10
	DST3563 (No Connectors)	3m	9.8mm	15mm	10