



Non-Irritating



Absorbent
Inner



Protects
Eye

Sentry Medical Eye Pads are constructed with an absorbent inner and covered with non-woven outer layer with sealed edges.

SHAPED TO FIT

Oval shaped to specifically fit the eye socket

SEALED EDGES

Non-woven outer layer with sealed edges prevents linting

SOFT AND COMFORTABLE

Non abrasive outer layer soft on delicate skin



Suitable for:

Minor eye injuries

Post surgical dressing

First aid

Reorder	Size	Sterility	Pcs/Unit	Units/Ctn
EP001	6cm x 8cm	Sterile	50	24
EP002	6cm x 8cm	Non Sterile	2000	1

