



Low Linting



Highly Absorbent



X-Ray Detectable

ThreadX Non-Woven Balls are constructed from a polyester and rayon blend, with a black x-ray detectable thread heat bonded to the material.

X-RAY DETECTABLE

X-ray detectable thread for patient safety

HIGHLY ABSORBENT

Suitable for controlling saliva and bleeding



Suitable for:

Open surgical procedures

Absorption of exudate

Reorder	Pcs	Sterility	Pcs/Unit	Units/Ctn
XNB020	5pcs	Sterile	600	1

