

# MICRO-TOUCH<sup>®</sup> AFFINITY<sup>™</sup>

Performance Series

## Neoprene, powder-free examination glove with excellent elasticity and barrier protection

This non-sterile neoprene synthetic formulation provides superior protection, comfort, fit, elasticity and tactile sensitivity with excellent resistance to a broad range of chemicals.

### KEY FEATURES AND BENEFITS:

- Offers comfort level similar to latex.
- Textured finger surfaces provide a secure grip in both dry and wet conditions for enhanced instrument handling.
- Excellent tear and puncture resistance
- Protects from latex Type I allergy
- Tested for use with chemotherapy drugs\*

\* See product packaging or contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

**WARNING:** No glove completely protects against all chemicals. Users should test the suitability of this product against the specific chemicals and environment where used.

### Recommended For

- Tested for use with chemotherapy drugs
- Protection from Type I latex allergy in HCW's or patients
- Enhanced grip in wet environments
- Sustained, high-risk examination procedures

### FEATURED TECHNOLOGIES



## MICRO-TOUCH<sup>®</sup> Everyday, Everywhere Protection

MICRO-TOUCH<sup>®</sup> premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH<sup>®</sup> AFFINITY<sup>™</sup>

Performance Series

## PRODUCT INFORMATION

|   |  |  |   |
|---|--|--|---|
| Material                                | Neoprene (Polychloroprene)   | Not Made From Natural Rubber Latex     | Yes   |
| Color                                   | Green  | Cuff Length                            | Standard  |
| Powder Content                          | Powder-Free  | External Glove Surface                 | Textured Fingers                                      |
| Freedom from Holes (Inspection level I) | 1.5 AQL  | Palm Thickness (mm/mil)                | 0.13 / 5.1  |
| Finger Thickness (mm/mil)               | 0.16 / 6.3   | Allergy Prevention                     | Latex (Type I)  |
| Available Sizes                         | XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10) | Tested For Use With Chemotherapy Drugs | Yes, in accordance with ASTM D6978 and US FDA cleared |
| Sterile                                 | No   | Antistatic                             | Not Tested  |
| Glove Length (mm/inches)                | 240 / 9.5  | Product Segmentation                   | Performance Series                                    |

## ORDERING INFORMATION

| Size         | XS (5.5 - 6) | S (6.5 - 7) | M (7.5 - 8) | L (8.5 - 9) | XL (9.5 - 10) |
|--------------|--------------|-------------|-------------|-------------|---------------|
| Product Code | 3770         | 3771        | 3772        | 3773        | 3774          |

## PACKAGING AND STORAGE

|                      |  |
|----------------------|--|
| Packaging            | 100 gloves per dispenser; 10 dispensers/1000 gloves per case   |
| Shelf Life           | 5 Years  |
| Storage Instructions | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition. |

## STANDARDS AND CERTIFICATIONS

ASTM D6977, CA Prop 65, Category III, EN 455 Part 1, EN 455 Part 3, ISO 13485, ISO 10993-10:2002

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email:  
insidesalesus@ansell.com

Latin America & Caribbean  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email:  
luvas.medicas@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email:  
apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email:  
protection@ap.ansell.com