

# If there's something to discover, you'll be sure to see it. The new HEINE DELTA 30 Dermatoscope.

Not only would our new high-end dermatoscope win any beauty pageant; it stands out due to the best imaging we've ever developed. Combined with a field of view measuring a veritable 30 millimetres and unprecedented colour rendering, in LED<sup>HQ</sup>, of course. Polarisation provides glare-free and

non-reflective working conditions. Above all else, it is comfortable to operate with its angled, ergonomic design and intuitive control. Feel free to take a look – you're in for a real discovery.

Crystal-clear image with the new achromatic HEINE optical system

Largest field of view of real 30mm (lens 32mm)

**LED<sup>HQ</sup>**  
LED NOW IN HEINE QUALITY.

Excellent colour rendering due to LED<sup>HQ</sup>

Ergonomic handling thanks to the angled product design

Digital documentation with the exclusive HEINE DERM app or the HEINE Cube System\*



Table charger for DELTA 30



DELTA 30 with contact plate with scale, USB cord with medical approved plug-in power supply, case

**K-230.28.305**

Table charger for DELTA 30

**X-002.99.212**

∴ Accessories and further articles under [heine.com/en/DELTA30](http://heine.com/en/DELTA30)

\*Not available in all countries.

