


# Active Protection

Effective exudate management and effective antimicrobial action in one dressing

*An easy-to-use, convenient, highly cost-effective solution*

ALLEVYN<sup>®</sup> Ag is a unique mix combining the existing ALLEVYN<sup>®</sup> triple layer foam technology with enhanced silver sulfadiazine (SSD) providing long wear time with long-lasting antimicrobial protection.


## Effectiveness against MRSA over 24 hours<sup>2</sup>



2. Data on file, ref 0807055.

## Enhanced SSD is released into the wound bed over seven days:

- Effective against a broad range of pathogens, including bacteria, yeast, fungi and antibiotic-resistant strains
- Easy to use – silver is released when in contact with wound fluid
- Does not form pseudoeschar (unlike other SSD preparations)



ALLEVYN<sup>®</sup> Ag Adhesive



ALLEVYN<sup>®</sup> Ag Non-Adhesive

**ALLEVYN<sup>®</sup> Ag is indicated for all exuding wounds where the risk of infection is suspected or needs to be managed, including pressure ulcers, leg ulcers, diabetic foot ulcers, skin tears, post operative wounds and donor sites.,**

### ALLEVYN<sup>®</sup> Ag Adhesive

Code	Description	Qty	ARTG
66800073	7.5cm x 7.5cm	box/10	155735
66800075	10cm x 10cm	box/10	155735
66800078	12.5cm x 12.5cm	box/10	155735
66800081	17.5cm x 17.5cm	box/10	155735

### ALLEVYN<sup>®</sup> Ag Non-Adhesive

Code	Description	Qty	ARTG
66800083	5cm x 5cm	box/10	155652
66800086	10cm x 10cm	box/10	155652
66800089	15cm x 15cm	box/10	155652
66800091	20cm x 20cm	box/10	155652

Smith & Nephew Pty Ltd Healthcare Division


**Australia**  
 T 61 3 8540 6777 F 61 3 9544 5086  
[www.smith-nephew.com.au/healthcare](http://www.smith-nephew.com.au/healthcare)  
**Customer Service**  
 T 13 13 60 F 1800 671 000

**New Zealand**  
 T 64 9 828 4059 F 64 9 820 2867  
[www.smith-nephew.com/nz](http://www.smith-nephew.com/nz)  
**Customer Service**  
 T 0800 807 663 F 0800 263 222

<sup>®</sup>Trademark of Smith & Nephew <sup>®</sup>Registered Trademark  
<sup>©</sup>Copyright of Smith & Nephew <sup>™</sup>Trademark  
 SN8091A 02/10

# All-round versatility

The only foam dressing combining a bordered silicone adhesive with effective antimicrobial action



Reference:  
 1. Smith & Nephew Data on file report 0903005.

## ALLEVYN<sup>®</sup> Ag Gentle Border kindly protects sensitive patients

- The silicone gel adhesive reduces trauma on removal for patients with fragile skin
- The enhanced silver sulfadiazine (SSD) provides sustained antimicrobial protection
- The ALLEVYN<sup>®</sup> foam prevents leakage and reduces risk of maceration of the wound and surrounding skin – fewer painful dressing changes means less stress for your patients

ALLEVYN<sup>®</sup> Ag Gentle provides a soft gel alternative allowing easy application of secondary retention.

ALLEVYN<sup>®</sup> Ag is a unique mix combining the existing ALLEVYN<sup>®</sup> triple layer foam technology with enhanced silver sulfadiazine (SSD) providing long wear time with long-lasting antimicrobial protection.

## Enhanced SSD is released into the wound bed over seven days:

- Effective in vitro against a broad range of pathogens, including bacteria, yeast, fungi and antibiotic-resistant strains
- Easy to use – silver is released when in contact with wound fluid
- Does not form pseudoeschar (unlike other SSD preparations)



ALLEVYN<sup>®</sup> Ag Gentle



ALLEVYN<sup>®</sup> Ag Gentle Border

**ALLEVYN Ag Gentle and Ag Gentle Border are indicated for wounds on patients with fragile or sensitive skin where the risk of infection is suspected or needs to be managed. Wound types include skin tears, partial and full thickness burns, donor sites, pressure ulcers, leg ulcers and diabetic foot ulcers.**

### ALLEVYN<sup>®</sup> Ag Gentle

Code	Description	Qty	ARTG
66800464	5cm x 5cm	box/10	169183
66800465	10cm x 10cm	box/10	169183
66800466	10cm x 20cm	box/10	169183
66800467	15cm x 15cm	box/10	169183
66800468	20cm x 20cm	box/10	169183

### ALLEVYN<sup>®</sup> Ag Gentle Border

Code	Description	Qty	ARTG
66800460	7.5cm x 7.5cm	box/10	169184
66800461	10cm x 10cm	box/10	169184
66800462	12.5cm x 12.5cm	box/10	169184
66800463	17.5cm x 17.5cm	box/10	169184