



NINGBO YINGMED MEDICAL INSTRUMENTS CO., LTD.

Add: Yingsheng Building No.456 Tai'an Road, Southern Business Zone
Yinzhou District, Ningbo, China

TEL :0086 574 88226978 FAX:0086 574 89065812

Technical Data sheet

Item name: KN95 face mask

Material	Gram	50g/m2	25g/m2	25g/m2	50g/m2	40g/m2
	Material	SPP	Melt-blow	Melt-blow	Soft Cotton	SPP
	Color	White	white	white	white	white
	Size	15.5*10.5cm				
	Type	With ear loop				
	Ply	5ply				
Technique	Technology	Semi-Automatic Ultrasonic machine				
	Cutting	Neat				
	connect	Spot welding				
Cleanliness	Overall	Clean				
Packing	Carton	Corrugate case				
	Inner package	White box				
	Packing size	57*40*32cm				
	Quantity per ctn	10pcs/box,40box/ctn				
Result	All test results comply with the requirements of Sandard GB2626-2006					