



# Warts.

Safely Treated.  
Kid's Hands Approved.

## Histofreezer® Portable Cryosurgical System

- **Safe**  
Low risk of scarring, safely used on children 4+ years old<sup>1,2</sup>
- **Clinically proven**  
93.4% resolution rates compared to liquid nitrogen<sup>3</sup>
- **Easy-to-use**  
Portable, 1-minute, one-handed procedure<sup>4</sup>




## Proven Performance\*



\*Occasionally multiple treatments spaced 2 weeks apart may be required. Verruca vulgaris shown.

## Easy 3-Step Procedure

1. Charge 2. Invert 3. Apply



# Histofreezer®

Tough on warts. Gentle on skin.

References 1) Shepherd, J.P and Dawber, R.P.R. Wound Healing and Scarring After Cryosurgery. *Cryobiology* 1984;Vol. 21:157-169. 2) Ricketti, J.C. and Niedbala, R.S. A Study of a New Method of Cryosurgical Treatment of Verrucae Plantaris. *The Lower Extremity* 1995;Vol. 2, No. 3. 3) Martinez, F.C. et al. Dermatological Cryosurgery in Primary Care with Dimethyl Ether Propane Spray in Comparison with Liquid Nitrogen. *Atención Primaria* 1996; Vol. 18, No. 5:211-216. 4) Histofreezer Portable Cryosurgical System Package Insert. ©2017 OraSure Technologies, Inc. Histofreezer is a registered trademark of OraSure Technologies, Inc. HISTO0376 02/17



220 East First Street  
Bethlehem, PA 18015 USA  
www.OraSure.com  
www.Histofreezer.com



## Histofreezer®

### Advanced Measurement Technology

The Histofreezer® Portable Cryosurgical System is a low cost, highly effective cryosurgical agent. This therapeutic treatment makes use of local freezing for the controlled destruction or removal of living tissue.

#### Key Points:

- Ultra-portable, one-handed operation, treatment takes less than 60 seconds to perform
- Results comparable to those achieved with liquid nitrogen
- Requires no bulky equipment or large investment, unlike Liquid Nitrogen Systems
- High patient acceptance, relatively painless and bloodless with minimal threat of scarring
- Highly efficient cryogen delivery (low waste) and no messy or leaky cones
- 2 x 80mL complete Histofreezer (canister + applicators)

## Histofreezer®



|           |                 |                 |
|-----------|-----------------|-----------------|
|           | 2x 80mL         | 5+              |
| HM1008301 | <b>\$199.00</b> | <b>\$169.00</b> |

## CryoPen B Kit

### Ideal for General practice and veterinary

- Now includes 3 Micro applicator tips: gold dot 1-4mm, silver dot 2-8mm & bronze dot 4-12mm micro-applicators.
- To be used with 23.5g cartridges



# \$2,200

CPBPLUS

|           |  |                 |
|-----------|--|-----------------|
| CPBP-C235 | Cartridge for Cryopen BPlus 6 x 23.5g Suitable for Cryopen-B model | <b>\$299.00</b> |
|-----------|--|-----------------|

## CryoPen O Kit

### Excellence in cryotherapy made affordable

- Now includes 3 Micro applicator tips: Blue dot 1-3 mm, white dot 2-5 mm & green dot 4-8 mm micro-applicators.
- To be used with 8g & 16g cartridges.



# \$2,395

CRYOPEN-O

|       |                                 |                 |
|-------|---------------------------------|-----------------|
| CPC8  | Cartridge for Cryopen O 24 x 8g | <b>\$195.00</b> |
| CPC16 | Cartridge for Cryopen O 6 x 16g | <b>\$95.00</b>  |

## Silver Nitrate Applicators

- 75% Silver Nitrate
- Single Use 15cm Bendable Applicator



|         |         |         |    |                |
|---------|---------|---------|----|----------------|
| SNA7484 | Pkt/100 | \$75.00 | 6+ | <b>\$69.00</b> |
|---------|---------|---------|----|----------------|



MEDICAL PRODUCTS  
**WALMAY**

# CRYOSPRAY

MODEL - WM161

The Walmay Cryospray has given over two decades of well proven reliable service to Medical Professionals.

It is the instrument of choice for all Dermatologists and Practitioners.



### CONTENTS

- 0.3L Flask capacity.
- Set of four cone-shaped spray tips (with varying size orifice) for accurate visual aiming during treatment.
- One curved spray tip for body orifice surface treatment.



- Includes USB containing a video for setting up and using the Walmay Cryospray WM161.

**ONLY**  
**\$875+GST** WMWAL

|         |                                 |                |
|---------|---------------------------------|----------------|
| WMCP121 | Extension Probe for Model WH121 | <b>\$65.00</b> |
| WMWP    | Extension Probe for Model WH099 | <b>\$65.00</b> |

NEW & IMPROVED MODEL

FIVE YEAR WARRANTY

## CryoPro Mini Flask 350mL

Allows Primary Care Doctors to quickly and easily treat a range of dermal problems including Keratoses, Condaylomata, warts, skin tags and other skin lesions.

Danish design and manufactured, the CryoPro has provided over 2 decades of proven reliable service to Dermatologists and Medical Practitioners.

- Suitable for use by both GP and Specialists
- 350ml Capacity
- Comes with Applicator Set (5 Tips and 1 Cleaning Adaptor)
- 3 Year Guarantee
- Simple & Reliable
- 12 - 14 Hour Static Holding Time



|        |             |          |                 |
|--------|-------------|----------|-----------------|
|        |             | Each     | 3+              |
| ELCPRO | 350ml Flask | \$849.00 | <b>\$819.00</b> |

## BRYMILL'S CRY-AC® Hand-Held Liquid Nitrogen Delivery System

### Features:


- Stainless steel and solid brass constructions for long life
- Slender, fingertip trigger for maximum control and ease of viewing
- Wide opening neck for easy filling
- Streamlined relief valve maintains a constant operating pressure to ensure consistent and accurate freezing
- Patented Safety Autovent maximises safety by allowing internal pressure to gradually vent as the cap is unscrewed
- Sturdy base for extra stability
- 3 year warranty

### Spray kit includes:

- Flask 300ml or 500ml
- 5 spray apertures with varying orifice size
- 1 x 78mm curved spray extension
- Instruction and technique booklet



|       |             |                   |
|-------|-------------|-------------------|
| DSBRY | 300ml Flask | <b>\$1,295.00</b> |
| DSBRZ | 500ml Flask |                   |





## Liquid Nitrogen Storage Container

### Features:

- Comes with 5 year vacuum warranty, 2 year container warranty
- Made of high-strength aluminum alloy.
- Adopting advanced powder coating technique
- Each unit comes complete with safety lock lid, neck plug, leather (PU) bag and ladle.

Freight charges apply



| Ladle outer diameter                    |   | Unit | 10L      | 20L       |
|---|---|------|----------|-----------|
| Ladle height                            |   | L    | 10       | 20        |
| Number of ladles that can fit per dewar |   | kg   | 6.0      | 9.5       |
| Comes with one ladle per dewar          |   | mm   | 50       | 50        |
| Ladle outer diameter                    |   | mm   | 300      | 394       |
| Ladle height                            |   | mm   | 552      | 672       |
| Number of ladles that can fit per dewar |   | L/d  | 0.12     | 0.12      |
| Comes with one ladle per dewar          |   | days | 86       | 168       |
| Ladle                                   | Ladle outer diameter                    | mm   | 38       | 38        |
|   | Ladle height                            | mm   | 120      | 120/276   |
|   | Number of ladles that can fit per dewar | pc   | 3~6      | 3~6       |
|   | Comes with one ladle per dewar          | pc   | 1        | 1         |
|   | 0.5ml                                   | pc   | 792      | 792       |
|   | 0.25ml                                  | pc   | 1788     | 1788      |
|   | 0.5ml                                   | pc   | -        | 1284      |
|   | 0.25ml                                  | pc   | -        | 2832      |
| Code                                    |   |      | TW10LD   | TW20XTO   |
| Price                                   |   |      | \$995.00 | \$1195.00 |

## Roller Base

This sturdy, cast-aluminium, five-wheel roller base enables you to move dewars easily and safely. To fit 396mm diameter.

# \$475

FCRB400



## Ladles

Ladle 20 Litre with handle

# \$159

TWDIP3

Ladle 10 Litre with handle

# \$159

BIOLD4



## Dewar Measuring Rod

Low conducting measuring rods that clearly show the LN2 level of a dewar in inches and centimetres. (92cm - 36inch)

# \$19

FCMR033



## Universal Extractor

Adjustable to suit dewars from 10L - 25L.

# \$195

WMDEXT

