




# DERM A942™ High Frequency Desiccator

The world's most user-friendly diathermy



The DERM 942 and its patent pending handpiece are designed and engineered to be the most reliable and durable desiccator available today. The unit's life is prolonged by incorporating the latest in chip technology and reducing the number of internal parts, thereby generating less internal heat.

- 40 watt high frequency desiccator
- Adjustable in 1/10th watt increments in the first 10 watts, in monopolar and bipolar
- More accurate power with Universal Power Control 100-240 VAC
- Less to do and less to buy makes this the most cost effective and user friendly unit on the market
- The Bovie DERM 942 is a ground reference unit which means no grounding pad is required.



**\$2,095**

DSA942



**Bovie Stand**

Table-top stainless steel stand

**\$89**

DCA813

## Bipolar Accessories

Bipolar power is intrinsically safer than monopolar HF power

DC701BA	Bipolar Adson Forcep 1mm wide serrated jaw	<b>\$369</b>
DSA827EU	Bipolar Cable	<b>\$215</b>
DSA803	Footswitch	<b>\$345</b>



## Aaron Disposable Non-Sterile Electrodes

(Fits Hyfrecator)		Box/100
DSA804	Sharp	<b>\$115</b>
DSA806	Blunt	<b>\$115</b>



## Aaron Disposable Non-Sterile Pencil Sheaths

Box/100

**\$62**

DCDA910



## Aaron Disposable Cautery

- Complete battery cautery handle with fine point tip
- Individually sterile packaged

DCAA01	Box/10	<b>\$265</b>
--------	--------	--------------



High temperature cauteries are popular for many types of procedures including pinpoint hemostasis



## The Hyfrecator® 2000 Electrosurgical System



### Features:

- Easy to use and affordable
- 3 programmable power settings
- Remote control hand piece
- The medical industry's leading device for over 60 years
- Cost effective alternative to various laser procedures
- 12 months warranty

### There is only one Hyfrecator®

#### Special features include:

Subtle control to 1/10 of a watt in low mode. Dual microprocessors: one to control the function and output of the Hyfrecator®; the second to terminate activation if a hazard or safety issue is detected.


#### Hyfrecator® 2000 Contents:

- Pack of 5 Electrolase Sharp Tip disposable electrodes, non sterile
- Pack of 5 Electrolase Blunt Tip disposable electrodes, non sterile
- 1 x Hyfrecator® Electrosurgical Handpiece
- 2 x Hyfrecator® Pencil Sheath, non sterile
- 2 x Handpiece Sheath, sterile
- Inservice Training DVD
- Operator's Manual
- Power Cord




# \$2,195

MQ25-5011



## Hyfrecator® 2000 Accessories

Below are listed most common lines. Full range available.



MQ25-5070	Foot Control for Hyfrecator® 2000	7-900-6	Each	<b>\$275</b>
MQ25-5072	Pencil Up/Down 2000	7-900-5	Each	<b>\$545</b>
MQ25-5076	Foot Switch	7-900-4	Each	<b>\$325</b>
MQ25-5200	General Purpose Electrode	711B	Each	<b>\$69</b>
MQ25-5250	Electrodes 714	714	Box/6	<b>\$299</b>
MQ25-5313	Bipolar Forceps 10cm	7-809-1	Each	<b>\$725</b>
MQ25-5316	Bipolar Forceps 13cm	7-809-4	Each	<b>\$695</b>
MS255317	HYF Forcep Adson Serrated Tip	7-809-7	Each	<b>\$450</b>
MQ25-5320	Bipolar Forecp Lead	7-809-12	Set	<b>\$295</b>

## Electrolase Disposable Electrodes



**TOP  
SELLER**



Non-Sterile			Box/100
MQ25-5350	Sharp	7-100-12	<b>\$110</b>
MQ25-5351	Blunt	7-101-12	<b>\$110</b>
Sterile			Box/50
MQ25-5352	Sharp	7-100-8	<b>\$110</b>
MQ25-5353	Blunt	7-101-8	<b>\$110</b>

## Disposable Pencil Sheaths

MQ25-5105	Pencil Sheaths Non-Sterile	Box/100	<b>\$55</b>
MQ25-5106	Pencil Sheath Sterile	Box/25	<b>\$189</b>

