

# Whiter, Brighter AND Greener



## Welch Allyn Green Series™ LED Medical Exam Lights

- Lower cost of ownership due to elimination of bulb replacement and lower energy consumption
- Helps you see true colour tissue rendition, and provides a uniform light spot
- Reduced risk of cross-contamination for your patients with touchless on/off and optional disposable sheaths
- 5-year limited warranty†



### GS 300

**44456** GS 300 General Exam Light with Mobile Stand

**44416** GS 300 General Exam Light with Table/Wall Mount

### GS 600

**44606** GS 600 Minor Procedure Light with Mobile Stand

**44616** GS 600 Minor Procedure Light with Table/Wall Mount

### GS 900

**44906** GS 900 Procedure Light with Mobile Stand

**44900-C** GS 900 Procedure Light with Ceiling Mount

**44900-W** GS 900 Procedure Light with Wall Mount

### EXAM LIGHT

**48816** GS Exam Light IV with Table/Wall Mount






## Green Series™ 300 General Exam Light

Utilising a single white LED, the GS 300 Exam Light provides a whiter, brighter light than halogen and consumes less energy.

### Features

- Perfect combination of intense light output and cool operation in a compact design
- Affordable price and rugged performance
- Touchless On/Off
- Multiple mounting options for floor space maximisation
- Lower cost of ownership and energy consumption due to elimination of bulb replacement
- Easy assembly
- 5 year limited warranty to include the LED light source and power supply only



		Each	3+
WL44416	GS 300 Exam Light with wall/table mount	\$599.00	<b>\$499.00</b>
WL44456	GS 300 Exam Light with mobile stand	\$699.00	<b>\$649.00</b>




## Green Series™ 600 Minor Procedure Light

Featuring three LEDs and a wider head, the GS 600 is designed specifically for demanding viewing conditions.

### Features

- Three white LEDs for greater intensity
- Touchless power controls and adjustable light output
- Disposable sheath available to help reduce the risk of cross-contamination
- Rugged, ergonomic design allows for improved performance and long-lasting durability
- Low cost of ownership due to elimination of bulb replacement and lower energy consumption
- Easy to assemble and clean
- 5 year limited warranty to include the LED light source and power supply only



		Each	3+
WL44616	GS 600 Procedure Light with wall/table mount	\$1,099.00	<b>\$899.00</b>
WL44606	GS 600 Procedure Light with mobile stand	\$1,199.00	<b>\$999.00</b>



## Green Series™ 900 Procedure Light

This procedure light, featuring six LEDs, combines outstanding durability, extraordinary maneuverability, unparalleled spot quality and multiple mounting options.

### Features

- Six white LEDs for greater intensity
- Touchless On/Off with dimming option
- Multiple mounting options for floor space maximisation
- Disposable sheath available to help reduce the risk of cross-contamination
- Rugged, ergonomic design allows for improved performance and long-lasting durability
- Lower cost of ownership due to elimination of bulb replacement and lower energy consumption
- 5 year limited warranty to include the LED light source and power supply only



		Each	3+
WL44900-C	GS 900 Procedure Light with ceiling mount	\$2,899.00	<b>\$2,499.00</b>
WL44900-W	GS 900 Procedure Light with wall mount	\$2,899.00	<b>\$2,499.00</b>
WL44906	GS 900 Procedure Light with mobile stand	\$2,899.00	<b>\$2,499.00</b>



## Green Series™ Exam Light IV

Perfect for the OB/GYN office or for specialties that require very focused light in specific areas. The GS Exam Light IV provides you with 3x the typical light output and features a highly flexible light pipe that is twice as long as conventional pipes. With a compact, unobtrusive head, it ensures pinpoint accuracy and provides true tissue colour rendition during an exam.

### Features

- One LED - 3x the typical light output; pinpoint accuracy; true tissue colour rendition
- Highly flexible light pipe for easy maneuverability
- Lightweight pipe is twice as long as conventional pipes
- Touchless On/Off
- Multiple mounting options for floor space maximisation
- Disposable sheath available to help reduce the risk of cross-contamination
- Lower cost of ownership due to elimination of bulb replacement and lower energy consumption
- Easy to assemble
- 5 year limited warranty to include the LED light source and power supply only



	Each	3+
WL48816	\$1,099.00	<b>\$899.00</b>



**HEINE** QUALITY  
MADE IN GERMANY

# HEINE EL10 Examination Light

A real highlight in all respects with  
LED in Heine Quality Illumination

Experience The Difference

LED NOW

LED NOW IN HEINE QUALITY.

- Continuous adjustment
- One-hand operation
- Compact design





## HEINE EL10 LED® Examination Light

Only a few lights are as bright, none are as homogeneous as this. Illumination without compromises for every application. Stationary and yet flexible. Position, brightness, and spot size can be adjusted with one hand.

LED –economical with high colour fidelity. The new examination light combines many characteristics that facilitate everyday working life.

- 5 year guarantee
- Bright and homogenous light
- Adjustable spot size
- Brightness control
- Easy to position
- One handed operation
- Compact design
- Hermetically sealed housing
- Various mounting options available with Stand, Wall and Clamp mounts

		Each	3+
HEJ0827001	EL10 LED Examination light with Wall Mount	\$2,215	<b>\$2,015</b>
HEJ0827002	EL10 LED Examination light with Clamp Mount for Table Top Mounting	\$2,365	<b>\$2,165</b>
HEJ0827003	EL10 LED Examination light with Wheeled Stand & Metal Base	\$2,575	<b>\$2,375</b>



LED NOW IN HEINE QUALITY.



## Ri-magic® Examination LED Light

The flexible examination light ri-magic® LED illuminates the examination and operation area in the best possible way.

- The most modern LED technology, similar to daylight and high-performance; ideal for examinations and smaller procedures
- Lifetime of the lamp up to 50,000 hours, replacements no longer necessary
- Significantly less heating and energy loss
- The focusing device at the front of the light head allows exact adjustment of the size of the illuminated examination area
- Highly flexible, long special arm can be bent in all directions, the light remains fixed in the chosen position
- Simple to use handle with rheostat for infinitely variable regulation of light intensity
- Wide-range power adaptor for 230V/120V



RIE1010235	Wall mount	<b>\$945.00</b>
RIE1010236	Mobile mount	<b>\$1,099.00</b>
RIE11343	Mobile stand and adaptor	<b>\$385.00</b>





## HEINE EL3 LED Examination Light

### Simply bright LED Illumination

The EL 3 LED provides bright (> 42,000 lux / > 350 lumen), shadow-free illumination precisely at the point of interest. at a working distance of 30 cm.

- Shadow-Free Illumination. The compact illumination head (Ø approx. 60 mm) allows convenient, precise positioning, particularly in difficult application situations. Optimal positioning eliminates unwanted shadows in the illuminated field.
- Maintenance-free and easy to clean.
- The light weight head and stable goose neck allow easy and intuitive adjustment and keep the light exactly in the required position.
- LED in HQ. The new standard in LED illumination defining optimal light intensity, homogeneity and colour rendering for each individual examination situation – Colour temperature: 5,000 K, Colour Rendering Index (CRI) > 90 on a scale of 100. This ensures true tissue colour for the most accurate diagnosis.
- 134cm Goose Neck allowing for head to toe examination
- Hermetically Sealed Housing
- No dust traps
- Easy to clean



		Each	3+
HEJ0827011	EL 3 LED Examination Light with wall mount	\$1,029.00	<b>\$949.00</b>
HEJ0827013	EL 3 LED Examination Light with clamp for table-top mounting	\$1,175.00	<b>\$1,075.00</b>
HEJ0827014	EL 3 LED Examination Light with wheeled stand, metal-base	\$1,385.00	<b>\$1,285.00</b>



## Halogen Examination Light

Cost-effective solution for surgery examination lighting

- 35 watt globe, 12V (approx 4000hrs)
- 70 cm long gooseneck
- Easy grip handle with on/off switch
- Additional glass filter reduces heat

	Each	3+
MGL308	\$279.00	<b>\$269.00</b>



Brackets to be purchased separately



## LED Examination Light

- 700mm Gooseneck
- Sure grip handle with on/off switch
- Cool white light, 5500°K colour,
- Concise beam angle of 30°
- 3 x 1 watt high power LED giving 16,000 lux at 30cm
- Low energy consumption, only 5 watt
- Low heat output
- Wall or mobile mounted
- 40,000 hour bulb life
- 2 year warranty

	Each	3+
MGL308LED	\$489.00	<b>\$479.00</b>





## Ordisi Iglux Examination Lights

Ultimate day procedure light

- Compact design with smooth edges that make it easy to clean
- 12 month warranty
- Light intensity adjustment through the touch panel controls at front of lamp
- Handle is removable and safe to sterilise
- 20,000 - 65,000 lux adjustable brightness

		Each	3+
ODIGLUXM	Light Ordisi Iglux Mobile		
ODIGLUXC	Light Ordisi Iglux Ceiling Standard 80cm Pole	\$3490.00	<b>\$3,140.00</b>
ODIGLUXCSH	Light Ordisi Iglux Ceiling Short 34.5cm Pole		



## Ordisi Mobile Examination Light

Cool light source

- Expected LED lifetime 50,000 hours
- Very power power consumption (10w)
- Strong, white light
- 42,000 lux/0.5m
- 12,000 lux/1m
- lightweight mobile base
- 12 month warranty

	Each	2+
ODLEDMEXAM	\$585.00	<b>\$569.00</b>



## Ordisi Wall Examination Light

Includes Swivel Bracket

- Cool light source
- Expected LED lifetime 50,000 hours
- Very power power consumption (10w)
- Strong, white light
- 42,000 lux/0.5m
- 12,000 lux/1m
- Goose style neck for easy maneuverability
- 12 month warranty

	Each	2+
ODFLH2731W	\$585.00	<b>\$569.00</b>





## Green Series Mobile Stand Examination Light IV

Examination Light IV

**\$349**

WL48950

Examination Light Tall

**\$349**

WL48960



Freight charges will apply

**WelchAlllyn**

## Mobile Base for LHH Examination Light

Freight charges will apply



**\$345**

LLTS

**LUXO**

## Mobile Base

- Stable 5-star solid base
- Black 50mm twin wheel castors
- Chrome plated stem with key wheel
- Height 850mm stem
- Diameter 12.5mm

Freight charges will apply



**\$149**

VIDCD

## Extended Wall Mount for GS300/600 Examination Light



**\$255**

WL44215

**WelchAlllyn**

## Ceiling Mount Bracket for LHH 60cm



**\$355**

LU2529

**LUXO**

## Wall Bracket and Base



**\$75**

LUX899

**LUXO**

## Maggylamp Wall Bracket



**\$65**

MGL307

**MAGGYLAMP**  
medical

## Maggylamp Mobile Base Bracket



**\$52**

MGL309

**MAGGYLAMP**  
medical

## ri-magic Mobile Stand + Adaptor

**Riester**

**\$385**

RIE11343



## Disposable Sheaths for Green Series™ 600/900

**\$205**

WL52630  
Box 250



**WelchAlllyn**

## Disposable Sheaths for Green Series™ Exam Light IV

**\$39**

WL52640  
Box 25



**WelchAlllyn**

## 35 Watt Halogen Bulb for Maggylamp

**\$24**

MGL306



**MAGGYLAMP**  
medical



## MAGGY-VUE Magnifying Lamp

- 4 Diopter circular lens, Ø 120mm
- 2 x magnification
- Shadow-free diffused lighting
- 22 Watt fluorescent lamp (included)
- Spring balanced joints
- Optional clip-on lens for increased magnification (3x)



MGL150WM	Wall Model	<b>\$425.00</b>
MGL150MM	Mobile Model	<b>\$425.00</b>
MGL151	Add-on Lens	<b>\$155.00</b>

**MAGGY LAMP**  
medical

## LAMP EXAMINER LFM30-101

- Magnification and cool lamp
- Round 4 diopter lens
- Internal spring-balance arm with parallel movement
- Magnification ratio: x 100% actual size



**\$225**

LULKMAGW

**LUXO**

## HEINE mini-c® Clip Lamp



Compact, long-life pocket diagnostic light  
Sturdy, long-life pocket diagnostic light in a compact, modern design with chrome head. Bright, concentrated light thanks to novel XHL Xenon Halogen lens-lamp. Mini-c clip lamp only fits a mini-c handle. Not compatible with mini 3000 instruments.

- Compact, oval, attractive light
- Fits ergonomically in the hand
- High-quality handle: Shockproof, sturdy, non-slip
- Switches off automatically when put back in pocket
- Replaceable batteries. Size AAA



**\$94**

HED0173109

## Ear Light for mini-c Clip Lamp

6 pieces in blister pack

HED0073105	Pkt/6	<b>\$114.00</b>
------------	-------	-----------------



## MAG-LITE®

High Quality Examination Torches

### Mini Maglite

- High-intensity light beam
- 1/2 turn, twist focus, spot-to-flood
- Spare lamp safely secured in the tailcap
- High-strength aluminium alloy case
- O-ring sealed for water resistance



AAA

**\$29.95**

SH85200

AA

**\$39.95**

SH85700

The Mini Maglite® 2-Cell AAA flashlight is a smaller and lighter version for ease of use and storage.

### MagLite 3 Cell

A very popular D cell Mag-Lite. Featuring a long handle, extending its length to 31cm.



**\$49.95**


SH86230





## Diagnostic Penlight

- Modern satin finish aluminium body – delux model
- Pocket clip doubles as on/off switch
- Useful 5cm rule and dilation guide
- 2 x AAA batteries included



**\$24.<sup>50</sup>**

SAS233

ParagonCare

## WOODS LIGHT Ultra Violet

For the detection of various skin disorders

- Convenient hand size
- Rugged plastic construction
- Plastic lens over tube to protect patients' eyes
- Includes torch



**\$25**

LRUVL

Tube for Woods Light

LRUVT

**\$9.50**

## Promar Torches

ParagonCare

- Sealed unit with 2xAAA size batteries included - non rechargeable, non replaceable batteries
- Concentric beam bulb with pocket clip/switch
- Colour - White



	Each	+6
LRDDT	\$3.95	<b>\$3.50</b>

## Handheld LED Illuminated Magnifier

- 85 x 50mm lens
- Crystal clear acrylic lens with
- 2.0 X mag and 4.0 X mag bubble
- Complete with batteries and pouch
- Ergonomic shaped handle fits comfortably in the hand



**\$29**

MGL485

## MAGGYLAMP UV Woods 4x4 WATT

- Fitted with 4 x 4 Watt UV light
- UV Wavelength Lens: black light blue, 365nm biconvex glass, rectangular 100 x 50mm
- Magnification: 2.25 x 5 dioptre

**\$549**

RB215



MAGGYLAMP  
medical

This Catalogue has been printed on paper sourced from Sustainable Forests under the PEFC accreditation.





## LoupeLight 2 with mPack mini for i-View loupe mount (HR / HRP)

Compact LED HQ Headlight – now even brighter

The combination of light and loupe delivers the ultimate in performance as high quality illumination with LED HQ maximises the image quality of the loupe. For a perfect view and a precise diagnosis, the LoupeLight 2 is now 25% brighter than its predecessor model, with a homogeneous edge-to-edge illumination of typ. 55,000 lux.

For intensive use, the LoupeLight 2 now features a replaceable connecting cord which connects the device to the mPack mini – the new and compact power source for the LoupeLight 2.

- Compact and Lightweight
- Bright, homogeneous and coaxial illumination (typ. 55,000 lux)

### HEINE mPack Mini

The new, high performance power source HEINE mPack mini for LoupeLight 2 and MicroLight 2 offers freedom from cables that tie the users to desk, wall or static power sources. With a weight of only 95 g, this mobile power source is now approx. 70 % smaller and lighter than its predecessor.

HEC0832265	LoupeLight 2 F HR/HRP + mPack Mini	<b>\$2,275.00</b>
HEX0099007	mPack Mini Belt Clip	<b>\$34.00</b>
HEC0032240	Connecting Cord SC1 (1,5m/Ø 2,4mm)	<b>\$197.00</b>
HEC0032243	Connecting Cord SC2 1.5m (1,5m/Ø 3,2mm)	<b>\$197.00</b>



LED HQ  
LED NOW IN HEINE QUALITY.



## HEINE LoupeLight 2 Sets with HR Binocular Loupes

For a perfect view and a precise diagnosis, the LoupeLight 2 is now 25% brighter than its predecessor model, with a homogeneous edge-to-edge illumination of typ. 55,000 lux.

For intensive use, the LoupeLight 2 now features a replaceable connecting cord which connects the device to the mPack mini – the new and compact power source for the LoupeLight 2.

- Compact and lightweight
- Bright, homogeneous and coaxial illumination (typ. 55,000 lux)
- Optional polarisation filter system

#### Includes:

- LoupeLight2
- mPack mini
- mPack mini belt clip
- Yellow Filter
- E4 USB Plug-In Transformer
- 2 x Protective Lenses
- 2 x Sterilisable Swivel Levers
- Retaining Cord
- Cleaning Fluid
- Cleaning Cloth
- Carry Case



HEC0832450	HR 2.5x/340mm, LL2, I-View A.S-Frame	<b>\$4,655.00</b>
HEC0832451	HR 2.5x/420mm, LL2, I-View A.S-Frame	
HEC0832452	HR 2.5x/520mm, LL2, I-View A.S-Frame	

LED HQ  
LED NOW IN HEINE QUALITY.





## ML4 LED HeadLight with Headband Professional L and Power Sources

**HEINE** QUALITY  
MADE IN GERMANY

Ultra Bright Coaxial Illumination – Now with LED <sup>HQ</sup>

- **LED in HQ** - the new standard in LED illumination defining optimal light intensity, homogeneity and colour rendering for each individual examination situation.
- **Adjustable illumination spot size** - 30 mm to 80 mm spot size range (at 420 mm working distance) to fit any examination situation.
- Bright and homogenous Illumination. Absolutely bright light spot that is uniform from edge-to-edge for the perfect illumination in all examination situations.
- **Stepless light intensity control** - Optimal brightness setting prevents reflexes.
- **Flexibility**. Flexible power source options: 100% mobility with a choice of cable-free mPack UNPLUGGED head-worn battery or mPack belt-worn battery pack.
- **Coaxial Illumination** - Coaxial design ensures a completely shadow-free image and allows for excellent illumination of difficult to see areas.
- **Comfortable and secure fit** - The ML4 Professional headlight has multiple adjustment points and soft padding for comfort and stability – even during long examinations.
- **Optional Filters** - Polarisation filter P 2 improves contrast and yellow filter reduces blue light.



HEJ0831413	ML4 Headlight with mPack	<b>\$2,795.00</b>
HEX0799666	mPack Unplugged + Plug Transformer	<b>\$1,095.00</b>

**\$2,950**

HEJ0831415



## HEINE ML4 LED Headlight with Binocular Loupes

Combining HEINE high performance binocular loupes with the ML4 headlight to provide ideal magnification and illumination for surgical and examination applications.



HEJ0831530	Headlight ML4 Unplugged with 2.5x 340mm Loupes	<b>\$4720</b>
HEJ0831531	Headlight ML4 Unplugged with 2.5x 420mm Loupes	<b>\$4720</b>
HEJ0831532	Headlight ML4 Unplugged with 2.5x 520mm Loupes	<b>\$4720</b>



## HEINE mPack

- **Portable power pack** - Does not tie the user to desk, wall or static power supplies.
- **Li-ion Technology** - Triple capacity compared with current NiCd systems. No memory effect, autonomy for 8.5 hours of continuous ON time at full power.
- **Fast Charge** - Charges to full capacity in 2 hours.
- **Charge Status Display** - Always lets you know how much operating charge you have left.
- **Stepless light intensity control** - Optimum brightness setting prevents reflexes.
- **Flexible Charging Options** - EN 50 Desk and wall charger, or plug-in transformer.
- **Intelligent** - Can be used as mains power supply even while in "charge" mode.
- **100 – 240 V Compatible Transformer** - incl. 4 plugs for world wide compatibility.



10% GST will apply to all items unless they are marked GST Free



## HEINE S-Frame + HR<sup>®</sup> 2.5 x Binocular Loupe with i-View

### HEINE S-Frame<sup>®</sup>

- Innovative spectacle frame for HR Binocular Loupe
- Appealing and functional design
- Can be customised to fit any shape of head and face

### HEINE HR<sup>®</sup> Binocular Loupe 2.5X

- The brightest image – Heine's exclusive coating technology reduces reflection to 0.3%, providing the brightest image of any comparable loupe
- Crispest image with colour correction
- Super-wide, super-deep field of view
- Dust and water resistant



HR Binocular Loupe Set with i-View Optics, S-Frame with i-View loupe mount for HR 2.5 x loupes, detachable protective lenses, 2 sterilisable swivel levers, retaining cord, cleaning fluid, carrying case.



HEC0032355	340mm working distance	<b>\$2,190.00</b>
HEC0032356	420mm working distance	
HEC0032567	520mm working distance	

## Mag-eyes The hands-free, head mounted magnifier

Mag-Eyes with 2 Lenses  
(1.6/2.0 X Magnification)

**\$79**

MGL400

Working distances (Lens to Subject)	
1.6 X (30-35cm)	2.0 X (25-30cm)
2.25 X (18-25cm)	2.75 X (10-18cm)

Mag-Eyes with 2 Lenses  
(2.25/2.75 X Magnification)

**\$95**

MGL410



## Guthrie Loupe

ParagonCare

- Adjustable head band padded for user comfort
- Plastic, Precision Lens
- Retractable visor
- 2.7x Magnification

**\$95**

LRM2001



Here's a bright idea!  
Why not visit  
[www.vitalmedicalsupplies.com.au](http://www.vitalmedicalsupplies.com.au)  
to place your lighting  
order today.