

MEDICAL REFRIGERATION

VACCINE RANGE 2019

- Operating temperature range: +2°C to +8°C
- Low energy consumption
- Fixed lock standard
- Visual and audio alarms
- Dual temperature probes - recording both cabinet air and simulated product temperatures for a more accurate reporting
- 2 year warranty (4 year extended warranty package also available)



Cloud based temperature monitoring



Remote access



Text and email alerts



Reporting



Battery backup



Medisafe Plus 371

First 12 months
cloud
monitoring
FREE!



Compliant Pharmacy Refrigeration System

Features

- Superior Quality; No domestic fridge components used
- Electronic temperature control and monitoring system ensures consistent temperatures are maintained on all shelves
- Dual temperature probes - recording both cabinet air and simulated product temperate for a more accurate reporting.
- Built-In cloud based temperature monitoring system (system contains its own independent wireless connection and does not rely on stores having their own WIFI connection)
- Provides continuous temperature monitoring.
- Email and SMS alerts notifying of temperature breach and power outage
- Min/Max Temperature Recorder logger for ease of daily temperature recordings Visual and audible auto resetting alarms and a power failure alarm
- VI311MS and VI371MS models include an internal evaporator guard and door fan switch
- All fridges come standard with an integral lock and key
- 24 hour technical support & service
- 24 months warranty
- Australian design
- First 12 months cloud monitoring is free, after which you will have the option to renew for a paid service of 12 or 36 months cloud monitoring.




1500L

\$9,995

VI1365MS

Height	210cm
Width	150cm
Depth	82cm
Litres	1500



CLOUD BASED REPORTING INSTALLED

*Requires a stable 4G mobile data signal to operate AQ box cloud monitoring



First 12 months cloud monitoring FREE!



80L

\$2,595

VI80MS

Height	83cm
Width	50cm
Depth	56cm
Litres	80



306L

\$3,395

VI311MS

Height	163cm
Width	59.5cm
Depth	59.5cm
Litres	306



381L

\$3,695

VI371MS

Height	193cm
Width	59.5cm
Depth	59.5cm
Litres	381



700L

\$6,995

VI625MS

Height	210cm
Width	75cm
Depth	82cm
Litres	700

Freight charges will apply



Medisafe AKG 157

Features

- Electronic temperature control & monitoring system
- Visual and audible auto resetting alarms,
- Min/Max temperature Recorder
- Dual temperature probes - recording both cabinet air and simulated product temperature for a more accurate reporting.
- Noise level of 35.5 db
- Low energy consumption
- Australian design
- 24 months warranty
- Power failure alarm
- Power supply 10 amp
- Built-In cloud based temperature monitoring system (system contains its own independent wireless connection and does not rely on stores having their own WIFI connection)



\$2,995

VI157MS



Height	824mm
Width	595mm
Depth	626mm
Litres	138

Bromic MediFridge™ Display Chiller 145L

- Integrated temperature Alarm System ensures sensitive stock is maintained and managed according to requirements
- Fan Forced air circulation for even temperature throughout the unit
- Increased air circulation around the unit thanks to unique custom wire 'open rung' shelves



\$1,695

BR140GD



Model	VI80MS	VI157MS	VI311MS	VI371MS	VI625MS	VI1365MS
Dimensions External (W x D x Hmm)	500x560x830	595x626x824	595x595x1630	595x595x1930	750x820x2100	1500x820x2100
Gross Litres	80	138	306	381	700	1500
Weight	53	76	82	92	185	280
Lock	Y	Y	Y	Y	Y	Y
Alarm	Y	Y	Y	Y	Y	Y
Power	10amp	10amp	10amp	10amp	10amp	10amp
Cloud Monitoring	Y	Y	Y	Y	Y	Y
Refrigerant	R600	R134A	R134A	R134A	R404A	R404A
Temp Probes	2	2	2	2	2	2
BMS Outlet	Y	Y	Y	Y	Y	Y

Medi-Temp Fridge/Freezer Digital Thermometer

Features

- 2 temperature sensors: 1 inbuilt inside the unit and another sensor on a 3 meter wire.
- Inbuilt alarms that can tell you that the temp limits are exceeded.
- Inbuilt Memory for temperature history readings.
- Inbuilt sensor range measures from: -10 up to +50C and,
- Outer sensor (on the end of the cable) measures temperatures from -50 C up to +70C.

\$29

HLT009



Medilog Fridge/Freezer Thermometer

Features

- Pre set alarm set points for the vaccine requirements
- Inbuilt Memory for temperature history readings.
- Includes a 3 meter long small wire probe and a stick on probe holder. Temp range -20C to + 70C. (C/F Selectable).
- Battery powered - Immune to power outs.

\$74

HLML0G

