



Warts.

Safely Treated.
Kid's Hands Approved.

Histofreezer® Portable Cryosurgical System

- **Safe**
Low risk of scarring, safely used on children 4+ years old^{1,2}
- **Clinically proven**
93.4% resolution rates compared to liquid nitrogen³
- **Easy-to-use**
Portable, 1-minute, one-handed procedure⁴


Proven Performance*



*Occasionally multiple treatments spaced 2 weeks apart may be required. Verruca vulgaris shown.

Easy 3-Step Procedure

1. Charge 2. Invert 3. Apply



Histofreezer®

Tough on warts. Gentle on skin.

References 1) Shepherd, J.P and Dawber, R.P.R. Wound Healing and Scarring After Cryosurgery. *Cryobiology* 1984;Vol. 21:157-169. 2) Ricketti, J.C. and Niedbala, R.S. A Study of a New Method of Cryosurgical Treatment of Verrucae Plantaris. *The Lower Extremity* 1995;Vol. 2, No. 3. 3) Martínez, F.C. et al. Dermatological Cryosurgery in Primary Care with Dimethyl Ether Propane Spray in Comparison with Liquid Nitrogen. *Atención Primaria* 1996; Vol. 18, No. 5:211-216. 4) Histofreezer Portable Cryosurgical System Package Insert. ©2017 OraSure Technologies, Inc. Histofreezer is a registered trademark of OraSure Technologies, Inc. HISTO0376 02/17



220 East First Street
Bethlehem, PA 18015 USA
www.OraSure.com
www.Histofreezer.com



CRYOPRO MINI FLASK 350ML

**OVER 2 DECADES OF PROVEN RELIABLE SERVICE TO
DERMATOLOGISTS AND MEDICAL PRACTITIONERS.**

- Suitable for use by both GP and Specialists
- Comes with Applicator Set (5 Tips and 1 Cleaning Adaptor)
- 3 Year Guarantee
- Simple & Reliable
- 12 - 14 Hour Static Holding Time



CryoPro Mini Flask 350mL

Allows Primary Care Doctors to quickly and easily treat a range of dermal problems including Keratoses, Condaylomata, warts, skin tags and other skin lesions.

Danish design and manufactured, the CryoPro has provided over 2 decades of proven reliable service to Dermatologists and Medical Practitioners.

- Suitable for use by both GP and Specialists
- 350ml Capacity
- Comes with Applicator Set (5 Tips and 1 Cleaning Adaptor)
- 3 Year Guarantee
- Simple & Reliable
- 12 - 14 Hour Static Holding Time



	Each	3+
ELCPRO 350ml Flask	\$934	\$853

BRYMILL'S CRY-AC® Hand-Held Liquid Nitrogen Delivery System

Features:


- Stainless steel and solid brass constructions for long life
- Slender, fingertip trigger for maximum control and ease of viewing
- Wide opening neck for easy filling
- Streamlined relief valve maintains a constant operating pressure to ensure consistent and accurate freezing
- Patented Safety Autovent maximises safety by allowing internal pressure to gradually vent as the cap is unscrewed
- Sturdy base for extra stability
- 3 year warranty

Spray kit includes:

- Flask 300ml or 500ml
- 5 spray apertures with varying orifice size
- 1 x 78mm curved spray extension
- Instruction and technique booklet



DSBRY 300ml Flask	\$1,395.00
DSBRZ 500ml Flask	\$1,495.00





Histofreezer® Advanced Measurement Technology

The Histofreezer® Portable Cryosurgical System is a low cost, highly effective cryosurgical agent. This therapeutic treatment makes use of local freezing for the controlled destruction or removal of living tissue.

Key Points:

- Ultra-portable, one-handed operation, treatment takes less than 60 seconds to perform
- Results comparable to those achieved with liquid nitrogen
- Requires no bulky equipment or large investment, unlike Liquid Nitrogen Systems
- High patient acceptance, relatively painless and bloodless with minimal threat of scarring
- Highly efficient cryogen delivery (low waste) and no messy or leaky cones
- 2 x 80mL complete Histofreezer (canister + applicators)

Histofreezer®



	2x 80mL	5+
HM1008301	\$199	\$169

Liquid Nitrogen Storage Container

Features:

- Comes with 5 year vacuum warranty, 2 year container warranty
- Made of high-strength aluminum alloy.
- Adopting advanced powder coating technique
- Each unit comes complete with safety lock lid, neck plug, leather (PU) bag and ladle.



Freight charges apply

Ladle outer diameter	Unit	10L	20L
Ladle height	L	10	20
Number of ladles that can fit per dewar	kg	6.0	9.5
Comes with one ladle per dewar	mm	50	50
Ladle outer diameter	mm	300	394
Ladle height	mm	552	672
Number of ladles that can fit per dewar	L/d	0.12	0.12
Comes with one ladle per dewar	days	86	168
Ladle	Ladle outer diameter	mm	38
	Ladle height	mm	120
	Number of ladles that can fit per dewar	pc	3~6
	Comes with one ladle per dewar	pc	1
0.5ml	pc	792	792
0.25ml	pc	1788	1788
0.5ml	pc	-	1284
0.25ml	pc	-	2832
Code		TW10LD	TW20XTD
Price		\$1,121	\$1,320

CryoPen B Kit

Ideal for General practice and veterinary

- Now includes 3 Micro applicator tips: gold dot 1-4mm, silver dot 2-8mm & bronze dot 4-12mm micro-applicators.
- To be used with 23.5g cartridges



\$2,365

CPBPLUS

CPBP-C235	Cartridge for Cryopen BPlus 6 x 23.5g Suitable for Cryopen-B model	\$322
-----------	--	--------------

Cryopen-O Cartridges

- For Cryopen-O 16g

CPC16	Pkt/6	\$106.43
-------	-------	-----------------

- For Cryopen-O 8g

CPC8	Pkt/24	\$210
------	--------	--------------

Silver Nitrate Applicators


- 75% Silver Nitrate
- Single Use 15cm Bendable Applicator

	Pkt/100	6+
SNA7484	\$86.90	\$75.90



Roller Base

This sturdy, cast-aluminium, five-wheel roller base enables you to move dewars easily and safely. To fit 396mm diameter.



\$479

FCRB400

Ladles

Ladle 20 Litre with handle

\$190.⁸⁰

TWDIP3

Ladle 10 Litre with handle

\$190.⁸⁰

BIOLD4



Dewar Measuring Rod

Low conducting measuring rods that clearly show the LN2 level of a dewar in inches and centimetres. (92cm - 36inch)



\$19

FCMR033

Universal Extractor

Adjustable to suit dewars from 10L - 25L.



\$199

WMDEXT

