

We invented the dermatoscope 30 years ago. Now we've done it again.

The new Dermatoscopes: DELTA 30 and DELTAone.



30
YEARS
1st DERMATOSCOPE



THE NEW HEINE DELTA 30 DERMATOSCOPE.

**If there's something to discover,
you'll be sure to see it.**



30
YEARS
1st DERMATOSCOPE





HEINE DELTA 30

If there's something to discover, you'll be sure to see it.

30
YEARS
1st DERMATOSCOPE

Not only would our new high-end dermatoscope win any beauty pageant; it stands out due to the best imaging we've ever developed. Combined with a field of view measuring a veritable 30 millimetres and unprecedented colour rendering, in LEDHQ of course. Polarisation provides glare-free and non-reflective working conditions. Above all else, it is comfortable to operate with its angled, ergonomic design and intuitive control. Feel free to take a look – you're in for a real discovery.

- Excellent colour rendering due to LEDHQ
- Digital documentation with the exclusive HEINE DERM app or the HEINE Cube System**
- Ergonomic handling thanks to the angled product design
- Crystal-clear image with the new achromatic HEINE optical system
- Largest field of view of real 30 mm (lens 32 mm)
- Glare-free and non-reflective working conditions thanks to polarisation
- Custom light intensity: 3 degrees of luminous intensity
- 10x magnification
- Switch between polarised and non-polarised light at the touch of a button = toggle function. For more information, please visit heine.com/en/Toggle
- No cross-contamination – autoclavable contact plate with scale
- Maintenance-free due to our dustproof design
- Dioptre compensation
- Rechargeable battery with Li-ion technology: no memory effect; you can safely recharge the battery regardless of its charge state
- Wireless connection via iPhone* and WLAN with the HEINE Cube System**



\$2,499

HEK2302830

NEW

* The Apple hardware is not included in the scope of delivery. For Apple hardware, please contact an authorized Apple distributor, an Apple Premium Reseller or Apple Inc. directly. Apple, the Apple Logo, iPhone and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries.**Not available in all countries.



Mobility and quality have now become *one*. The new HEINE DELTA*one* Dermatoscope.



30
YEARS
1st DERMATOSCOPE

This mobile masterpiece of dermatoscopy is taking medical professionals by storm. The crystal-clear image was accomplished using the new achromatic HEINE optical system. What's more, our LED^{HQ} produces extraordinary colour rendering. Plus the polarisation provides glare-free and non-reflective working conditions. The new DELTA*one* is a tool you'll gladly carry in your pocket.
www.heine.com





HEINE DELTAone

Mobility and quality have now become one.



This mobile masterpiece is taking medical professionals by storm. The crystal-clear image was accomplished using the new achromatic HEINE optical system. What's more, it produces extraordinary colour rendering, due to the quality of our LEDHQ. Examinations are now particularly comfortable due to the polarisation allowing doctors to work glare-free and without disrupting reflections. These three features – clarity, LEDHQ and non-reflective working – are what turn the DELTAone into a vital working tool wherever you go.



- Field of view of real 21mm (lens 22 mm)
- Custom light intensity: 3 degrees of luminous intensity
- 10x magnification
- Switch between polarised and non-polarised light at the touch of a button = toggle function. For more information, please visit heine.com/en/Toggle
- No cross-contamination > autoclavable contact plate with scale
- Maintenance-free due to our dustproof design
- Dioptre compensation
- Rechargeable battery with Li-ion technology: no memory effect; you can safely recharge the battery regardless of its charge state
- Wireless connection via iPhone* and WLAN with the HEINE Cube System**
- Excellent colour rendering due to LEDHQ
- Fits into any pocket
- Crystal-clear image with the new achromatic HEINE optical system
- Glare-free and non-reflective working conditions thanks to polarisation
- Digital documentation with the exclusive HEINE DERM app or the HEINE Cube System



HEK2102830	DELTAOne Standard	\$1,799
HEK2112830	DELTAOne White	\$1,859

* The Apple hardware is not included in the scope of delivery. For Apple hardware, please contact an authorized Apple distributor, an Apple Premium Reseller or Apple Inc. directly. Apple, the Apple Logo, iPhone and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries.**Not available in all countries.





DELTA 20T Rechargeable Dermatology Set

The Specialist

The HEINE Delta 20T offers unrivaled German optics and the latest LEDHQ illumination creating the most precise diagnostic environment. Confidently diagnose even the most intricate lesions in both polarised and non-polarised modes.

- >16x Magnification
- Polarised and non-polarised
- SLR digital documentation available
- Optional small contact plate available
- 5 year warranty

Includes:

- HEINE DELTA 20T Dermatoscope Head
- Contact plate with scale
- BETA 4 NT rechargeable handle
- NT 4 table charger
- Dermatoscopy Oil 10 ml
- Hard case

\$3,279

HEK6224420



DELTA 20T USB Rechargeable Set

Includes:

- HEINE DELTA 20T Dermatoscope Head
- Contact plate with scale
- BETA 4 USB rechargeable handle
- USB Cord and plug In power puppy
- Dermatoscopy Oil 10 ml
- Hard case

\$2,866

HEK6228388



NC2 Dermatoscope

LEDHQ
LED NOW IN HEINE QUALITY.

The All-Rounder

The new HEINE NC2 is the all-in-one Dermatoscope which adapts perfectly to a dynamic and varied clinical environment. Polarised, contact, non-contact and digital documentation are all covered with this device.

- 6x - 10x Magnification
- Polarised
- Contact & Non-Contact
- iPhone 6/7/8 Digital Documentation
- 5 year warranty



		Each
HEH8417812	NC2 Dermatoscope with C/Plate and zipper case	\$1,704
HEH8427812	NC2 Dermatoscope without C/Plate or zipper case	\$1,557
HEH0078126	iPhone 6 case for NC2 Dermatoscope	\$101
HEH0078127	Mobile Phone Case NC2/7 iPhone 7	\$101

iC1 Dermatoscope

LEDHQ
LED NOW IN HEINE QUALITY.

The Digital Doctor

The HEINE iC1 is a completely digital Dermatoscope suitable for clinicians where digital documentation and record keeping is of paramount importance. Simple, wireless and great for patient education.

- 40 x Digital Magnification
- Polarised & Non-Polarised
- Contact
- iPhone 6/7/8 Digital Documentation
- 5 year warranty



		Each
HEK2702830	iC1 Dermatology Set (6/6S)	\$1,689
HEK2732830	iC1 Digital Dermatology Set (7/8)	\$1,689

mini3000 LED Dermatoscope

The Locum

For the medical professional who values maximum quality and minimum size. Conveniently compact, lightweight and with a unique on/off switch, the HEINE mini3000 series makes an ideal portable set.

- High quality achromatic optics with 10x Magnification
- Non-polarised
- Solid metal construction
- Superior LEDHQ illumination
- Integrated pocket on/off switch
- 5 year warranty

LEDHQ
LED NOW IN HEINE QUALITY.

\$836

HED8778021



Vital Medical Supplies New Website

The top 7 reasons for change

#1

The Shopping Cart has an auto save function and will retain your order even if you log out.

#2

The website is built on a very stable platform, which means minimal down time.

#3

The website will no longer time you out.

#4

There is more space on the website to include better images and product content.

#5

You can learn more about each product with our easy to find search engine.

#6

We can assist you in simplifying and streamlining your orders through an order template custom for you.

#7

The new website has a fresh new look.

www.vitalmedicalsupplies.com.au



Why choose **Vital Medical Supplies**

**Price Match
Guarantee**



Efficient Delivery



**Over 17,000
Products**



**5 Star
Customer Care**



and so much more...

 Vital Medical
Supplies

www.vitalmedicalsupsplies.com.au



HEINE DELTA 30 and DELTAone Accessories



		Each
HEK0034260	AD. Smartphone Delta 30/One/6	\$119
HEK0034261	AD. Smartphone Delta 30/One/7+8	\$119
HEK0034262	AD. Smartphone Delta 30/One/X+XS	\$119
HEK0034263	AD. Smartphone Delta 30/one/XR	\$119
HEX0029921	Charger Delta 30	\$445
HEK0034220	Contact Plate Delta 30	\$289
HEK0034222	Contact Plate Delta One	\$269
HEK0034206	Small Contact Plate Delta 30	\$289
HEK0024207	Small Contact Plate Delta One	\$269
HEK0034270	Universal Smart Phone Charger	\$99



HEINE iC1 Accessories



		Each
HEX0099303	USB cord + E4 Medically approved Plug-In power supply	\$101
HEK0034201	Contact Plate Heine iC1 with scale	\$374
HEK0034250	Mobile phone case for iPhone 6 Series	\$101
HEK0034253	Mobile phone case for iPhone 7/8 Series	

DELTA 20 T Accessories



		Each
HEK0034205	Contact plate small, 8mm dia. (for the examination of inaccessible lesions)	\$233
HEX0099008	BETA Belt Clip	\$61
HEX0009923	Adaptor Cord (1.3m)	\$99
HEK0034186	Nikon SLR Photo adaptor set for Delta 20/Plus/T	\$895
HEK0034185	Cannon SLR Photo adaptor set for Delta 20/Plus/T	\$895

DELTA 20T Dermatoscopy Head Only



HEK0834221	Contact Plate with Scale	\$2,282
------------	--------------------------	----------------