

Drapes

Sterile Drapes



Fenestrated Poly-lined Towel

		Box/100
31120079	Absorbent, 18x26 in. (45x66cm), 2.5 in aperture, sterile	\$139

Drape Sheets

		Box/20
31120123	Tiburon® Medium Sheet, 40 x 71 in. (101 x 180cm)	\$89
31120124	Tiburon® Half Sheet, 40 x 57 in. (101 x 145cm)	\$79
31120143	Tiburon® Utility Drape Sheet, 44 x 60 in. (112 x 152cm), 4 x 6 in. (10 x 15cm) poly-reinforced fenestration	\$122

Drape Sheet, Square-Folded

		Box/30
31120142	Tiburon® Drape Sheet, square-folded, 38.5 x 38.5 in. (98 x 98cm)	\$119

Poly U-Drapes

31120108	Blue poly U-drape, 60 x 84 (152 x 213cm)	Box/22	\$259
31120109	Blue poly U-drape, 60 x 70 in. (152 x 183cm) split with adhesive	Box/24	\$209

Leggings

		Box/20
31120104	Leggings, folded, 6 in. cuff, 31 x 48 in. (79 x 123cm), 2/pkg.	\$159

Under-Buttocks Drapes

		Box/20
31120114	Under-buttocks drape with graduations, fluid control pouch with filtering screen and drainage port, plastic backing	\$215

Aperture Drapes

		Box/10
31120156	Drape Aperture Sml Surgical	\$46

Minor Procedure Drapes

		Box/25
31120163	Drape Minor Procedure	\$169

Stockinettes

		Box/10
31120004	Stockinette, X-large 12 x 48 in. (31 x 122cm)	\$89



A FULL RANGE OF CARDINAL HEALTH SURGICAL DRAPES ARE AVAILABLE THROUGH VITAL MEDICAL SUPPLIES. FOR INFORMATION PLEASE CALL OUR CUSTOMER SERVICE ON 1300 557 651.

Equipment, Drapes and Covers




Mayo Stand Covers

		Box/30
31120095	Mayo stand cover, reinforced poly, 23 x 55 in. (58 x 140cm)	\$119

Heavy-Duty Back Table Covers

		Box/24
31120094	Heavy-duty back table cover, 60 x 90 in. (152 x 229cm), thick SMS reinforcement, 29 x 90 in. (74 x 229cm)	\$249



Sterile Drapes



Soft, drapable DuPont Sontara® laminated with PE to give it 100% (impervious) fluid-spray rating on one side and absorbent on the other.



Surgical Drape Plain		Each	100+
Waterproof			
GM29-220	60x60cm	\$2.26	\$1.93
GM29-223	90x90cm	\$3.91	\$3.47

Surgical Fenestrated Drape		Each	120+
Waterproof with 7.5 x 6cm fenestration			
GM29-612EE	45x60cm	\$1.82	\$1.49
GM29-614EE	45x75cm	\$1.98	\$1.76
		Each	100+
GM29-620EE	60x60cm	\$2.09	\$1.76
GM29-634EE	90x90cm	\$4.24	\$4.13

Minor Procedure Drape		Each	100+
Waterproof with 7.5 x 6cm fenestration and 10 x 10cm adhesive surround			
GM29-780EE	60x60cm	\$3.41	\$3.03

Plastic Waterproof Sheets		Each	100+
GM22-200EE	45x60cm	\$0.61	\$0.50
GM22-201EE	61x91cm	\$0.77	\$0.55
GM22-203EE	90x120cm	\$0.99	\$0.77
GM22-205EE	120x140cm	\$1.38	\$1.10

Dressing Towel Sterile


		Each	100+
GM30-666EE	40x40cm	\$0.63	\$0.50

3M™ Steri-drape™



Small Absorbent Prevention Fabric drape with adhesive area and aperture

Fabric is impervious to fluid strikethrough over the entire drape surface. It is also absorbent, has superior strength, and is soft and conformable.



		Box/25	2+
3M1092	56 x 64 cm	\$130.91	\$123.21

Large range of drapes available on request