



Projet 101 Pulsatile Ear Irrigator

Designed by Doctors with innovative pulsatile irrigation

Utilising the most modern production methods the Guardian Projet 101 was designed to fulfil the most demanding surgeries and hospital ward requirements for a reliable, leak free and easy to use machine that is not expensive to run and does not require constant maintenance and subsequent down time.

- Extra long mains lead and handle hose so machine can be sited where it is most convenient
- Integral wall mount frees up valuable work space by wall mounting the Projet™ however easily removed for cleaning etc
- Sliding flow regulator conveniently mounted in handle. Water flow can be temporarily stopped by sliding regulator to off position or flow regulated by sliding control to desired position
- Ultra thin atraumatic tips help visualisation and direction of water flow
- Digital control of water pressure and flow is easy to operate and remains constant and reproducible

\$445

PROJET101

**2
Year
Warranty**



PROJET-TIP Tips Box/100 **\$165.00**

Projet




Bionix Ear Safe Curettes

The Lighted Ear Curette™ with Magnification

Bionix Lighted Ear Curettes combine illumination and magnification to help healthcare providers take less time to see and safely remove cerumen from the ear canal.

Includes:

- All tip styles (except CeraPik™)
- 1 Light Source
- 1 Mag. Lens
- Includes 50 pieces



NEW

\$147

EC2245

Safe Ear Curettes

Bionix Safe Ear Curettes allow healthcare provides to safely and gently remove cerumen with styles to match different types of wax and personal preference

Includes:

- 10 of each style
- 15 white
- Includes 75 pieces



NEW

\$103.⁹⁵

EC5777



Ear Syringe NJ Phillips Precision Instrument

www.njphillips.com.au

For use in the removal of waxy build-up in the ear

- Ergonomic design
- Extremely accurate
- Robust, reliable and durable
- Easily adjusted
- Calibrated from 0.5ml to 5ml
- Silicone tips

\$295

MCS1043



ASSEMBLING METAL TIP UNIT:

1. Loosen needle nut until reasonably loose.
2. Slide needle tip base into the lower slot of the needle nut.
3. Screw the needle nut tightly back onto the main threaded mount.

ASSEMBLING RUBBER TIP UNIT:

4. The rubber tip **MUST** always be sterilised prior to use.
5. Remove the metal tip and needle nut.
6. Push the rubber tip onto the main threaded mount. Turn clockwise as you push to ensure it is firmly in place.

Spares

MCNITIS020	Nozzle only	\$39.95
MCNIWX1187	Repair kit	\$39.95
MCS137	Rubber tips Pkt/25	\$52

Ear Syringe With O Ring

ParagonCare

- Ring fitted on plunger giving smooth movement
- Unit supplied with appropriate nozzles (2x)



Spares

LRRESORT	O-ring for nozzles	\$6.75
LRLESOR	O-ring for plunger	\$8.95
SAS018	Nozzle Set (Pair)	\$10.50
SAS001	Shield	\$5.50

\$75

LRLESWOR

Ear Syringe

With Dual O Rings, 2 Nozzles & Shield



Dual Tips

\$52

SY2010

SAS010R	Rubber O ring for tip	\$3.75
---------	-----------------------	---------------

Ear Syringe

Metal Premium German Quality



\$189

HTESYR



Propulse® Irrigator Accessories

KIT0099	Propulse Cleaning Tablets	Box/200	\$79
SWCU08036	Noots Tank Disposable	Box/40	\$23.95
KIT0029	Propulse Ear Wax Removal Scoop	Box/40	\$139
KIT1012	Propulse Protect Cape	Box/36	\$39.50
KIT5002B	Propulse QrX™ Disposable Tips	Box/100	\$195



You can place your order via fax on **1300 557 631**.