



DERM A942™ High Frequency Desiccator

The world's most user-friendly diathermy



The DERM 942 and its patent pending handpiece are designed and engineered to be the most reliable and durable desiccator available today. The unit's life is prolonged by incorporating the latest in chip technology and reducing the number of internal parts, thereby generating less internal heat.

- 40 watt high frequency desiccator
- Adjustable in 1/10th watt increments in the first 10 watts, in monopolar and bipolar
- More accurate power with Universal Power Control 100-240 VAC
- Less to do and less to buy makes this the most cost effective and user friendly unit on the market
- The Bovie DERM 942 is a ground reference unit which means no grounding pad is required.



\$1,995

DSA942



Bovie Stand

Table-top stainless steel stand

\$99

DCA813

Bipolar Accessories

Bipolar power is intrinsically safer than monopolar HF power

| | | |
|----------|---|--------------|
| DC701BA | Bipolar Adson Forcep 1mm wide serrated jaw | \$379 |
| DSA827EU | Bipolar Cable | \$225 |
| DSA803 | Footswitch | \$355 |



Aaron Disposable Non-Sterile Electrodes

| (Fits Hyfrecator) | | Box/100 |
|-------------------|-------|--------------|
| DSA804 | Sharp | \$125 |
| DSA806 | Blunt | \$125 |




Aaron Disposable Non-Sterile Pencil Sheaths

Box/100

\$65

DCDA910



Aaron Disposable Cautery

- Complete battery cautery handle with fine point tip
- Individually sterile packaged

| | | |
|--------|--------|--------------|
| DCAA01 | Box/10 | \$329 |
|--------|--------|--------------|



High temperature cauteries are popular for many types of procedures including pinpoint hemostasis



The Hyfrecator® 2000 Electrosurgical System



Features:

- Easy to use and affordable
- 3 programmable power settings
- Remote control hand piece
- The medical industry's leading device for over 60 years
- Cost effective alternative to various laser procedures
- 12 months warranty

There is only one Hyfrecator®

Special features include:

Subtle control to 1/10 of a watt in low mode. Dual microprocessors: one to control the function and output of the Hyfrecator®; the second to terminate activation if a hazard or safety issue is detected.


Hyfrecator® 2000 Contents:

- Pack of 5 Electrolase Sharp Tip disposable electrodes, non sterile
- Pack of 5 Electrolase Blunt Tip disposable electrodes, non sterile
- 1 x Hyfrecator® Electrosurgical Handpiece
- 2 x Hyfrecator® Pencil Sheath, non sterile
- 2 x Handpiece Sheath, sterile
- Inservice Training DVD
- Operator's Manual
- Power Cord



\$1,995

MQ25-5011



Hyfrecator® 2000 Accessories

Below are listed most common lines. Full range available.



| | | | | |
|-----------|-----------------------------------|----------|-------|--------------|
| MQ25-5070 | Foot Control for Hyfrecator® 2000 | 7-900-6 | Each | \$275 |
| MQ25-5072 | Pencil Up/Down 2000 | 7-900-5 | Each | \$555 |
| MQ25-5076 | Foot Switch | 7-900-4 | Each | \$325 |
| MQ25-5200 | General Purpose Electrode | 711B | Each | \$69 |
| MQ25-5250 | Electrodes 714 | 714 | Box/6 | \$299 |
| MQ25-5313 | Bipolar Forceps 10cm | 7-809-1 | Each | \$725 |
| MQ25-5316 | Bipolar Forceps 13cm | 7-809-4 | Each | \$695 |
| MS255317 | HYF Forcep Adson Serrated Tip | 7-809-7 | Each | \$450 |
| MQ25-5320 | Bipolar Forecp Lead | 7-809-12 | Set | \$295 |

Electrolase Disposable Electrodes



**TOP
SELLER**



| Non-Sterile | | | Box/100 |
|-------------|-------|----------|--------------|
| MQ25-5350 | Sharp | 7-100-12 | \$110 |
| MQ25-5351 | Blunt | 7-101-12 | \$110 |
| Sterile | | | Box/50 |
| MQ25-5352 | Sharp | 7-100-8 | \$110 |
| MQ25-5353 | Blunt | 7-101-8 | \$110 |

Disposable Pencil Sheaths

| | | | |
|-----------|----------------------------|---------|--------------|
| MQ25-5105 | Pencil Sheaths Non-Sterile | Box/100 | \$55 |
| MQ25-5106 | Pencil Sheath Sterile | Box/25 | \$189 |

