

# SWITCH TO INTELLIGENT INFECTION PREVENTION

## THE Vernacare COMPACT RANGE

Consider a disposal unit for your single use / single patient environment

### COMPACT



- Energy efficient
- Quiet operation
- Multiple installation options
- EBOS Code: 22441073

### COMPACT+



- Prevents aerosolisation
- Minimise spillage and splashback
- Easy to clean
- Saves space

#### Other Optional features include:

- Reports in real-time
- Access Control/Key
- Multiple Cycle Options (Economy/ Heavy weight)
- EBOS Code: 22441086

## Make the right choice

### COMPACT



### COMPACT+

Small wards or single patient environments

Ward size & type

Large wards with multiple patient beds



Max 1 item per cycle

Items per cycle

Max 4 items per cycle



Low number of immobile patients who are toileted or washed in bed

Patient type

High number of immobile patients who are toileted or washed in bed



Patient bathroom, dirty utility room or small sluice room  
*Wall-mounted to save space.*

Installation location

Large dirty utility room or sluice room  
*Compact+ is floor-standing and needs additional space for installation.*

Contact our Account Managers for your Quote today!  
Call 1300 557 651



# The Vernacare Disposable System

Exclusive to Vital Medical Supplies, the Vernacare disposable system has led the way in modern, hygienic and cost effective human waste management for more than five decades. Designed to help protect the nurse and patient from the risk of infection, the comprehensive range of pulp receptacles are manufactured in the UK and carry the B.S.I Kitemark certification to guarantee effective liquid retention and maceration.

## Pulp and Pulp Support Products

 <b>101AA100</b> Traditional Bedpan Liners 100s	 <b>111AA200</b> Medium Bowl	 <b>112AA100</b> Large Bowl	 <b>105AA260</b> Kidney Dish	 <b>116AA200</b> Multi Cup
 <b>106AA300</b> CS Receiver	 <b>102AA100</b> Male Urinal Traditional	 <b>103AA120</b> Male Urine Bottle	 <b>109AA100</b> Midi Slipper Pan	 <b>114AA100</b> Maxi Slipper Pan
 <b>110AA100</b> Measuring Jug	 <b>118AA100</b> Washbowl 4L	 <b>121AA075</b> Bed Pan	 <b>122AA050</b> Disposable Bed Pan Support	 <b>108AA032</b> Vernafem Urinal
 <b>119AA060</b> Vernatray	 <b>104AA200</b> General Purpose/Vomit Bowl	 <b>113BA100</b> Mid Stream Sampler	 Plastic Assorted Pan Supports	

## Pan Room Essentials

### Code

### Description

Great support and storage solutions. Our comprehensive range of racks, holders and supports play a key role in helping you maintain a tidy efficiently organised ward and pan room environment.

### VERNAGEL



450MA100	VERNAGEL POWDER 7g 100's
456MA400	VERNAGEL POWDER 4kg
456MA475	VERNAGEL POWDER 475g

### PAN ROOM ESSENTIALS



501ZS001	BEDPAN SUPPORT RACK - STEEL
503ZS001	BEDPAN SUPPORT HOLDER - FLAT PLATE
504ZS001	MALE URINAL BOTTLE HOLDER -BEDSIDE HANGER
BSRE632	BEDPAN SUPPORT RACK WH 400mm
MUH1697	URINAL HOLDER TRADITIONAL
UDE635	UNIVERSAL DISPENSER - POWDER COATED
505ZS001	GP BOWL/ SMALL BOWL/ JUG DISPENSER RACK
506ZS001	COMMUNE ATTACHMENT
508ZS001	MALE URINAL -TRAD. CARTON HOLDER
509ZS001	MALE URINAL CARTON HOLDER
510ZS001	WASHBOWL RACK

# Do your share for the environment

## Replace Plastics

The Vernacare Single-Use products are produced using water and recycled over-issued newsprint and the wide range of disposable single use containers are designed to replace re-usable plastic or stainless steel alternatives. Ensuring a clean bedpan, bowl or container every time, proven to help reduce the risk of cross infection and more importantly environmentally friendly.



Vernacare manufacture over **170m** Single-use containers each year

The single-use system is used in over **50** countries worldwide



**6 weeks:**

Vernacare's pulp containers have a similar composition to toilet paper and can fully decompose in just over a month



**450+ years:**

Plastic toiletting products have a similar decomposition timeline to plastic bottles making them a significant environmental threat



1. Hospitalhealth.com.au. (2019). Clinical Waste - A Major Environmental Burden.

## Detergent Proof Wash Bowl

Maximum capacity 4 litres




Increased patient safety

Holds warm water and detergent for up to 2 hours

Fill level 2 litres

Case Quantity: 100

Vital Code: 118AA100



## Kidney Dish

General purpose receptacle

Ideal for holding a variety of medical items including dressings, wipes or surgical instruments

Case Quantity: 260

Vital Code: 105AA260



## Quality Guaranteed

Vernacare's medical pulp products carry the BSI Kitemark and are manufactured in line with BS ENISO 9001: 2000 and PAS29:1999, using 100% over-issued clean newspaper. As the biggest medical pulp moulder in the world, producing around 150 million units per year. Vernacare have invested heavily in automation and quality control measures to ensure that the product quality is of a consistently high standard.





## Bedpan

Regular hospital Green

# \$26.50

DOSBP  
GST Free



## Urinals

Male-natural

# \$13.95

DB040  
GST Free



## Sterile and Non-Sterile Single Use Bowls



		Each	10+
MGTMG1010	120ml Iodine bowl graduated clear	\$0.54	<b>\$0.43</b>
MGT06790	150ml Bowl Sterile Blue	\$7.65	<b>\$5.45</b>
		Each	50+
MGTDP1559	150ml Bowl Non-Sterile Blue	\$1.65	<b>\$0.83</b>



## Bowls - Stainless Steel

ParagonCare

Iodine Bowls		Each	3+
LRM042	65x33mm	\$5.20	<b>\$4.50</b>
SL142	77x43mm	\$8.20	<b>\$6.75</b>
SL144	100x50mm	\$9.20	<b>\$8.10</b>
LRM048	120x60mm	\$11.50	<b>\$10.25</b>
Bowls			
LRM082	140x60mm	\$11.30	<b>\$9.95</b>
LRM084	180x80mm	\$23	<b>\$19.95</b>
LRM128	3.7 Litre	\$37	<b>\$29</b>



## Bowls

Iodine Bowls		Each	10+
DOB08	80mm	\$5.20	<b>\$3.95</b>
DOB11GR	110mm	\$5.40	<b>\$4.75</b>
Bowls			
DOB14	140mm	\$9.30	<b>\$7.95</b>
DOB18	185mm	\$4.90	<b>\$4.30</b>
DOB24	240mm	\$16.20	<b>\$14.95</b>



ParagonCare



## Kidney Dishes - Stainless Steel

		Each	3+
LRM018	172x93x34mm	\$13.90	<b>\$12.50</b>
SY409P20	200x95x38mm	\$18.50	<b>\$16.75</b>
SY410P25	250x110x45mm	\$27	<b>\$22</b>
LRM034	300x145x57mm	\$34	<b>\$27</b>

## Kidney Dishes

		Each	10+
DOK16	160mm	\$6.95	<b>\$4.95</b>
DOK30	320mm	\$14.95	<b>\$13.95</b>
Yellow			
DOK21Y	220mm	\$8.95	<b>\$7.95</b>
DOK25YEE	255mm	\$8.35	<b>\$6.50</b>



## Sterile and Non-Sterile Single Use Kidney Dishes

		Each	10+
GM05-800	700ml Kidney Dish Non-Sterile Clear	\$0.94	<b>\$0.66</b>
GM06-887EE	700ml Kidney Dish Sterile Clear	\$2.75	<b>\$1.93</b>
MGT05801	700ml Kidney Dish Non-Sterile Yellow	\$0.87	<b>\$0.65</b>
MGT06883	700ml Kidney Dish Sterile Yellow	\$3.36	<b>\$2.75</b>



10% GST will apply to all items unless they are marked GST Free



## Autoplas Plastic Autoclavable Instrument Trays



Compartment Trays		Each	10+
DB033	3C 160mmx140mmx20mm	\$16	<b>\$14</b>
DOT3C	3C 205mmx165mmx30mm	\$15	<b>\$12</b>



Flat Instrument Tray		Each	10+
DOT22	220mmx140mmx20mm	\$10.20	<b>\$8.95</b>
DOT24	245mmx195mmx40mm	\$13.95	<b>\$12.95</b>



Perforated Trays		Each	10+
DOTP1815	180mmx150mmx30mm	\$9.95	<b>\$8.95</b>
DOTP2715	270mmx150mmx30mm	\$13.50	<b>\$11.95</b>
DOTP2730	270mmx300mmx60mm	\$29.95	<b>\$27.95</b>
DOTP2740	270mmx400mmx60mm	\$49.95	<b>\$47.95</b>

## Instrument Trays - Stainless Steel

ParagonCare



			Each	3+
LRM002	w/lid	200x80x38mm	\$37	<b>\$32</b>
SL231	w/lid	200x130x50mm	\$39	<b>\$36</b>
LRM028	w/lid	300x200x65mm	\$73	<b>\$65</b>

## Non-Sterile Injection Tray



		Each	10+
GM05-111B		\$0.32	<b>\$0.22</b>

## Dressing Jars

ParagonCare



STAINLESS STEEL

		Each	3+
LRM029	100x125mm	\$34	<b>\$28</b>
LRM026	127x165mm	\$47	<b>\$39</b>

## JUGS Stainless Steel

ParagonCare



		Each	3+
LRM033	500CC	\$46	<b>\$39</b>
LRM030	1000CC	\$49	<b>\$43</b>
LRM025	2000CC	\$55	<b>\$45</b>

## Forceps Jar Stainless Steel

ParagonCare



		Each	3+
LRM031	55x114mm	\$37	<b>\$33.50</b>



## Customer Value Promise

Our experienced team is committed to bringing you the best value.

Vital will always offer the best deal.

If you think you've found a better one, why not give us a call and we'll ensure you get the best solution for your practice