




Green Series™ 300 General Exam Light

Utilising a single white LED, the GS 300 Exam Light provides a whiter, brighter light than halogen and consumes less energy.

Features

- Perfect combination of intense light output and cool operation in a compact design
- Affordable price and rugged performance
- Touchless On/Off
- Multiple mounting options for floor space maximisation
- Lower cost of ownership and energy consumption due to elimination of bulb replacement
- Easy assembly
- 5 year limited warranty to include the LED light source and power supply only



		Each	3+
WL44416	GS 300 Exam Light with wall/table mount	\$599	\$524
WL44456	GS 300 Exam Light with mobile stand	\$803	\$720




Green Series™ 600 Minor Procedure Light

Featuring three LEDs and a wider head, the GS 600 is designed specifically for demanding viewing conditions.

Features

- Three white LEDs for greater intensity
- Touchless power controls and adjustable light output
- Disposable sheath available to help reduce the risk of cross-contamination
- Rugged, ergonomic design allows for improved performance and long-lasting durability
- Low cost of ownership due to elimination of bulb replacement and lower energy consumption
- Easy to assemble and clean
- 5 year limited warranty to include the LED light source and power supply only



		Each	3+
WL44616	GS 600 Procedure Light with wall/table mount	\$1,255	\$986
WL44606	GS 600 Procedure Light with mobile stand	\$1,523	\$1,245



Green Series™ 900 Procedure Light

This procedure light, featuring six LEDs, combines outstanding durability, extraordinary maneuverability, unparalleled spot quality and multiple mounting options.

Features

- Six white LEDs for greater intensity
- Touchless On/Off with dimming option
- Multiple mounting options for floor space maximisation
- Disposable sheath available to help reduce the risk of cross-contamination
- Rugged, ergonomic design allows for improved performance and long-lasting durability
- Lower cost of ownership due to elimination of bulb replacement and lower energy consumption
- 5 year limited warranty to include the LED light source and power supply only



		Each
WL44900-C	GS 900 Procedure Light with ceiling mount	\$3,044
WL44900-W	GS 900 Procedure Light with wall mount	
WL44906	GS 900 Procedure Light with mobile stand	



Green Series™ Exam Light IV

Perfect for the OB/GYN office or for specialties that require very focused light in specific areas. The GS Exam Light IV provides you with 3x the typical light output and features a highly flexible light pipe that is twice as long as conventional pipes. With a compact, unobtrusive head, it ensures pinpoint accuracy and provides true tissue colour rendition during an exam.

Features

- One LED - 3x the typical light output; pinpoint accuracy; true tissue colour rendition
- Highly flexible light pipe for easy maneuverability
- Lightweight pipe is twice as long as conventional pipes
- Touchless On/Off
- Multiple mounting options for floor space maximisation
- Disposable sheath available to help reduce the risk of cross-contamination
- Lower cost of ownership due to elimination of bulb replacement and lower energy consumption
- Easy to assemble
- 5 year limited warranty to include the LED light source and power supply only



		Each	3+
WL48816		\$1,360	\$1,133





HEINE EL10 LED® Examination Light

The EL10 LED examination light distinguishes itself through an absolutely homogeneous field of light, ergonomic design for intuitive and fast adjustment of brightness and spot size. Combine this with the compact illumination head and you have an almost coaxial illumination perfect for difficult situations.

- Bright and homogeneous illumination
- Coiled goose neck for perfect positioning
- Up to 45,000 LUX at 30cm
- Solid construction and compact design
- Spot size diameter between 8.5 and 14cm
- Easy one-hand operation
- Colour temperature 3000k
- Hygienic design for efficient cleaning and disinfection
- Colour Rendering Index >90
- 3 mounting options: Wheeled stand, wall or table mount
- 5 year warranty



		Each	3+
HEJ0827001	EL10 LED Examination light with Wall Mount	\$2,600	\$2,356
HEJ0827002	EL10 LED Examination light with Clamp Mount for Table Top Mounting	\$2,774	\$2,534
HEJ0827003	EL10 LED Examination light with Wheeled Stand & Metal Base	\$3,021	\$2,777



LED HQ

LED NOW IN HEINE QUALITY.



ri-magic® Examination LED Light

The flexible examination light ri-magic® LED illuminates the examination and operation area in the best possible way.

- The most modern LED technology, similar to daylight and high-performance; ideal for examinations and smaller procedures
- Lifetime of the lamp up to 50,000 hours, replacements no longer necessary
- Significantly less heating and energy loss
- The focusing device at the front of the light head allows exact adjustment of the size of the illuminated examination area
- Highly flexible, long special arm can be bent in all directions, the light remains fixed in the chosen position
- Simple to use handle with rheostat for infinitely variable regulation of light intensity
- Wide-range power adaptor for 230V/120V
- 2 Year warranty



RIE1010235	Wall mount	\$993
RIE1010236	Mobile mount	\$1,045
RIE11343	Mobile stand and adaptor	\$295



HEINE EL3 LED Examination Light

The light from LEDHQ is homogeneous and evenly distributed across the entire field of view, resulting in even illumination up to the outer areas without disturbing reflections, shadows, bright or dark spots.

- Maintenance-free and easy to clean
- Optimal positioning with lightweight head and solid goose-neck
- Up to 42,000 LUX at 30cm
- Bright and homogeneous illumination field between 11-13cm
- Colour temperature 5000k
- Colour Rendering Index >80
- 3 mounting options: Wheeled stand, wall or table mount



		Each	3+
HEJ0827011	EL 3 LED Examination Light with wall mount	\$1,203	\$1,109
HEJ0827013	EL 3 LED Examination Light with clamp for table-top mounting	\$1,398	\$1,281
HEJ0827014	EL 3 LED Examination Light with wheeled stand, metal-base	\$1,695	\$1,619



Halogen Examination Light

Cost-effective solution for surgery examination lighting

- 35 watt globe, 12V (approx 4000hrs)
- 70 cm long gooseneck
- Easy grip handle with on/off switch
- Additional glass filter reduces heat

	Each	3+
MGL308	\$279	\$269



Brackets to be purchased separately



LED Examination Light

- 700mm Gooseneck
- Sure grip handle with on/off switch
- Cool white light, 5500°K colour,
- Concise beam angle of 30°
- 3 x 1 watt high power LED giving 16,000 lux at 30cm
- Low energy consumption, only 5 watt
- Low heat output
- Wall or mobile mounted
- 40,000 hour bulb life
- 2 year warranty



	Each	3+
MGL308LED	\$489	\$479





Ordisi Iglux Examination Lights

Ultimate day procedure light

- Compact design with smooth edges that make it easy to clean
- 12 month warranty
- Light intensity adjustment through the touch panel controls at front of lamp
- Handle is removable and safe to sterilise
- 20,000 - 65,000 lux adjustable brightness

		Each	3+
ODIGLUXM	Light Ordisi Iglux Mobile	\$3,295	\$3,185
ODIGLUXC	Light Ordisi Iglux Ceiling Standard 80cm Pole	\$3,454	\$3,179
ODIGLUXCSH	Light Ordisi Iglux Ceiling Short 34.5cm Pole	\$3,839	\$3,454



Ordisi Mobile Examination Light

Cool light source

- Expected LED lifetime 50,000 hours
- Very power power consumption (10w)
- Strong, white light
- 42,000 lux/0.5m
- 12,000 lux/1m
- lightweight mobile base
- 12 month warranty

	Each	2+
ODLEDMEXAM	\$644	\$626



Ordisi Wall Examination Light

Includes Swivel Bracket

- Cool light source
- Expected LED lifetime 50,000 hours
- Very power power consumption (10w)
- Strong, white light
- 42,000 lux/0.5m
- 12,000 lux/1m
- Goose style neck for easy maneuverability
- 12 month warranty

	Each	2+
ODFLH2731W	\$644	\$626





Green Series Mobile Stand Examination Light IV

Examination Light IV

\$409

WL48950

Examination Light Tall

\$409

WL48960



Freight charges will apply

WelchAlllyn

Mobile Base for LHH Examination Light

Freight charges will apply



\$345

LLTS

LUXO

Mobile Base

- Stable 5-star solid base
- Black 50mm twin wheel castors
- Chrome plated stem with key wheel
- Height 850mm stem
- Diameter 12.5mm

Freight charges will apply



\$169

VIDCD

Extended Wall Mount for GS300/600 Examination Light



\$293

WL44215

WelchAlllyn

Ceiling Mount Bracket for LHH 60cm



\$355

LU2529

LUXO

Wall Bracket and Base



\$75

LUX899

LUXO

Maggylamp Wall Bracket



\$65

MGL307

MAGGYLAMP
medical

Maggylamp Mobile Base Bracket



\$52

MGL309

MAGGYLAMP
medical

ri-magic Mobile Stand + Adaptor

Riester

\$310

RIE11343



Disposable Sheaths for Green Series™ 600/900



	Each	3+
WL52630	\$230	\$199

WelchAlllyn

Disposable Sheaths for Green Series™ Exam Light IV



\$40.95

WL52640
Box 25

WelchAlllyn

35 Watt Halogen Bulb for Maggylamp



\$24

MGL306

MAGGYLAMP
medical



MAGGY-VUE Magnifying Lamp

- 4 Diopter circular lens, Ø 120mm
- 2 x magnification
- Shadow-free diffused lighting
- 22 Watt fluorescent lamp (included)
- Spring balanced joints
- Optional clip-on lens for increased magnification (3x)



MGL150WM	Wall Model	\$425
MGL150MM	Mobile Model	\$425
MGL151	Add-on Lens	\$155

MAGGY LAMP
medical

LAMP EXAMINER LFM30-101

- Magnification and cool lamp
- Round 4 diopter lens
- Internal spring-balance arm with parallel movement
- Magnification ratio: x 100% actual size



\$225

LULKMAGW

LUXO

HEINE mini-c® Clip Lamp



Compact, long-life pocket diagnostic light
Sturdy, long-life pocket diagnostic light in a compact, modern design with chrome head. Bright, concentrated light thanks to novel XHL Xenon Halogen lens-lamp. Mini-c clip lamp only fits a mini-c handle. Not compatible with mini 3000 instruments.

- Compact, oval, attractive light
- Fits ergonomically in the hand
- High-quality handle: Shockproof, sturdy, non-slip
- Switches off automatically when put back in pocket
- Replaceable batteries. Size AAA



\$110.²⁰

HED0173109

Ear Light for mini-c Clip Lamp

6 pieces in blister pack

HED0073105	Pkt/6	\$134
------------	-------	--------------



MAG-LITE®

High Quality Examination Torches

Mini Maglite

- High-intensity light beam
- 1/2 turn, twist focus, spot-to-flood
- Spare lamp safely secured in the tailcap
- High-strength aluminium alloy case
- O-ring sealed for water resistance



AAA

\$29.⁹⁵

SH85200

AA

\$39.⁹⁵

SH85700

The Mini Maglite® 2-Cell AAA flashlight is a smaller and lighter version for ease of use and storage.

MagLite 3 Cell

A very popular D cell Mag-Lite. Featuring a long handle, extending its length to 31cm.



\$49.⁹⁵

SH86230




Diagnostic Penlight

- Modern satin finish aluminium body – deluxe model
- Pocket clip doubles as on/off switch
- Useful 5cm rule and dilation guide
- 2 x AAA batteries included

\$24.⁵⁰

SAS233



ParagonCare

WOODS LIGHT Ultra Violet

For the detection of various skin disorders

- Convenient hand size
- Rugged plastic construction
- Plastic lens over tube to protect patients' eyes
- Includes torch

\$25

LRUVL

Tube for Woods Light
LRUVT \$9.50



Promar Torches

ParagonCare

- Sealed unit with 2xAAA size batteries included - non rechargeable, non replaceable batteries
- Concentric beam bulb with pocket clip/switch
- Colour - White



	Each	+6
LRDDT	\$3.95	\$3.50

Handheld LED Illuminated Magnifier

- 85 x 50mm lens
- Crystal clear acrylic lens with
- 2.0 X mag and 4.0 X mag bubble
- Complete with batteries and pouch
- Ergonomic shaped handle fits comfortably in the hand

\$29

MGL485



MAGGYLAMP UV Woods 4x4 WATT

- Fitted with 4 x 4 Watt UV light
- UV Wavelength Lens: black light blue, 365nm biconvex glass, rectangular 100 x 50mm
- Magnification: 2.25 x 5 diopre

\$549

RB215



Mag-eyes The hands-free, head mounted magnifier

Mag-Eyes with 2 Lenses (1.6/2.0 X Magnification)

\$79

MGL400

Mag-Eyes with 2 Lenses (2.25/2.75 X Magnification)

\$95

MGL410



Working distances (Lens to Subject)	
1.6 X (30-35cm)	2.0 X (25-30cm)
2.25 X (18-25cm)	2.75 X (10-18cm)

Guthrie Loupe

ParagonCare

- Adjustable head band padded for user comfort
- Plastic, precision lens
- Retractable visor
- 2.7x Magnification

\$95

LRM2001





HEINE Loupes & Headlights

Perfect view for highest accuracy during all procedures.

❖ Flexible Flip-Up Design

The i-view bar allows for a steep angle of declination to be achieved and is as individual as your daily work.

❖ Large Field of View

2.5x High-resolution binocular loupes gives a crisp and clear view (up to 180mm) and a long depth of field up to 260mm.

❖ Lightweight Design

Combination of the Carbon-Kynetium frame, light weight i-view bar, adjustable head strap, frame arms and nose pad allows for the perfect fit.

❖ Perfect Illumination

The Loupe Light 2 in LED HQ offers optimum brightness at 55k Lux, edge-to-edge homogeneity, true colour rendering >90 CRI, ideal temperature management and an unlimited working life.



HEINE 2.5x HR Loupes & LoupeLight 2

Most trusted by ENT specialists: The combination of light and loupe delivers the ultimate in performance as high quality illumination with LED HQ maximises the image quality of the loupe.

For intensive use, the LoupeLight2 now features a replaceable connecting cord which connects the device to the mPack mini – the new and compact power source for the LoupeLight2.

- Compact and lightweight
- Flexible Flip-up design
- Large Field of View
- Perfect Illumination

Includes:

- 2.5x HR Loupes
- LoupeLight2
- mPack mini
- mPack mini belt clip
- E4 USB Plug-In Transformer
- 2 x Protective Lenses
- 2 x Sterilisable Swivel Levers
- Retaining Cord
- Cleaning Fluid
- Cleaning Cloth
- Carry Case



		Each
HEC0832450	HR 2.5x/340mm, LL2, I-View A.S-Frame	\$5,447
HEC0832451	HR 2.5x/420mm, LL2, I-View A.S-Frame	
HEC0832452	HR 2.5x/520mm, LL2, I-View A.S-Frame	
HEC0832265	LoupeLight 2 F HR/HRP + mPack Mini	\$2,863
HEX0099007	mPack Mini Belt Clip	\$23.40
HEC0032240	Connecting Cord SC1 (1,5m/Ø 2,4mm)	\$231
HEC0032243	Connecting Cord SC2 1.5m (1,5m/Ø 3,2mm)	\$231



HEINE 2.5x HR Loupes on S-Frame

HEINE S-Frame®

- Innovative spectacle frame for HR Binocular Loupe
- Appealing and functional design
- Can be customised to fit any shape of head and face

HEINE HR® Binocular Loupe 2.5X

- The brightest image – Heine's exclusive coating technology reduces reflection to 0.3%, providing the brightest image of any comparable loupe
- Crispest image with colour correction
- Super-wide, super-deep field of view
- Dust and water resistant

		Each
HEC0032355	340mm working distance	\$2,575
HEC0032356	420mm working distance	
HEC0032567	520mm working distance	

HR Binocular Loupe Set with i-View Optics, S-Frame with i-View loupe mount for HR 2.5 x loupes, detachable protective lenses, 2 sterilisable swivel levers, retaining cord, cleaning fluid, carrying case.





HEINE ML4 LED Headlight

HEINE QUALITY
MADE IN GERMANY

Ultra Bright Coaxial Illumination – Now with LED ^{HQ}

- **LED in HQ** - the new standard in LED illumination defining optimal light intensity, homogeneity and colour rendering for each individual examination situation.
- **Adjustable illumination spot size** - 30 mm to 80 mm spot size range (at 420 mm working distance) to fit any examination situation.
- Bright and homogenous illumination. Absolutely bright light spot that is uniform from edge-to-edge for the perfect illumination in all examination situations.
- **Stepless light intensity control** - Optimal brightness setting prevents reflexes.
- **Flexibility**. Flexible power source options: 100% mobility with a choice of cable-free mPack UNPLUGGED head-worn battery or mPack belt-worn battery pack.
- **Coaxial Illumination** - Coaxial design ensures a completely shadow-free image and allows for excellent illumination of difficult to see areas.
- **Comfortable and secure fit** - The ML4 Professional headlight has multiple adjustment points and soft padding for comfort and stability – even during long examinations.
- **Optional Filters** - Polarisation filter P 2 improves contrast and yellow filter reduces blue light.



		Each
HEJ0831413	ML4 Headlight with mPack	\$3,322
HEJ0831415	ML4 LED HeadLight with Headband Professional L and Power Sources	\$3,452



HEINE ML4 LED Headlight with Binocular Loupes

Combining HEINE high performance binocular loupes with the ML4 headlight to provide ideal magnification and illumination for surgical and examination applications.



		Each
HEJ0831530	Headlight ML4 Unplugged with 2.5x 340mm Loupes	\$5,637
HEJ0831531	Headlight ML4 Unplugged with 2.5x 420mm Loupes	
HEJ0831532	Headlight ML4 Unplugged with 2.5x 520mm Loupes	



HEINE mPack

- **Portable power pack** - Does not tie the user to desk, wall or static power supplies.
- **Li-ion Technology** - Triple capacity compared with current NiCd systems. No memory effect, autonomy for 8.5 hours of continuous ON time at full power.
- **Fast Charge** - Charges to full capacity in 2 hours.
- **Charge Status Display** - Always lets you know how much operating charge you have left.
- **Stepless light intensity control** - Optimum brightness setting prevents reflexes.
- **Flexible Charging Options** - EN 50 Desk and wall charger, or plug-in transformer.
- **Intelligent** - Can be used as mains power supply even while in “charge” mode.
- **100 – 240 V Compatible Transformer** - incl. 4 plugs for world wide compatibility.

