



TENS HV-F128 - High Frequency TENS Therapy

Relieve pain and muscular tension at the touch of a button

The HV-F128 is the premium TENS device from Omron. With 12 automatic therapy modes the Omron HV-F128 can easily be set to offer either pain relief or muscle stimulation at the touch of a button. The LCD screen displays the time remaining of the therapy and the frequency of the electronic stimulation being delivered.

Features:

- Management of muscular injuries
- 12 Automatic therapy modes
- LCD panel shows frequency and time remaining of therapy
- High frequency stimulation to relieve pain quickly
- Single channel device



\$225

DAHVF128
GST Free

HVF-128 Long Life Pads (pair)

\$50

DALLPAD

OMRON

ProTens Machine

- Dual Channel
- Hard Carry Case
- Battery Powered
- Adjustable pulse width, pulse rate plus intensity
- Timer of 15, 30 and 60 minutes
- Will stop output at pre-set time

\$49

ACPROTENS



allcare

AC35050	ProTens Electrodes 5 x 5cm	Pkt/4 \$10
---------	----------------------------	----------------------

Allcare Brand Exercise Band Rolls

- Steady 30% increase from colour to colour
- To use, simply cut a length from the roll and tie a knot to form a loop. The knot can be undone easily to change the length or size of the loop. Allband Exercise Band may also be used as a single length
- Available in 6 colour coded resistances, providing a wide range of treatments
- 99.9% latex allergen free
- Made under the international quality standard ISO 9002



ACBAND25YE	Yellow Light 25m	\$37.00
ACBAND25RD	Red Medium 25m	\$45.00
ACBAND25GN	Green Firm 25m	\$48.00
ACBAND25BU	Blue Extra Firm 25m	\$52.00
ACBAND25BK	Black X-Extra Firm 25m	\$63.00
ACBAND25SL	Silver XX-Extra Firm 25m	\$75.00

allcare

Premium Rigid Sports Tape

The highest grade tape available, combines excellent adhesion and maximum tensile strength to provide ultimate support to joints under stress. Cerated edge for easy hand tearing.



ACTAPE12	12.5mm X 13.7m	\$4.50
ACTAPE25	25mm X 13.7m	\$7.20
ACTAPE38	38mm X 13.7m	\$8.50
ACTAPE50	50mm X 13.7m	\$13.00

allcare

K Tape

AllCare K Tape is a tape that has 140% elasticity which matches the same characteristics of the human skin. It gives non-restricted support that permits muscles to perform at full range of motion.

ACKT31BK	Black 50mm x 31m	\$55.00
ACKT31BU	Blue 50mm x 31m	
ACKT31FL	Flesh 50mm x 31m	
ACKT31PK	Pink 50mm x 31m	
ACKT31RD	Red 50mm x 31m	
ACKT5BK	Black 50mm x 5m	\$11
ACKT5BU	Blue 50mm x 5m	
ACKT5FL	Flesh 50mm x 5m	
ACKT5PK	Pink 50mm x 5m	
ACKT5RD	Red 50mm x 5m	



allcare



Ortho Walker

New contoured strut design allows the walker frame to conform to your patients anatomy ensuring a better fit. Unique shock absorbing sole reduces impact from the heel strike to improve comfort during ambulation. Rocker bottom is specially engineered to be low & wide to help promote a natural gait & increase stability.



WK001R30	Regular 17" small	\$47.00
WK001R40	Regular 17" medium	
WK001R50	Regular 17" large	
WK001R60	Regular 17" XLarge	
WK001S30	Short 11" small	
WK001S40	Short 11" medium	
WK001S50	Short 11" large	
WK001S60	Short 11" XLarge	



Ortho Air Walker

The AllCare Air Walker designed for the treatment of stable fractures and ankle sprains. Rocker bottom and low profile heels facilitate more nature gait pattern. Adjustable air cells increase stability while decreasing pain and swelling.



WK007R30	Regular 17" small	\$65.00
WK007R40	Regular 17" medium	
WK007R50	Regular 17" large	
WK007R60	Regular 17" XLarge	
WK007S30	Short 11" small	
WK007S40	Short 11" medium	
WK007S50	Short 11" large	
WK007S60	Short 11" XLarge	



Sam Splint

The Sam Splint Original is an extremely pliable multi-purpose splint and is built from a thin core of aluminium alloy, sandwiched between two layers of closed cell foam. Bent into various shapes and curves, it becomes extremely strong and supportive for any fractured or injured limb.

Available in a range of sizes, it is lightweight, X-ray penetrable, reusable, waterproof and very portable – almost any bone in the body can be splinted.

FNSVPSS36	Sam Splint 36" roll	90L x 10W cm	\$19.50
FNVPS36XL	Sam Splint 36" XL	90L x 14W cm	\$29.95
FNSVPSS18	Sam Splint 18" Junior	45L x 10W cm	\$21.95
FNSVPSS9	Sam Splint 9"	23L x 10W cm	\$12.95
FNSVPSSF	Sam Splint Finger		\$3.50

The Sam Splint stores small yet fits all! It is a simple splint with many features.



Actimove® Mitella Comfort

Arm Sling

After trauma injuries at shoulder, elbow, wrist or arm

- Comfortable and secure support
- Thumb loop prevents sling migration
- Soft neck pad provides extra comfort
- Adjustable neck strap



SN7281961	Small	\$58.12
SN7281962	Medium	
SN7281963	Large	



Actimove® Sling

Collar and Cuff Support

- Versatile concept of strap and Y-tab fasteners allows for various application techniques
- Non-stretch foam strap holds the arm in the desired position
- Blue colour for improved aesthetics and clean look
- Economical dispenser box for multiple applications. Contains 12 meter roll for cutting to size
- Includes pack of 14 Y-tab fasteners

SN72859-18	12m Roll	Box/2	\$50.88
------------	----------	-------	----------------



Actimove® Mitella Universal Arm Sling for Kids

- Actimove® Mitella for Kids
- Colorful arm support
- Actimove® Mitella for Kids is a paediatric arm support with a fun and colourful design.

BE7456113	Each	\$64.58
-----------	------	----------------





Philly Patriot Collar

The Philly Patriot Collar is moulded from high density polyethylene with foam padding and a patented chin support. This collar restricts cervical spine flexion, extension and rotation to promote patient recovery. The one-piece design enables efficient storage where space is limited. The multi-height adjustment fits almost all adults (four positions).

The one-piece design can easily be used; simply measure patient, adjust and lock selected height and apply. For either a supine or seated position, application is quick and easy.

Key Features

- Adjustable hook and loop fasteners ensure accurate fit
- Anterior and posterior flex tabs
- Oversized trachea opening
- Easily customised for individual requirements
- Latex-free

Adult
\$35
FNCCPCPA

GST Free



STAX Finger Splints



MUS01	Size 1	\$4.60
MUS02	Size 2	
MUS03	Size 3	
MUS04	Size 4	
MUS05	Size 5	
MUS06	Size 6	
MUS07	Size 7	

Collar & Cuff Arm Sling

Contents

- 5mx6m of sling material
- 8 plastic fasteners

\$12.95
RMCC06-B



Finger Splints



DD4813	13mmx45cm	\$4.95
DD4819	19mmx45cm	\$4.95
DD4825	25mmx45mm	\$4.95
DM4832	32x480mm	\$4.50

Actimove® Genu Tri-Panel Knee Immobiliser

Actimove® Genu Tri-Panel provides effective immobilisation of the knee and helps to reduce pain following trauma or surgery. The 'tri-panel' design allows splint adjustments for a broad spectrum of patient leg circumferences and the hook and loop sensitive closure system permits easy application and removal. Removable and adjustable aluminium stays (2 posterior and 2 on each side panel) allow good anatomical fit.



SN7281914	50cm Small	\$91.62
SN7281915	50cm Large	\$99.02
SN7281916	60cm Small	\$106
SN7281917	60cm Large	\$113.93



BSN medical

Actimove® Umerus Comfort Shoulder Immobiliser

- Waist strap provides secure, anatomical fit
- Optimal sling position maintained by thumb loop
- Customised straps for accurate, adjustable fitting
- Easy injury inspection



SN7281939	Small 34-38cm	\$70.25
SN7281940	Medium 37-43cm	
SN7281941	Large 42-48cm	
SN7281942	X-Large 48+cm	

BSN medical

Order Online

- Save time
- Compare products
- Review Order history

www.vitalmedicalonline.com.au





Instant Ice Packs



B013	Instant Cold Compress 15.3 x 22.9cm	\$2.20
------	-------------------------------------	---------------

Reusable HOT/COLD Gel Pack

- Can be reused for either hot or cold
- Remains soft and pliable after freezing
- Cleans easily with soap and water
- 100mm x 150mm



Product A

- To heat, microwave or place in hot water, please follow
- the instructions on the back of the compress
- To obtain cold therapy, freeze for at least 1 hour
- This gel pack can reach a high of 60°C and a low of -12°C
- Maintains therapy for approximately 30 minute



Product B



Product C

		Each	10+
LR3057	Product A 10 x15cm	\$2.95	\$2.50
LR3510	Product B 12.5 x 25cm	\$3.95	\$3.75
SM10130312	Product C 25cm x 17cm	\$5.95	\$4.20

ParagonCare

BodyWARMER

After opening, the body warmer gradually heats up. Once temperature is reached it will last for around 12-15 hours, providing warmth for the user. It comes with adhesive to apply to body or clothing. It provides long lasting warmth to assist users in cold weather.

135mm x 100mm

\$2.50

LRBODYW



TheraBeads® Professional Packs

DMI® TheraBeads® Cervical Professional Pack

- Designed for use in professional practices, hospitals and clinics
- Wraps around neck
- Color-coded zones indicate temperature range
- Durable, long-life fabric emits heat through 1 side only maximizing moist heat therapy
- Latex Free
- 9" x 24"

\$95

BR7105




DMI® TheraBeads® Standard Professional Pack

- Designed for use in professional practices, hospitals clinics and physiotherapy
- Comes with extra long life durable fabric that's used with a special skin temperature range
- Monitor uses bold, colour coded zones to indicate at a glance when pack is in therapeutic temperature range
- Helps ensure the patient receives prescribed treatment accurately and consistently
- 12" x 12"

\$95

BR7103





Standard Wheelchair

- 110kg user weight capacity
- Seat 45cm W x 40cm D
- Frame: Single cross brace, powder coating
- Front wheel: 8"/6", rubber tires
- Rear wheel: 24" spoke wheel, pneumatic tires
- Armrest: Flip back with desk length
- Footrest: Swing away footrest
- Wheelchair is foldable



\$302

PM18WCS

Patient Mover

- 110kg user weight capacity
- Seat 46cm W x 41cm D
- Frame: Single cross brace, powder coating
- Front wheel: 8", PVC tires
- Rear wheel: 12" PU tires
- Armrest: Flip back with desk length
- Footrest: Swing away footrest



\$316

PM18WCM

Bariatric Wheelchair

- 160kg user weight capacity
- Seat 55cm W x 40cm D
- Frame: Single cross brace, powder coating
- Front wheel: 8"x2", PU tires
- Rear wheel: 24" spoke wheel, PU tires
- Armrest: Flip back with desk length
- Footrest: Swing away footrest
- Seat and backrest: Nylon with sponge



\$375

PM22WCB

Combi Transit Wheelchair

- Complies to Australian Standards
- Swing-away removable leg rests
- Solid 200mm front castors
- Removable arm rests
- 46cm seat width
- 100kg capacity



\$325

LRCMWCT GST Free

ParagonCare

Folding Wheelchair

- Complies to Australian Standards
- Removable arms
- Heel Loops
- Swing away detachable footrests
- Breathable padded upholstery
- Rear backrest storage pocket
- 46cm seat width
- 19.5kg gross weight
- 120kg SWL
- Solid Tyres



Top Seller!

\$335

AXWC604 Freight charges will apply

ParagonCare

Crutches Aluminium Elbow Crutch Double Adjustable



Code	Type	Height to Handgrip	Cuff to Handgrip	Patient Height	Max user Weight	
463400M	Med	550-800mm	200-275mm	1470-1720mm	125kg	\$65
463400L	Large	680-930mm	200-275mm	1580-1830mm		\$65
463400XL	X-Lrg	850-1110mm	210-285mm	1760-2010mm		\$69

Aluminium Under Arm

				+2
41110020	Youth	3' - 3.7'	\$49	\$45
41110021	Medium	4.7' - 5.3'		
41110022	Adult	5.2' - 5.10'		
41110023	Tall Adult	5.10' - 6.6'		

