

# seca 286

## Wireless ultrasonic measuring station for height and weight with voice guidance

### EMR-integrated

The seca 286 ultrasonic measuring station offers fully automatic measuring and was designed to be user-friendly for hands-free operation.

Voice guidance can be optionally selected, combined with the included poster and sticker on the measuring station illustrates the correct stance and posture. This makes it easy for the patient to use without any assistance from the nursing staff.



### Intelligently precise

- Capacity: 300 kg
- Graduation (g): 50 g
- Measuring range in cm: 60-210 cm
- Graduation [Measuring Rod]: 1 mm
- Power supply: Power adapter
- Measure (W x H x D): 434 x 2270 x 466 mm
- Approval class: III (3)
- Wireless connection to EMR or Direct print to any printer via Seca SE456 dongle
- Optional printer (shown)

2 Year Warranty

**\$3,775**

PN286



seca 467



**seca**  
Precision for health



# seca mBCA 515

## Medical Body Composition Analyser



The seca 515 is the **ONLY** Medical Body Composition Analyser in the world that is medically validated, verified and measured against all the international gold standards.

Package includes Windows-based software package (seca software 115) for detailed patient records and report printing

**Gold standards include:**

- Dexa standard
- MRI standard
- Bodpod (air displacement plethysmography) standard
- D-20 Standard
- NABR Standard

**This unit will measure:**

- Fat Free Mass
- Muscle Mass
- Body Composition Chart
- Bioelectrical impedance vector analysis
- Total Body Water
- Skeletal Muscle Mass
- Visceral fat



**\$18,895**

PN515

**seca**  
Precision for health

**Specifications:**

Dimensions (W x H x D)	976 x 1,251 x 828 mm
Capacity	300 kg
Graduation	50 g < 150 kg > 100 g
Weight	36 kg
Display	8.4" Touch-screen display, can be rotated 360°
Power Supply	Power adapter
Voltage	100 V – 240 V
Power Frequency	50 Hz – 60 Hz
Calibration Class	3
Interfaces	Ethernet
Compatible Printers	Laser printer and inkjet printer via PC software seca analytics 115


### seca 777

#### Digital column scale with eye-level display and height rod

The price-performance champion of measuring and weighing.

The seca 777 digital column scale with eye-level display is the medical-grade solution for measuring and weighing. It features a tilt-proof full metal platform and robust single piece column, that won't give way when pushing the keypad. All patients (especially ones affected by obesity) will feel secure on the platform even when stepping on the edge of the base. The scale is easy to move by staff due to the large rubber coated castors and a wide wheelbase.

- Capacity: 250kg
- Medical-grade, robust scale designed for years of use with high patient volume.
- Built-in measuring rod with 4 - 90" measuring range.
- USB interface for EMR integration solutions or PC connectivity.
- Large rubber coated castors for excellent movability.
- Wide and robust tilt-proof metal platform, no plastic
- 2 Year warranty



**\$1,045**

SEC777



## seca 700

Mechanical Column Scales With Eye-Level Beam

- Capacity: 220kg
- Graduation: 50g
- Reset-to-zero function
- Height rod included
- 2 year warranty

Measuring rod:

- Measuring range: 60-200cm
- Graduation: 1mm
- 2 year warranty

### \$730

PN700



seca

## seca 703

Hospital Grade Scale

When weight, height and BMI have to be obtained, the new seca 703 can do it all for toddlers to obese adults. The seca 703 with optional integrated stadiometer is a complete measuring station that can measure both weight and height in just one time-saving step.

- Capacity: 300kg
- Graduation: 50g < 150kg > 100 g
- Dimensions (WxHxD): 360 x 930 x 520 mm
- Platform (WxHxD): 335 x 80 x 345 mm
- Optional Measuring range: 7-230 cm
- Weight: 17.6 kg
- Power supply: Batteries / power adapter (optional SE400)
- Data transmission: seca 360° wireless technology to PC or seca printer
- Functions: TARE, pre-TARE, mother/child-function,
- HOLD, auto-HOLD, BMI, auto-CLEAR, damping, automatic switch-off, kg/lbs switch-over, SEND/PRINT
- 2 year warranty
- optional height rod (PN220)

### \$1,570

PN703



seca

## seca 284

Wireless measure station for body weight, height and auto BMI functions



Scale:

- Capacity: 300kg
- Graduation: 50g
- Functions: TARE, pre-TARE, mother/child-function, HOLD, auto-HOLD, auto-BMI, CLEAR, auto-CLEAR, auto weighing range switch-over, kg/lbs switch-over, damping, SEND/PRINT
- EMR Compatible
- Optional: printer seca 466

Measuring rod:

- Measuring range: 30-220cm
- Graduation: 1mm
- Functions headpiece: HOLD, user-defined zero setting, automatic switch-off, SEND/PRINT
- 2 year warranty

### \$2,830

PN284



seca

## seca 769

Electronic Column Scale

- Capacity: 200kg
- Graduation: 100g
- Functions: TARE, HOLD, auto-HOLD, BMI, kg/lbs switch-over, automatic switch-off
- Height Measure rod optional (SEC220)
- Batteries included; Optional mains power adapter (code: SE400)
- 2 year warranty

### \$783

PN769



seca



**seca 813**  
Electronic Flat Scale with Very High Capacity

**\$251**

PN813

- Capacity: 200kg
- Graduation: 100g
- Functions: auto-HOLD, step-off, kg/lb switch over, tip-on automatic switch-on, automatic switch-off.
- Optional carry bag Seca 414
- 2 year warranty



**seca sensa 804**

With chrome electrodes for precise measuring

**\$156.45**

PN804

- Capacity: 150kg
- Graduation: 100g
- 24 electrodes. For analysis of body fat and body water (BF/BW)
- Automatic user recognition and automatic BMI calculation. For fast measuring results
- Functions: auto-HOLD, BMI, BF/BW-measurement, user memory, automatic user recognition, tip-on automatic switch-on, automatic switch-off
- 2 year warranty
- Optional carry bag: Seca 415



**seca aura 807**

Advanced technology for enhanced comfort.

**\$124.95**

SEC807

- Capacity: 150kg
- Graduation: 100g
- Functions: auto-HOLD, step-off, tip-on automatic switch-on, automatic switch-off
- Stable construction for users up to 150kg
- Roughened glass surface for maximum slip resistance
- Large LCD figures for easy reading
- Optional carry bag Seca 415
- 2 year warranty



**seca clara 803**

Clarity in shape and measurements

PN803W	White	<b>\$103.95</b>
PN803B	Black	

- Capacity: 150kg
- Graduation: 100g
- Functions: auto-HOLD, tip-on automatic switch-on, automatic switch-off
- Two-component rubber coating. For pleasant footing
- Four stable rubber feet. For safe standing
- 28-millimeter LCD figures. For easy reading
- 2 year warranty
- Optional carry bag: Seca 415



**seca 876**

Flat Scale For Mobile Use

**\$577**

SEC876

- Capacity: 250kg
- Graduation: 100g < 150kg > 200g
- Functions: mother/child function, kg/lbs switch-over, tap-on automatic switch-on, automatic switch-off.
- 2 year warranty
- Optional carry bag: Seca 415

**seca**

## seca 874Dr

User-friendly and specially designed for use in medical practices.

- Capacity: 200kg
- Graduation: 50 g < 150 kg > 100 g
- Dimensions: 321(w) x 61(H) x 361(D) mm
- 6 batteries, type AA
- 2 year warranty
- Dual display for patient and doctor
- Mother/child-function, Auto-HOLD
- Optional carry bag: Seca 415



**\$545**

SEC874

**seca**

## SOEHNLE Digital Personal Scale

Simple use for mobile weighing

- Max. Capacity 150 kg
- Less weight, 1.4kg
- Easy-to-read LCD display, 26 mm digit height
- Hold function: The weight value will be held for approx. 6 seconds
- Automatic zeroing at switch-on
- Units in kg, lb or st, switchable at the bottom of the scale
- Energy-saving auto-off function after approx. 8 seconds



	Each	3+
SO7810	\$185.90	<b>\$141.90</b>

**SOEHNLE**  
PROFESSIONAL

# SOEHNLE

## PROFESSIONAL

# Superior Quality and Reliability



SOEHNLE Industrial Solutions is one of the leading European companies in professional weighing and measuring technology.

Founded in 1868, SOEHNLE Industrial Solution has 150 years of experience developing accurate measuring devices in the private and commercial sector, and their strong, innovative spirit continues to drive the development of new concepts and designs.

SOEHNLE scales combine renowned German technological excellence with sophisticated design.

Distributed by





## seca 644

Digital Bariatric scale with handrail and BMI function

- Capacity: 360kg
- Graduation: 50 g
- Dimensions:  
640(w) x 1,150(H) x 660(D) mm
- Rechargeable batteries
- 2 year warranty
- Optional measuring rod: Seca 223

**\$1,675**

PN644



seca

## seca culina 852

Digital portion and diet scales

- Capacity: 6.6 lbs / 3 kg
- Graduation: 0.05 oz / 1 g
- Dimensions (WxHxD): 210 x 140 x 210 mm
- Dimensions/scale platform (WxHxD): 140 x 53 x 140 mm
- Weight: 710 g
- Functions: TARE, automatic switch-off
- Power supply: Batteries
- 2 year warranty
- Includes Bowl with food safe surface

**\$131.25**

ECO852



seca


## seca 719

With four load cells for increased stability

- Capacity: 180kg
- Graduation: 100g
- Display on a high column for easy reading
- High stability for demanding users.
- Functions: auto-HOLD, step-off, tip-on automatic switch-on, automatic switch-off
- Battery operated
- 2 year warranty

**\$346**

PN719



seca

## seca 954

Digital chair scale with wireless transmission

- Capacity: 300 kg
- Graduation: 50 g < 150 kg > 100 g
- Dimensions (WxHxD): 563 x 897 x 978 mm
- Weight: 24.5 kg
- Power supply: Power adapter/rechargeable batteries
- Data transmission: seca 360° wireless technology
- Functions: TARE, pre-TARE, HOLD, auto-HOLD, BMI, auto-CLEAR, lbs/kg switch-over, damping, automatic switch-off, SEND/PRINT
- 2 year warranty

**\$2,515**

SE954



seca

## seca 952

Chair scale for weighing while seated

- Capacity: 200kg
- Graduation: 100g
- Functions: TARE, auto-HOLD, kg/lbs switch over, automatic switch-off
- Battery operated
- 2 year warranty

**\$1,570**

PN952



seca



## SOEHNLE Mechanical Personal Scale

- Large and clear circle pointer display
- Easy to read in kg-steps
- Manual zero resetting possible
- Attractive chrome trim ring
- Sturdy steel housing
- Weighing platform with easy-to-clean plastic covering
- Capacity: 160kg



**\$141.<sup>90</sup>**

SO6163

**SOEHNLE**  
PROFESSIONAL

## seca 786

Mechanical Column Scales With Large Round Dial

- Capacity: 150kg
- Graduation: 500g
- Functions: reset-to-zero function
- 2 year warranty
- Optional: measuring rod seca 224
- 2 year warranty



**\$545**

PN786

**seca**

## seca 750

Robust Mechanical Floor Scales

Pretty and stable with their pure, light design to make a good impression in any clinic. Powder coated steel casing, circular display and an easy to clean tread area of imitation leather.

- Capacity: 150kg
- Graduation: 1kg
- Dimensions (mm): 317W x 117H x 477D
- Functions: reset to zero function
- 2 year warranty



		Each	3+
PN750WW	White	\$177.45	<b>\$156.45</b>
PN750WB	Black		

**seca**

## seca 762

Robust Mechanical Personal Scales with Fine Graduation

- Capacity: 150kg
- Graduation: 500g
- Reset-to-zero function
- 2 year warranty



**\$209**

PN761

**seca**

## seca 760

seca 760 mechanical flat scales, with large dial face.

seca 760 mechanical flat scales, with large dial face. Traditional weighing mechanics in classic colours. With chrome accents, the seca 760 is available in burgundy, midnight blue, ecru, as well as the traditional black and white – to suit most decors. The seca 760 has a platform of skin friendly grooved rubber with the classic scale dial and robust, precise weighing mechanism.

- Capacity: 150 kg
- Graduation: 1 kg
- Available in burgundy red/chrome or midnight blue/chrome
- Functions: reset-to-zero function



PN760B	Midnight Blue	<b>\$237</b>
PN760R	Burgundy Red	

**seca**



## seca 310 Hanging Baby Scale

Mechanical circular dial scale for use far from any hospital

- Capacity: 25 kg
- Graduation: 50 g
- Dimensions (WxHxD): 220 x 225 x 52 mm
- Reset-to-zero function
- 2 year warranty



**\$136.50**

SEC310

**seca**

## seca 333i baby scale and optional measuring rod with digital display.

The seca 333i makes a good impression with its outstanding quality and simple operation.

- Completely EMR ready with Wi-Fi technology.
- Capacity: 20 kg
- Graduation: 5 g
- Dimensions (WxHxD): 650 x 110 x 308 mm
- Weight: 3.7 kg
- Power supply: Power adapter, batteries
- Suits most USB barcode scanners for patient ID



### Digital measuring rod:

- Measuring results are easy to read because of the illuminated LCD display.
- Includes head and foot positioners.
- Ergonomically shaped adjustment which prevents the rod from tilting

PN331L	Baby Scale with digital display	<b>\$1,308</b>
PN234	Measuring rod with digital display	<b>\$583</b>

**seca**

## Carts for Baby Scales

Sturdy mobile support carts with customised indentations to suit all seca baby scales

### seca 402:

- Dimensions (WxHxD): 715 x 800 x 435 mm
- Weight: 11.5 kg
- Max. load: 100 kg
- 4 bumpers
- 2 brakes
- Optional: paper roll holder code SE408

### seca 403:

- Dimensions (WxHxD): 715 x 800 x 435 mm
- Weight: 27.5 kg
- Max. load: 100 kg
- 4 bumpers
- 2 brakes
- Integrated drawer and shelf
- 2 year warranty

SE402	Cart Only	<b>\$762</b>
SE403	Cart with shelf and drawer	<b>\$1,045</b>
SE408	Paper towel holder for scales cart	<b>\$177.45</b>

seca 402



seca 403

**seca**

## Cases for Baby Scales



SE431	431 Backpack for baby scale	<b>\$152.25</b>
PN428	428 Carry Case for baby scale	<b>\$183.75</b>

**seca**

## SOEHNLE Easy - Fold Baby Scale

Foldable baby scale – easy to transport

- Patented folding mechanism with magnetic lock
- Low space requirement, easy to transport
- Easy to read LCD display
- Tare function: Covers/paddings can be tared
- Hold function: The weight value will be held for a few seconds after unloading
- Switchable to kg/lb units
- Energy-saving auto-off function
- Battery operation with 2 x 1.5 V AAA alkaline batteries (included in scope of delivery)



**\$329**

SO8320

**SOEHNLE**  
PROFESSIONAL





## seca 374

Wireless Baby Scale With Extra Large Weighing Tray



- Capacity: 20kg
- Graduation: 5g < 10kg > 10g
- Functions: TARE, BMIF, auto-HOLD, auto-CLEAR, kg/lbs switch-over, damping, automatic switch-off, SEND/PRINT
- Optional: measuring rod seca 233, head and foot positioners seca 418, switch-mode power adapter seca 400
- 2 year warranty

### \$913

PN374



SEC233	Optional measuring rod for 374	<b>\$451</b>
--------	--------------------------------	--------------

## seca 354

Electronic baby scales with fine graduation, also useable as flat scales for children.

- Capacity: 20kg
- Graduation: 5g < 10kg > 10g
- Functions: TARE, BMIF, auto-HOLD, kg/lbs switch-over, automatic switch-off
- 2 year warranty
- Optional: measuring board seca 417, carrying case seca 414 (for combination with measuring board) and seca 413

### \$205

PN354



## seca 334

Mobile Electronic Baby Scale

### Scale

- Capacity: 20kg
- Graduation: 5g < 10kg > 10g
- Functions: TARE, BMIF, auto-HOLD, kg/lbs switch-over, automatic switch-off
- 2 year warranty
- Optional: measuring rod seca 232, head positioner seca 419, switch-mode power adapter seca 400, carry case seca 428 for scale only. (Or larger carry case seca 414 for scale and measuring rod)

### \$703

PN334



## seca 232

Measuring rod for SECA 334 baby scale.

- Measuring range: 35-80cm
- Graduation 1mm
- 2 year warranty

### \$415

SEC232



## seca 413

Carrying case for seca 354

### \$99.75

PN413





## Soehnle Height Rod - Electric

Handy height rod for mobile use

- Ultrasound measuring technology, temperature-compensated
- Tilt sensor for correction if measuring instrument is not horizontal, resulting in an error message
- Low space requirement due to fold function
- Plastic housing
- LCD display, 9 mm digit height
- Acoustic acknowledgment at end of measurement and reading displayed
- Units in cm or inch
- Battery operation with 2 x 1.5 V AA batteries (included in scope of delivery instruments)



**\$274**

SO5003

**SOEHNLE**  
PROFESSIONAL

## Soehnle Measuring Tape - Mechanical

Small and handy

- Measuring tape for wall-mounting
- No wear parts and maintenance-free
- Units in mm

**\$31.90**

SO5002

**SOEHNLE**  
PROFESSIONAL



## seca 210

Mobile Measuring Mat For Babies And Toddlers

- Measuring range: 10-99 cm
- Graduation: 5mm
- 2 year warranty

**\$120.75**

PN210



## seca 218

Disposable Measuring Tape

- Made of high-quality paper to prevent cross infection
- Simple removal of measuring tape for immediate use
- Always at your fingertips thanks to the exclusive wall dispenser that can be simply installed anywhere
- Exact results read easily
- Scale from 0-100 cm
- Wall dispenser included
- 500 disposable tapes included
- Refills (Pkt/1000) also available

**\$124.95**

PN218



## seca 211

Disposable Head Circumference Tapes

- Made of sturdy material to prevent cross infection
- Always at your fingertips thanks to the exclusive wall dispenser that can be simply installed anywhere
- Exact results read easily
- Scale from 12-59 cm
- Wall dispenser included
- 500 disposable tapes included
- Refills (Pkt/1000) also available

**\$78.75**

PN211



## seca 212

Measuring Tape For Head Circumference of Babies and Toddlers

- Measuring range: 5-59cm
- Graduation: 1mm
- 2 year warranty

**\$9.40**

SEC212



**seca**

## Retractable Measuring Tape - 1.5 metres

**\$3.95**

SAS423



## seca 206

Mechanical Measuring Tape

- Measuring range: 0-220 cm
- Graduation: 1mm
- 2 year warranty

**\$36.75**

PN206



**seca**

## seca 201

Retractable circumference measure tape

- Measuring range: 0-205cm
- Graduation: 1mm
- 2 year warranty

**\$21**

SEC201



**seca**