

## OnTrak, the **New Evolution** of ABPM from Spacelabs Healthcare



**Includes:  
2 x BP cuffs**

1 x standard adult size  
BP cuff and 1 x extra  
large adult size  
nylon BP cuff

### Vigilance with Confidence

OnTrak is the result of over forty years of expertise in the research and development of ambulatory blood pressure recording.

Spacelabs OnTrak Ambulatory Blood Pressure monitor has been comprehensively tested and validated against all three of the recognised international protocols, to give you reassurance and confidence when assessing a patient's hypertension.

### Ease of Use

- Large 2.4" color screen makes it easy to program giving full control of functionality.
- Protocol, Child Mode and Maximum initial inflation Pressure can be set on the unit.
- Monitor your patient for up to seven days or up to 300 readings on two AA batteries.
- The 'Comfort Mode' setting controls initial inflation pressure reducing patient anxiety and thus enhancing compliance.



## OnTrak ABP System

The New Evolution of ABPM

- Uses 2 standard AA batteries
- Record up to 300 readings on single set of batteries for multi-day study use
- Robust reusable cuffs
- Cuff covers for optimal patient comfort and cuff protection
- Error code reporting with reason descriptions
- Maintenance contracts available for maximum up time
- BHS: A/A rating
- ESH: ESH-IP2 recommended for clinical use on adult population
- AHA: Meets criterion 1 and 2 of ANSI/AAMI/ISO 81060-2:2009 for both adults and children
- Includes Sentinel ABP Reporting Software
- Optional - connect ECG's & Holters to Sentinel
- Optional - Sentinel Network Capabilities
- Supplied with 2 cuffs: Adult and Large Cuffs

**Includes:  
2 x BP cuffs**  
1 x standard adult size BP cuff and 1 x extra large adult size nylon BP cuff

**\$3,145**

EC90227



## ABPM-7100 Complete System

Including Software

- BHS A/A Intelligent inflation provides excellent patient comfort and improves compliance
- Fully flexible program options for various measurement periods and inflation
- Includes Welch Allyn CardioPerfect Workstation Software - the professional's choice
- Easily create, print and customise detailed reports
- 100% latex-free and PVC cuffs reduces risk of allergic reactions
- Can be integrated with other CardioPerfect devices such as ECG and Spirometry to streamline diagnostics
- Supplied with 2 cuffs - adult and large adult

	Each	3+
WL7100S	\$3,670	<b>\$3,276</b>



## HEM-907 Blood Pressure Monitor

The HEM-907 is a top-of-the-range digital blood pressure monitor. Featuring Intellisense™ technology, it can function via both Automatic Oscillometric and Automatic Auscultation measurement mode. Designed for hospital, general practitioners and other medical professionals.

- Dual functions – Automatic Oscillometric and Automatic Auscultation
- Accuracy and precision, complying with AAMI standards
- Quiet operating system, ideal for clinical environment
- Easy to operate using intellisense technology
- Consecutive measurement option with variable waiting and interval time

	Each	5+
DA907	\$895	<b>\$625</b>



Includes small/medium/large cuff & adaptor

## HEM7322 Premium

Automatic Blood Pressure Monitor

- Simple onetouch operation
- 90 reading memory with date and time
- Enhanced Intellisense technology
- Dual Check System
- Displays the average of the last 3 readings
- Detects irregular heart beat & body movement
- Blood Pressure Level Indicator
- Flexible fan shaped cuff
- Carry bag included & Cuff size (Med –Lrg) included
- Long Battery Life – 1,500 readings with new alkaline batteries
- 4 years extended warranty when registered online
- Optional AC adapter



	Each	5+
DAHEM7322	\$195	<b>\$175</b>

## Revolutionized office Blood Pressure




MESI Dual BP mTABLET MODULE

- First tubeless professional BP cuff supporting 3 cuff sizes
- The right pressure obtained by controlled inflation and deflation
- Advanced analysis options with pulse waveforms
- Advanced error detection system
- Module must be used with the mTABLET (refer to page 100)

**\$750**

CUFFSYS

ADDITIONAL APPS AVAILABLE:



# OMRON

# Global No. 1 Blood Pressure Monitors

Professional  
Blood Pressure Monitor  
**DA907**

Professional  
Blood Pressure Monitor  
**DAHBP1320**



## Irregular Heartbeat Detection

OMRON Blood Pressure Monitors can detect an irregular heartbeat (arrhythmia) which can assist with early diagnosis\*

*\*N.B. Only a doctor can diagnose arrhythmia*



Smart Elite +  
Blood Pressure Monitor  
**DAHEM7600T**



Ultra-Premium  
Blood Pressure Monitor  
**DAHEM7320**



Premium  
Blood Pressure Monitor  
**DAHEM7322**



Standard  
Blood Pressure Monitor  
**DVSEM7121**



Premium Wrist  
Blood Pressure Monitor  
**DAHEM6221**

For more information visit [www.omronhealthcare.com.au](http://www.omronhealthcare.com.au)

\*For people with high blood pressure. Consult your doctor to evaluate the readings. Check your device periodically for accuracy.  
\*Do not use these devices in combination with electronic implants such as pacemakers.



## HBP1320 Professional Blood Pressure Monitor

- IntelliSense™ Technology
- Digital Blood Pressure Monitor
- Non-invasive
- Portable

Package includes the following:

- HBP1320 Professional Blood Pressure Monitor.
- HBP1320-CUFFM- Cuff Medium (22-32cm) inc bladder and cuff cover
- HBP1320-CUFFL-Cuff Large (32-42cm) inc bladder and cuff cover
- HBP1320-PBAT- Rechargeable battery pack
- HBP1320- ADAPT- AC Adapter / battery charger

**NEW**



	Each	5+
DAHBP1320	\$535	<b>\$515</b>



## Omron HEM6161 Basic Wrist Blood Pressure Monitor

The HEM-6161 makes it easy to monitor blood pressure while on-the-go. This basic model comes with cuff wrap guide and irregular heartbeat detection to ensure the monitor is used correctly.

**NEW**



	Each	5+
DAHEM6161	\$145	<b>\$135</b>



## HEM7121 Standard Automatic Blood Pressure Monitor

The Omron Standard blood pressure monitor provides simple and accurate blood pressure readings without compromising user comfort. Cuff size (Med –Lrg) included.

- **One Touch Use** Easy one touch operation allows home blood pressure monitoring that is both simple and accurate
- **30 Memory Storage** 30 memory storage allows the user to review the last 30 readings with a touch of a button
- **Reading time** Approximately 50 seconds
- 4 years extended warranty when registered online



**\$119**

DVSEM7121



## HEM7320 Ultra Premium Automatic Blood Pressure Monitor

- Cuff Wrapping Guide
- Enhanced Intellisense technology
- Simple one touch operation
- 90 reading memory with date and time
- Displays the average of the last 3 readings
- Detects irregular heart beat
- Detects body movement
- Blood Pressure Level Indicator
- Long Battery Life –1,000 readings with new alkaline batteries
- Carry bag included
- Easy Fit Cuff (Med –Lrg) included
- 4 years extended warranty when registered online
- Optional AC adapter

**\$230**

DAHEM7320



## Smart Elite+ Blood Pressure Monitor

- Tubeless Technology
- Bluetooth connectivity
- Automatically applies the right amount of pressure
- 100 Data memory with date and time via OMRON Connect App
- Snaps easily to the arm
- Indicator confirms the cuff has been wrapped correctly
- If an irregular heartbeat is detected an indicator icon will appear
- Alerts user if any body movement has occurred during the reading.



	Each	5+
DAHEM7600T	\$295	<b>\$265</b>





Hillrom™

Welch Allyn  
ProBP™ 2000

Designed and validated<sup>1</sup> for healthcare professionals

## Accurate blood pressure readings in ~20 seconds<sup>1</sup>

<sup>1</sup>ALPERT B, Validation of the Welch Allyn Pro BP 2000, a professional-grade inflation-based automated sphygmomanometer with arrhythmia detection in a combined pediatric and adult population by ANSI/AAMI/ISO standard testing, Blood Pressure Monitoring, 2018 Vol 23 (6) p. 315-317.



ACCURACY



COMFORT



SPEED

IF YOU KNOW WELCH ALLYN,  
YOU KNOW HILLROM

Hill-Rom Pty Ltd  
Suite 4.01, 2-4 Lyonpark Road, Macquarie Park NSW Australia

Always read the labels and instructions for use.

© 2020 Welch Allyn, Inc. ALL RIGHTS RESERVED. APR53201-AU\_r12/02/2020 ENG

hillrom.com.au

## Connex® ProBP™ 3400 Digital Blood Pressure Device

Accurate, reliable, connected blood pressure readings in the palm of your hand.

- Small, handheld design makes device easy to carry and capture readings on rounds or in triage settings
- Hospital-grade SureBP™ technology enables fast, accurate 15-second blood pressure readings, improving patient comfort
- USB-wired technology enables connectivity to select EMR/HIS systems
- Compatible with blood pressure cuff sizes ranging from small child (12 cm) to thigh (55 cm)
- Lithium ion battery technology delivers 100 readings on a single charge with rapid recharge time
- Mercury and latex free to ensure clinician and patient safety
- Wall mount, desk mount or mobile stand available
- Includes Flexiport Adult and Large Adult BP Cuffs



	Each	3+
WL34XFHT-6	\$1,045	<b>\$1,039</b>

WelchAllyn®

## Connex® ProBP™ 2000

Digital Blood Pressure Device

### Features

- Blood pressure readings are captured on inflation in about 20 seconds Meets AAMI standards for blood pressure accuracy
- Compatible with Welch Allyn FlexiPort Reusable Blood Pressure Cuffs, from child size 9 up to high size 13 (15 - 55 cm range)
- Device comes with one adult size 11 FlexiPort Reusable Blood Pressure cuff
- Meets AAMI standards for blood pressure accuracy



		Each	3+
WL2000	ProBP 2000 Battery Operated	\$262	<b>\$241</b>
WL2000P	ProBP 2000 Power Supply included	\$314	<b>\$304</b>

WelchAllyn®

10% GST will apply to all items unless they are marked GST Free



## GAMMA G7®

The top-of-the-line model with aluminium alloy construction

- Latex-Free
- Exclusive air-release valve. Provides fine and precise air-release rate adjustment, quick deflation and extreme durability. A rubber-stop prevents valve-lock
- Exceptionally strong housing construction with triple-mounted, shockproof movement
- Impact-resistant housings crafted in special, shock-absorbing materials incorporate the latest advances in injection-moulding technology
- Pressure-proof manometer. Ensures durability and reliability
- Micro-filter. Protects valve and manometer
- Easy-to-read 56mm diameter scale
- 50% larger bulb volume for fast inflation of the cuff
- Light alloy, adjustable spoon. For both right or left-handed operation

The HEINE G7 is unique in its advanced construction. The all-metal housing makes it virtually indestructible.



**\$252**  
HEM0009237



## GAMMA GP®

Combines all the design and quality features of the GAMMA G5, but with a push-button valve. Fine and precise air-release rate adjustment, ultra-fast air-release, practically wear-free. 5 year manufacturers warranty



	Each	4+
HEM0009240	\$214	<b>\$196</b>



## GAMMA G5®

The HEINE standard sphygmomanometer. Virtually the same as the G7 Model except the housing is moulded of 2 component thermo plastic. 5 year manufacturers warranty



	Each	4+
HEM0009235	\$173	<b>\$167</b>



## GAMMA XXL LF®

The versatile professional model

- Latex-free
- Durable polycarbonate construction
- Will hold up in the busiest hospital and practice environments
- Oversized 130mm diameter scale
- Can be pivoted for clear, easy reading from any angle
- Long spiral inflation tube for increased working distance
- Precisely calibrated and durable
- Precision air-release valve
- Pressure-proof manometer. Ensures durability and reliability
- Wall, desk, wheeled stand. Adjustable height



Freight charges will apply

	Each	4+
HEM0009310 Desk (Non-spiral tubing)	\$307	<b>\$277</b>
HEM0009311 Wall	\$327	<b>\$319</b>
HEM0009313 Mobile	\$654	<b>\$648</b>



## DuraShock DS44

This integrated-style gauge features a FlexiPort cuff-mounted design that's compact and lightweight. The gauge can be rotated to any orientation desired by the clinician.



**\$166.95**

WLDS44-11

WelchAllyn®

## DuraShock DS54 + Adult Cuff

The Welch Allyn DuraShock DS54 provides a traditional thumbscrew-style air release valve.

- Gear-free DuraShock™ technology
- Shock-resistant for accuracy and lower cost over time
- 5-year calibration warranty
- Not made with natural rubber latex



**\$145.95**

WL5401189

NEW

WelchAllyn®



## Aneroid Sphygmomanometer

SA0569	Economy Aneroid	<b>\$29</b>
SA0405	Professional Aneroid	<b>\$29</b>



## Palm Sphygmomanometer

The CM-3414 is a highly accurate. This unit can be converted for left or right hand use and it has an easy to read 65mm diameter dial.



**\$59**  
SA0401

## Deluxe Aneroid Sphygmomanometer Kit - Latex Free

- Dual head Stethoscope
- Nurse Utility Scissors
- Pen Light
- Tape Measure
- Stethoscope Name Tag
- Carry Case

**\$79**  
AD2000RED



## ABN Aneroid Sphygmomanometer

- Large 15cm square easy-to-read dial
- Dial face adjustable
- Cuff, bladder, bulb, valves and coiled tubing all incl.
- Practical half round storage basket

SA0580	Desk Model	<b>\$159</b>
SA0568	Wall Model	<b>\$159</b>
SA0579	Floor Model	<b>\$395</b>

\*Freight charges will apply



You can place your order via fax on **1300 557 631**.

## Shock Proof Palm CM-3412 Sphygmomanometer

The CM-3412 utilises shockproof technology to provide a cost effective rugged latex free device. Ideal for use in field situations it has a black ruggedized ABS plastic case with 55mm diameter dial.

Its single tube latex free inflation system allows smooth and accurate readings. The unit is supplied in a vinyl zipper carry case.

**\$69**  
SA0582



## Babyphon®

Child-friendly all in blue – that's Babyphon®! The palm-type sphygmomanometer specially designed for paediatricians. With a set of 3 blue, calibrated latex free cuffs for children, infants and new-born. 2 years warranty.

**\$177.45**  
RIE1010017





## SPECTRUM Professional with Stethoscope

**\$79**

SA0406



## Home blood pressure kit with Stethoscope

**\$49**

SA0570



## Diplomat® – Desk Model Mercurial Sphygmomanometer



Extremely precise, constantly improved and with impressive features. The 99.999% mercury purity guarantees maximum measuring accuracy.

- High quality
- Micro filter for long life
- Precision air-release valve, wear-free and with fine regulation
- Mercury lock
- Special seal against mercury contamination
- Made in Germany
- 2 year warranty
- Includes Adult cuff latex free

**LATEX FREE**



		Each	2+	5+
RIE1010027	Blue	\$209	\$199	<b>\$178</b>
RI1002LF	Silver	\$209	\$199	<b>\$178</b>

## Empire® N Mercury Sphygmomanometer with velcro cuff

Mercury sphygmomanometer with ABS plastic casing.

An attractive design combined with sturdy everyday-use quality. Available as desk, wall, or floor model. Reasonably priced, practical and safe. Attractive optics and tested everyday practicality to cover every situation in general practice. Its specially designed stand provides a very high level of stability.

Full range of spares available.

2 years warranty.

		Each	4+
RIE1102LF	Wall, adult size velcro cuff	\$251	<b>\$241</b>
RIE1122LF	Floor Ø55 cm, adult size velcro cuff	\$489	<b>\$478</b>

## Big Ben® Analogue Sphygmomanometer

- Generously dimensioned scale (Ø 147.2mm) and high-contrast scale layout for optimum readability
- Spacious basket on the reverse side for tubes and cuffs
- Micro-filter to protect the release valve and measurement system
- Specially hardened copper-beryllium membrane, which is resistant to aging, can withstand pressures of up to maximum 600mmHg
- Precision mechanism mounted on a special support, in order to ensure complete stability of the measuring system and the zero-point
- Available in wall, desk or floor model
- 2 years warranty



		Each	4+
RIE1010002	Wall LF	\$183.75	<b>\$173.25</b>
RIE1010005	Trolley Lf	\$346	<b>\$335</b>
RIE1010001	Desk	\$209	<b>\$199</b>
RIE10392	Basket	\$93.45	N/A

**LATEX FREE**







## GAMMA G7/G5/GP Latex-Free Cuff & Bladder



				Each
HEM0009637	Child	10x34cm	13-20cm	<b>\$74</b>
HEM0009617	Adult small	14x53cm	20-29cm	<b>\$74</b>
HEM0009613	Adult	14x58cm	29-41cm	<b>\$76</b>
HEM0009706	Adult Large	17x70cm	35-47cm	<b>\$101</b>
HEM0009615	Thigh	20x86cm	>41cm	<b>\$100</b>

## GAMMA XXL LF® Latex-Free Cuff & Bladder

				Each
HEM0009703	Child	10x34cm	13-20cm	<b>\$79</b>
HEM0009705	Adult small	14x53cm	20-29cm	<b>\$79</b>
HEM0009623	Adult	14x58cm	29-41cm	<b>\$81</b>
HEM0009631	Adult Large	17x70cm	35-47cm	<b>\$104</b>
HEM0009625	Thigh	20x86cm	>41cm	<b>\$104</b>

## Accessories for HEINE EN200 BP

			Each
HEM0009803	Quick Cuff Adult small	<b>\$78</b>	
HEM0009804	Quick Cuff Adult		
HEM0009805	Quick Cuff Adult large		
HEM0009806	Quick Cuff Adult XLarge		

Quick Cuff - Single Patient		Box/10
HEM0009813	Small	<b>\$309</b>
HEM0009814	Adult	
HEM0009815	Large Adult	



## Bladderless Cuffs

- Latex Free
- High quality, durable material
- Washable and easily maintained.

SA0464	Child	19 - 27cm	<b>\$20.25</b>
SA0465	Adult	26 - 41cm	<b>\$21</b>
SA0466	Large Adult	34 - 51cm	<b>\$33</b>



## Accessories for OnTrak BP

### Disposable Cuffs

		Box/20
SL-ODSA	Adult small	<b>\$329</b>
SL-ODAD	Adult	<b>\$398</b>
SL-ODLA	Adult large	<b>\$405</b>
SL-ODXL	Adult XLarge	<b>\$662</b>



### Ultracheck Reusable Cuffs

SL-ONCH	Child	<b>\$258</b>
SL-ONSA	Small adult	<b>\$268</b>
SL-ONAD	Adult	<b>\$294</b>
SL-ONLA	Adult large	<b>\$331</b>
SL-ONXL	Adult XLarge	<b>\$331</b>



### Reusable Pouch for OnTrak BP

SL-OPOU	Pouch	<b>\$107.10</b>
SL-OSST	Shoulder strap	<b>\$68.25</b>



### Disposable Cuff Cover

		Box/50
SL-OCCO		<b>\$216</b>



## Sphygmomanometer Spares

Cuffs			
SA0453	New Born	18 x 5.5cm	<b>\$10.35</b>
SA0474	Child	34 x 11cm	<b>\$12</b>
SA0577	Adult	53 x 15cm	<b>\$14</b>
SA0485	Large Adult	62 x 18cm	<b>\$22</b>
SA0524	Obese/Thigh	82 x 21cm	<b>\$28</b>
Blisters 2 Tube			
SA0454	New Born	9 x 4cm	<b>\$16</b>
SA0475	Child	18 x 9cm	<b>\$14</b>
SA0455	Adult	22 x 12cm	<b>\$8</b>
SA0487	Large Adult	32 x 15cm	<b>\$19</b>
SA0523	Obese/Thigh	38 x 18cm	<b>\$29</b>
Complete Set (bladder, cuff, bulb, valve and connector)			
SA0469	Child		<b>\$32</b>
SA0470	Adult		<b>\$36</b>
SA0471	Large Adult		<b>\$42</b>
SA0472	Obese		<b>\$48</b>
Spares			
SA0574	Control Valve (air release valve)		<b>\$4.95</b>
SA0572	Bulb		<b>\$5.95</b>
SA0576	Bulb Sphyg Large with Valve		<b>\$6.50</b>
SA0578	Tubing Coiled 2m		<b>\$15</b>
Connectors			
SA0525	Plastic black		<b>\$2.50</b>



## FlexiPort™ Reusable Blood Pressure Cuffs

**WelchAllyn**

### With Two-Tube Tri-Purpose Connector & Empty Tube

WL072TP	Infant	<b>\$54.60</b>
WL082TP	Small Child	<b>\$54.60</b>
WL092TP	Child	<b>\$61.95</b>
WL102TP	Small Adult	<b>\$58.80</b>
WL112TP	Adult	<b>\$61.95</b>
WL122TP	Large Adult	<b>\$82.95</b>
WLREUSE12T	Large Adult Long	<b>\$91.35</b>

### With Screw-Type Connector & One-Tube

WL091SC	Child	<b>\$57.75</b>
WL101SC	Small Adult	<b>\$53.55</b>
WL111SC	Adult	<b>\$61.95</b>
WL121SC	Large Adult	<b>\$79.80</b>
WL131SC	Adult Thigh	<b>\$98.70</b>

### with No Tubes or Connector (Raw Cuff)

WLREUSE08	Small Child	<b>\$46.20</b>
WLREUSE09	Child	<b>\$48.30</b>
WLREUSE10	Small Adult	<b>\$50.40</b>
WLREUSE11	Adult	<b>\$57.75</b>
WLREUSE11L	Adult Long	<b>\$55.65</b>
WLREUSE12	Large Adult	<b>\$82.95</b>
WLREUSE12L	Large Adult Long	<b>\$86.10</b>
WLREUSE13	Adult Thigh	<b>\$100.80</b>



## Sphygmomanometer Accessories

WL5088-01	Air Release Valve and Bulb	<b>\$65.10</b>
WL5082-168	Tri-Purpose Sphygmomanometer Connector, 10/Bag	<b>\$26.25</b>
WL5089-13	Coil Tubing	<b>\$74.55</b>
WL460061	Mobile Stand for Connex BP	<b>\$461</b>
WL5082-189	13" Flexiport Tubing Pack/10	<b>\$18.90</b>
WL17-31-90	Bulb for Chromed Valve	<b>\$7.88</b>

WL16-40-00	Sphyg Adaptor Male Luer	<b>\$5.78</b>
WL107195	ProBP 2000 Stand Kit	<b>\$367</b>
WL107194	ProBP 2000 Wall Mount Kit	<b>\$57.75</b>
WL107201	ProBP 2000 Holster Kit	<b>\$36.75</b>
WL107041	ProBP 2000 AC Adapter	<b>\$61.95</b>

**WelchAllyn**

## FlexiPort Reusable Blood Pressure Cuffs Multi-Cuff Kit

**WelchAllyn**

Long lasting, easy to use, and affordable, Welch Allyn cuffs are the top choice for hospitals and clinics around the world. Why? Because Welch Allyn can offer you virtually every type, size, and configuration of cuffs, along with the ability to simplify inventory and reduce confusion with our innovative FlexiPort connection, which makes every cuff a one or two-tube cuff.

- Primary sizes are colour coded to make it easy to find the proper cuff size
- Easy one-hand attachment/detachment of tubing with FlexiPort makes changing cuffs much faster and easier
- Latex-free material minimizes risk of allergic reactions
- Rotatable port reduces stress to cuff tubing and port, improving patient comfort and cuff durability
- Folded edge reduces risk of cuts and scrapes, maximizing patient comfort
- Meets all latest clinical guidelines for proper fit from AAMI and AHA
- 3 Year Warranty



**\$193.<sup>20</sup>**

WLREUSEMLT



## HEM-907 Accessories

DAHEM-907R	Bladder for Small Cuff	<b>\$17</b>
DAHEM-907D	Bladder for Medium Cuff	<b>\$20</b>
DAHEM-907H	Bladder for Large Cuff	<b>\$32</b>
DAHEM-907C	Cuff Small 17-22cm	<b>\$19</b>
DAHEM-907M	Cuff Medium 22-32cm	<b>\$19</b>
DAHEM-907K	Cuff Large 32-42cm	<b>\$29</b>
DAHEM-907X	Cuff Extra Large 42-50cm*	<b>\$159</b>
DVHEMTUBE	Air Tube 1.3m	<b>\$10.95</b>
DAHEM-907S	Mobile Stand	<b>\$385</b>
DAHEM-907W	Wall Mounting Kit	<b>\$265</b>
DAHEM-907B	Battery Pack	<b>\$42</b>
DAHEM-907A	AC Adaptor	<b>\$69</b>

\* Please note the DAHEM-907X is supplied with cuff cover, bladder and 120cm tube.



On Roller Stand



## ABPM Accessories

- Reusable and ultra soft ABPM cuffs.
- Suits most brands of holter ABP-monitors.
- Latex free/PVC free.
- Compliant with ANSI/AAMI SP10-2002, EN60601-01:1990

ECVA004	Adult Plus Single Tube PVC Bladder BP Cuff, 32 - 38 cm	<b>\$72.45</b>
ECVA003	Adult Single Tube PVC Bladder BP Cuff, 24 - 32 cm	<b>\$68.25</b>
ECVA001	Child Single Tube PVC Bladder BP Cuff, 14 - 20 cm	<b>\$57.75</b>
ECVA005	Large Adult Single Tube PVC Bladder BP Cuff, 38 - 55 cm	<b>\$75.60</b>
ECVMCNT2-1	Metal Female Bayonet, 5/32 inch. ID Barbed. For HP / Philips, Siemens, Spacelabs, Criticare, IVAC/AI, Aris and Datascope	<b>\$19.95</b>
ECVA002	Small Adult Single Tube PVC Bladder BP Cuff, 20 - 24 cm	<b>\$63</b>

## AC Adapter

Suitable with **OMRON** models: HEM - 7211, 7221, 6052, 7203 & T9 (old models: 1A2, SEM1, SEM2, T9P)

**\$31**  
DVADAPT

## Standard Arm Cuffs

DACUFFCHIL	Small Cuff (17 - 22cm)
DACUFFFSC	Standard Cuff (22 - 32cm)
DACUFFLARG	Large Cuff (32 - 42cm)

**\$39**

Suitable for **OMRON** models: HEM - 6052, 6221, 7211, 7221, 7203, 7320, 7322, 7121 & T9, M10IT

## Medium/Large Fit Cuffs

Suitable with **OMRON** models: HEM 7121/ HEM 7130/ HEM 7322/HEM 7320

**\$39**  
DACUFF7320

## Riester Sphygmomanometer Spares

Latex-Free Cuff & Bladder 2 Tube		
RIE1010050	Child 35.5 x 10cm	<b>\$51.45</b>
RIE1010045	Adult 54.5 x 14.5cm	<b>\$51.45</b>
RIE1010046	Small Adult 42 x 13cm	<b>\$61.95</b>
RIE1010047	Large Adult 70 x 15cm	<b>\$72.45</b>
Latex-Free Cuff & Bladder 1 Tube		
RIE1010040	Child 35.5 x 10cm	<b>\$51.45</b>
RIE1010035	Small Adult 42 x 13 cm	<b>\$51.45</b>
RIE1010037	Adult 54.5 x 14.5 cm	<b>\$72.45</b>

More Riester spares		
RIE10363	Control Valve	<b>\$25.20</b>
RIE10382LF	Tubing coiled Latex Free	<b>\$27.30</b>
RIE10373	Connector Male Metal	<b>\$8.93</b>
RIE10374	Connector Female Metal	<b>\$8.93</b>
RIE10372	Connector Set M/F Metal	<b>\$15.70</b>

