



3M™ Steri-Strip™ Adhesive Skin Closures

3M Steri-Strip Adhesive Skin Closures (Reinforced)

- Used in the closures of lacerations, surgical incisions and following early suture/staple removal
- Non-woven, rayon-backed material reinforced with filaments for strength



		Box/50	4+
3MR1540	3x75mm	\$42.10	\$39.25
3MR1541	6x75mm	\$42.10	\$39.25
3MR1542	6x38mm	\$42.10	\$39.25
3MR1546	6x100mm	\$47.50	\$44.95
3MR1547	12x100mm	\$47.50	\$44.95

3M Steri-Strip Elastic Skin Closures

- Designed to stretch to accommodate swelling, body movement and contours
- Soft, non-woven backing for patient comfort



		Box/50	4+
3M4541	6x75mm	\$45.95	\$43.80
3M4542	6x38mm	\$45.95	\$43.80
3M4549	12x50mm	\$45.95	\$43.80
3M4547	12x100mm	\$57.00	\$54.00

3M Steri-Strip Blend Tone Skin Closures

- Wound support with a less noticeable look
- Non-reinforced closures with a pressure-sensitive adhesive



		Box/50	4+
3M1551	6x75mm	\$59.00	\$56.00



3M™ Steri-Strip™ Wound Closure System

- Ideal for wounds when extra protection is required
- Combines 3M™ Steri-Strip™ Skin Closures with 3M™ Tegaderm™ Dressings into one product
- Provides water proof coverage



		Size	Box/25
3M8514	3M Steri-Strip Wound Closure System - 3 Closure/envelope	60x48mm	\$111.88
3M8516	3M Steri-Strip Wound Closure System - 5 Closure/envelope	60x48mm	\$147.91




LEUKOSTRIP®

Adhesive Hypoallergenic Wound Closure Strips

LEUKOSTRIP® is made of a polyamide material, which will allow excess exudate to escape without impacting on the adhesion of the hypoallergenic adhesive.

- Elastic polyamide textile
- Hypoallergenic adhesive
- Reliable adhesive strength



Popular Sizes		Box/50	2+	4+
BE2876	4x38mm	\$55.00	\$53.00	\$49.00
BE2903	4x76mm	\$55.00	\$53.00	\$49.00
BE2878	6.4x76mm	\$55.00	\$53.00	\$49.00

LEUKOSTRIP® combines flexibility with secure grip!

Longer Sizes		Box/50	2+	4+
BE2879	6.4x102mm	\$69.00	\$68.00	\$63.00
BE2880	13x102mm	\$89.00	\$86.00	\$79.00

Skin Tone LEUKOSTRIP® (Tan Colour)		Box/50	2+	4+
BE2882	4x38mm	\$59.00	\$56.00	\$52.00
BE2883	6.4x76mm	\$62.00	\$59.00	\$54.00





Leukosan® Strip - Sterile

Wound closure strips are an option for topical skin closure due to their atraumatic painless application, low infection risk and good cosmetic results. Thanks to its knitted elastic fabric, Leukosan Strip represents the latest state of development, adding yet more benefits to the list. The knitted elastic fibres of Leukosan Strip allow the wound to heal without additional trauma to the surrounding skin.

White

BE7262901	3x75mm	Box/50 Pouch of 5 Strips	\$47.46
BE7262902	6x75mm	Box/50 Pouch of 3 Strips	\$47.46
BE7262904	12x100mm	Box/50 Pouch of 6 Strips	\$63.21
BE7262905	25x100mm	Box/25 Pouch of 4 Strips	\$74.76
BE7262900	6x38mm	Box/50 Pouch of 6 Strips	\$47.46
BE7262903	6x100mm	Box/50 Pouch of 10 Strips	\$68.78


Natural

BE7262908	6x75mm	Box/50 Pouch of 3 Strips	\$60.27
BE7262909	12x100mm	Box/50 Pouch of 6 Strips	\$72.14
BE7262907	6x38mm	Box/50 Pouch of 6 Strips	\$48.41



Omnistrips® Skin Tone

Omnistrips® are flexible wound closure strips with rounded edges to prevent rolling or lifting prematurely. Omnistrips® are tan coloured making them cosmetically pleasing for many skin tones. The application is made easy with anchoring tabs and the acrylic adhesive reduces the possibility of allergic reactions.



Watch an application video

	Size	Box/50
4454681	3x76mm	\$49.00
4454682	6x38mm	\$49.00
4454683	6x76mm	\$49.00
4454684	6x101mm	\$65.00



CLOSURE TECHNOLOGY

CONFIDENCE IN CLOSURE



SUTURES



HISTOACRYL®



MANIPLER



Dafilon®

Dafilon is a non-absorbable monofilament suture made of Polyamide polymers available in blue, undyed and black. Dafilon sutures are indicated for use in skin closure, plastic and reconstructive surgery, microsurgery and ophthalmic surgery.



**TOP
SELLER**

Most common sizes:

		Box/36	2+
BNSH093221	3-0 19mm 45cm	\$149.35	\$143.17
BNSH093220	4-0 19mm 45cm		
BNSH093219	5-0 19mm 45cm		
BNSH093212	5-0 16mm 45cm		
BNSH093211	6-0 16mm 45cm		

		Box/12
BNSG093221	3-0 19mm 45cm	\$64.90
BNSG093220	4-0 19mm 45cm	
BNSG093219	5-0 19mm 45cm	
BNSG093212	5-0 16mm 45cm	
BNSG093211	6-0 16mm 45cm	

		Box/36	2+
BNSH093236	2-0 24mm 45cm	\$166.87	\$161.72
BNSH093235	3-0 24mm 45cm		
BNSH093207	5-0 12mm 45cm		

		Box/36	2+
BNSH093535	3-0 24mm 75cm	\$180.25	\$174.07
BNSH093534	4-0 24mm 75cm		

		Box/36	2+
BNSH093520	4-0 19mm 75cm	\$166.87	\$161.72
BNSH093519	5-0 19mm 75cm		

www.bbraun.com.au

B | BRAUN
SHARING EXPERTISE

Monosyn® Quick



Monosyn Quick is an undyed short term absorbable synthetic monofilament suture made of glyconate for the approximation of soft tissue with a fast healing process. Main indications include Plastic surgery, Gynaecology / Obstetrics, Urology and Paediatrics, Odontology, ENT and Maxillofacial.

		Box/36	2+
BNS0025003	3-0 19mm 45cm	\$243	\$236
BNS0025002	4-0 19mm 45cm		
BNS0025006	5-0 19mm 70cm	\$329	\$320

B | BRAUN
SHARING EXPERTISE

Monosyn®



Monosyn® is a mid-term absorbable synthetic monofilament suture made of glyconate. Monosyn® offers excellent handling properties and high knot security. Monosyn® is recommended for all soft tissue approximations.

Monosyn® Undyed		Box/36	2+
BBC0023205	3-0 19mm 45cm	\$243	\$236
BBC0023204	4-0 19mm 45cm	\$243	\$236
BNS023403	5-0 19mm 70cm	\$319	\$310

Monosyn® Violet		Box/12
BNSG022205	3-0 19mm 70cm	\$118.46
BNSG022504	4-0 19mm 70cm	\$123
BNSG023403	5-0 19mm 70cm	\$115.36

B | BRAUN
SHARING EXPERTISE

Novosyn Quick



- Polyglactin 910
- Polyglactin 370 + Calcium Stearate (smooth surface)
- Uniform coating surface provides optimal flexibility & pliability
- Compacted braiding structure results in less dead space, high tensile strength and optimal knot security

**Replaces
Safil® Quick**


		Box/36	2+
BNS0046235	3-0 24mm 70cm	\$359	\$339
BBC0046221	3-0 19mm 45cm	\$279	\$269
BBC0046220	4-0 19mm 45cm	\$255	\$249

B | BRAUN
SHARING EXPERTISE



Novosyn®

- Polyglactin 910
- Polyglactin 370 + Calcium Stearate (smooth surface)
- Uniform coating surface provides optimal flexibility & pliability
- Compacted braiding structure results in less dead space, high tensile strength and optimal knot security



B BRAUN
SHARING EXPERTISE

Replaces
Safil®

Novosyn® Violet		Box/36	2+
BBC0068235	3-0 24mm 70cm	\$237	\$232
BBC0068220	4-0 19mm 45cm	\$222	\$216
BBC0068219	5-0 19mm 45cm	\$222	\$216

		Box/12	
BBG1048235	3-0 24mm 70cm		\$98.86
BBG006842	4-0 24mm 70cm		\$91.88
BBG0068219	5-0 19mm 45cm		\$81.38

Novosyn® Undyed		Box/36	2+
BBC0069221	3-0 19mm 45cm	\$222	\$216
BBC0069235	3/0 24mm 70cm	\$222	\$216
BBC0069220	4-0 19mm 45cm	\$222	\$216
BBC0069219	5-0 19mm 45cm	\$222	\$216

Silkam®

Silkam is a coated non-absorbable suture made of braided natural silk filaments available in black and white. Silkam sutures are indicated for general surgery, skin closure, oral surgery, neurosurgery and ligatures.



		Box/12	
BNSG076221	3-0 19mm 45cm		\$65
BNSG076220	4-0 19mm 45cm		\$60.78
BNSG076219	5-0 19mm 45cm		\$65
BNSG076235	3-0 24mm 75cm		\$71.08

		Box/36	2+
BNSH076235	3-0 24mm 75cm	\$205	\$195
BNSH076221	3-0 19mm 45cm	\$159.66	\$154.51
BNSH076220	4-0 19mm 45cm	\$159.66	\$154.51
BNSH076219	5-0 19mm 45cm	\$159.66	\$154.51

B BRAUN
SHARING EXPERTISE

Optilene®

- Polypropylene + Polyethylene
- High Tensile Strength
- New dye formulation provides improved thread visibility
- High thread pig-tailing resistance and reduced curling
- Re-engineered elasticity providing controlled suture elongation
- Low memory effect

Replaces
Premilene®



		Box/36	2+
BBC3090206	6-0 12mm 45cm	\$324	\$315
BBC0090435	3-0 24mm 45cm	\$247	\$239
BB3090219	5-0 19mm 45cm	\$247	\$239
BB3090220	4-0 19mm 45cm	\$247	\$239

B BRAUN
SHARING EXPERTISE

Syneture® DERMALON

Non-absorbable, monofilament, nylon

		Box/36	2+
DG1756-41	3-0 19mm 45cm	\$249	\$242
DG1756-31	4-0 19mm 45cm		
DG1756-21	5-0 19mm 45cm		



MAXON

Polyglyconate, monofilament
synthetic absorbable

		Box/36	2+
DG6618-51	2-0 24mm 70cm	\$369	\$360
DG6608-41	3-0 19mm 45cm		
DG6608-31	4-0 19mm 45cm		
DG6618-31	4-0 24mm 70cm		
DG6608-21	5-0 19mm 45cm		





SharpPoint Nylon

SharpPoint Plus Nylon Black Monofilament		Box/12	3+	5+
CRA1629N	4-0 24mm 45cm	\$84.50	\$74.50	\$69.50
CRAC0159N	3-0 18mm 45cm			
		Box/36	3+	5+
CRD662D	4-0 19mm 45cm	\$169	\$149	\$139

SharpPoint Nylon Black Monofilament		Box/36	3+	5+
CRAC0182D	4-0 18mm 45cm			
CRAC0568D	5-0 18mm 45cm	\$169	\$149	\$139
CRAC0586D	3-0 24mm 70cm			



LOOK®

- Extra sharp cutting edges extending from the point to its widest portion reduce tissue drag
- Precise, smooth penetration through tissue minimises trauma
- Needle channel is fitted to the specific size of the suture material

Nylon Sutures
Non-absorbable, monofilament nylon



Silk Sutures
Non-absorbable, coated/braided silk



Polysyn Sutures
Absorbable braided polyglycolin acid



		Box/12	3+	5+
LOO925B	3-0 24mm 45cm			
LOO922B	4-0 19mm 45cm	\$45	\$42	\$38
LOO918B	5-0 18mm 45cm			

		Box/12	3+	5+
LOO770B	2-0 24mm 45cm			
LOO784B	3-0 18mm 45cm	\$69	\$65	\$59
LOO781B	4-0 19mm 45cm			


		Box/12	3+	5+
LOO443B	2-0 24mm 70cm			
LOO423B	3-0 19mm 70cm	\$129	\$125	\$105
LOO421B	5-0 19mm 70cm			

ETHICON

The clinical advantage A full range of Ethicon Sutures is available

MONOCRYL*

Absorbable, monofilament poliglecaprone 25



		Box/12	2+
JJ77757	4-0 19mm 45cm	\$220	\$210
JJ77756	4-0 13mm 45cm	\$259	\$249
JJ77755	5-0 13mm 45cm	\$259	\$249

ETHILON*

Non-absorbable, monofilament nylon



		Box/12	2+
JJ76271	2-0 26mm 45cm	\$99	\$94
JJ76264	3-0 24mm 45cm	\$99	\$94
JJ76258	4-0 19mm 45cm	\$106	\$99
JJ76274	5-0 19mm 45cm	\$115	\$109

PDS* II

Absorbable, monofilament polydioxanone



		Box/12	2+
JJ78604	3-0 19mm 45cm	\$229	\$219
JJZ824G	3-0 19mm 45cm Conventional Cutting	\$229	\$219
JJ78606	4-0 16mm 45cm	\$248	\$239
JJ78602	5-0 13mm 45cm	\$265	\$259

PROLENE*

Non-absorbable, monofilament polypropylene



		Box/12	2+
JJ76444	3-0 24mm 45cm	\$135	\$129
JJ76440	4-0 19mm 45cm	\$149	\$142
JJ76427	5-0 19mm 45cm	\$158	\$149

VICRYL*

Absorbable braided monofilament Polyglactin 910



		Box/12	2+
JJ76486	3-0 24mm 70cm	\$381	\$359
JJ76470	4-0 19mm 70cm	\$435	\$415

VICRYL RAPIDE*

Absorbable braided Polyglactin 910



		Box/12	2+
JJ48583	3-0 24mm 70cm	\$655	\$645
JJ48013	4-0 19mm 45cm	\$665	\$655



DemeLON™ Nylon Sutures



DemeTech's Nylon Suture is a non-absorbable, sterile surgical suture composed of the long chain aliphatic polymers Nylon 6 and Nylon 6.6.

DemeTech's Nylon sutures are available in either black or blue. While DemeTech's Nylon is not absorbed, progressive hydrolysis may cause a gradual loss of up to 20% in its tensile strength each year. DemeTech's Nylon sutures are indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmology.

Features and Benefits

- Small USP sizes are designed specifically for microsurgical and ophthalmic procedures
- Made in U.S.A.

Properties

Tensile Strength	Material	Absorption Profile
Gradual loss of tensile strength over time	MONOFILAMENT Black/Colourless/Blue	Non Absorbable Suture

Description	UOM	Price	2+	3+	Description	UOM	Price	2+	3+
DTN270 2/0 75cm 24mm Reverse Cutting Blue	12				DTN034 4/0 75cm 24mm Reverse Cutting Blue	12			
DTN033 3/0 75cm 19mm Reverse Cutting Blue	12				DTN614 5/0 45cm 12mm Reverse Cutting Blue	12			
DTN035 3/0 75cm 24mm Reverse Cutting Blue	12	\$34.95	\$32.95	\$29.95	DTN031 5/0 75cm 19mm Reverse Cutting Blue	12	\$34.95	\$32.95	\$29.95
DTN032 4/0 75cm 19mm Reverse Cutting Blue	12				DTN613 6/0 45cm 12mm Reverse Cutting Blue	12			
DTN030 4/0 75cm 17mm Reverse Cutting Blue	12				DTN637 6/0 45cm 15mm Reverse Cutting Blue	12			

DemeLENE™ Polypropylene Sutures



DemeTech's Polypropylene suture is a non-absorbable, sterile surgical suture composed of a synthetic linear polyolefin. DemeTech's Polypropylene is blue to enhance visibility in the body.

Features and Benefits

- Flexible
- Allows for secure knot tying
- Smooth textile
- Made in U.S.A.

Properties

Tensile Strength	Material	Absorption Profile
Indefinite	MONOFILAMENT Blue	Non Absorbable Suture

Description	UOM	Price	2+	3+	Description	UOM	Price	2+	3+
DTP694 3/0 45cm 19mm Reverse Cutting Blue	12				DTP686 5/0 45cm 17mm Reverse Cutting Blue	12			
DTP699 3/0 45cm 24mm Reverse Cutting Blue	12				DTP692 5/0 45cm 19mm Reverse Cutting Blue	12	\$42.95	\$38.95	\$34.95
DTP397 3/0 75cm 24mm Reverse Cutting Blue	12	\$42.95	\$38.95	\$34.95	DTP392 5/0 75cm 19mm Reverse Cutting Blue	12			
DTP693 4/0 45cm 19mm Reverse Cutting Blue	12				DTP681 6/0 45cm 15mm Reverse Cutting Blue	12			
DTP396 4/0 75cm 24mm Reverse Cutting Blue	12								

DemeSILK™ Silk Sutures



DemeTech's Silk suture is a non-absorbable, sterile, surgical suture composed of an organic protein called fibrin.

DemeTech's silk sutures are processed to remove the natural waxes and gums from its origin.

As Silk is a natural fibre, it does elicit an acute inflammatory reaction, which is followed by gradual encapsulation of the suture by fibrous connective tissue (occurring usually in 14-21 days). DemeTech's silk sutures are indicated for use in general soft tissue approximation and/or ligation.

DeNon-absorbable

- Multifilament in nature
- Made in U.S.A.

Properties

Tensile Strength	Material	Absorption Profile
Gradual loss of tensile strength over time	MULTIFILAMENT BLACK WHITE	Non Absorbable Suture


Description	UOM	Price	2+	3+
DTS074 2/0 75cm 24mm Reverse Cutting Black	12			
DTS071 3/0 75cm 19mm Reverse Cutting Black	12			
DTS070 4/0 75cm 19mm Reverse Cutting Black	12	\$34.95	\$32.95	\$29.95
DTS072 4/0 75cm 24mm Reverse Cutting Black	12			
DTS083 5/0 75cm 17mm Reverse Cutting Black	12			



Choose Histoacryl® for Basic Skin Closure

175 YEARS

B|BRAUN
SHARING EXPERTISE



SAVE TIME

- Quick to apply – Dries in less than a minute
- No need for a dressing as the anti-bacterial barrier protects the wound
- No need for a second visit to the doctor to remove stitches – excellent cosmetic results
- Water-resistant – showering is possible

SAVE MONEY

- No need for refrigeration
- No need for local anesthesia
- No need for sterile needle-holder/instruments

		Box/5	2+			Box/10	3+
BN105005	Blue	\$135	\$129	BN1050044	Blue	\$210.00	\$205
BN1050060	Clear	\$135	\$129				

Switch to Histoacryl® for the repair of small wounds (less than 7cm in length)

GluStitch® Twist Wound Closure

Using GluStitch® Twist improves cosmetic outcomes

GluStitch® Twist is a sterile, topical tissue adhesive.

GluStitch® Twist is provided in a single-use applicator packaged in a protective foil pouch. It remains liquid until it is applied to the skin and spread into a thin layer. GluStitch® Twist is intended for the topical closure of surgical skin incisions and trauma-induced lacerations in areas of low skin tension that are simple, thoroughly-cleansed, and have easily approximated skin edges.

- Topical tissue adhesive
- 6 x 0.5ml Sterile Ampoules
- Simple twist-off cap
- Sets in seconds
- Violet colour to improve visibility
- Great value



		Box/6	3+	5+
GLU-6	Violet	\$119.90	\$115.50	\$108.90
GLU-6CL	Clear			

NOW AVAILABLE IN CLEAR

Leukosan® Adhesive

- Secure closure flexes with the skin
- No heat generation after application, very child friendly
- Easy controlled application



Box/10

\$230

SN7254100

BSN medical

Sequer Staple Remover

Sterile



Each

\$2.95

19220606EE

Skin Staple Remover



Each

\$17.54

3MSR-3

3M

LEUKOCLIP® SD

Skin Stapler with Exchangeable Cartridges

Indications: Wound closure following small, medium and large incisions in all operative disciplines.



BE47111	Leukoclip® SD Stapler Handle	Each	\$82.00
BE47114	Leukoclip® SD Staple Cartridge	Each	\$19.35
BE47113	Leukoclip® SD Staple Cartridge 35 Staples	Pkt/35	\$20.00

SmithNephew