

# INTELLIGENT PUBLIC ACCESS DEFIBRILLATOR

## **iPAD** CU-SP1

### CPR DETECTION


CPR is vital to ensure that the casualty has the best chances of survival.

The iPAD CU-SP1 detects if CPR is being performed and when appropriate prompts and encourages the first responder. If CPR is not being performed, the iPAD CU-SP1 prompts the first responder to 'perform CPR'. If CPR is already being performed, the iPAD CU-SP1 encourages the first responder by prompting them to 'continue CPR'.



### SMART PADS - SMART STORAGE!

The pre-connected 'Smart' electrode pads are stored in a clear storage compartment on the underside of the unit. This ensures that the iPAD SP1 is always ready to be used in the shortest amount of time possible. The unique technology built in to the iPAD CU-SP1 and the 'Smart' electrode pads allows the unit to detect, when connected, if a set of pads has already been used. If they have, the user is advised of this. In addition, the visual indicator on the front of the iPAD CU-SP1 indicates if the connected pads are within their life expectancy. The indicator will change when the pads have only 3 months life left before their expiry date-giving you plenty of time to arrange for replacements. The indicator will change again when the expiry date is reached. Checking the life of your pads is as quick and simple as looking at the LCD display.



### AUTO VOLUME ADJUST

Public places can often be noisy. Crowded stations, airports with constant announcements and streets filled with cars and people often make it difficult to hear what's going on around you. Trying to use an AED in noisy places can be difficult if the first responder can't clearly hear the voice prompts and instructions from the AED.

The iPAD CU-SP1 listens to the ambient noise, and automatically adjusts the volume of its prompts to make them clear and easily heard over the background noise. When the correct volume has been reached the iPAD CU-SP1 keeps the prompts at this level.



### PATIENT MODE SWITCH

- Easily switch from Adult to Child mode without changing pads.
- Safety cover prevents accidental switching





## IPAD CU-SP1

Automatic External Defibrillator



- Charging time: Less than 10 secs
- Charging time after CPR: 6 secs
- Output Energy Adult-150Joules /pediatric-50Joules
- CPR guidance – voice prompts for adult and child patient
- On/Off button & Shock button
- Indicator LCD display
- Automatic volume adjusting
- Status display: Device, battery Pads, sensing Pads expiring date, Pads connection status and CPR monitoring
- Patient analysis shockable rhythms – Ventricular Fibrillation, Ventricular Tachycardia
- Multi recording and storage capacity
- Includes: Defibrillator, Pads, Battery, Manual
- Wall cabinet and bracket options available

**\$1,895**

BIOCUSP1 GST Free

Replacement Battery

**\$365**

BIOCUI103



[www.biomedex.com.au](http://www.biomedex.com.au)

## Defibtech Lifeline PRO Defibrillator

Designed for the Professional



- Three operating modes available – AED, ECG monitoring and manual override
- High ratings for water and dust protection – IP55
- Manufactured in the USA
- Charge time less than 4 seconds
- Free software upgrades throughout 8 year warranty
- Daily automatic self-tests
- Full colour interactive display, ECG display, shock control, step-by step instructions for performing CPR and external defibrillation
- Compact and lightweight – 1.4kg
- 8 year warranty
- Long battery life with 125 shocks or 8 hours of continuous operation
- Rechargeable battery available

**\$4,295**

DCF-2410EN  
GST Free



[www.defibtech.com.au](http://www.defibtech.com.au)



## Defibtech Lifeline VIEW Defibrillator



- Award winning design
- Simple to use – First responder to untrained bystanders
- High ratings for water and dust protection – IP55
- Manufactured in the USA
- Charge time less than 4 seconds
- Free software upgrades throughout 8 year warranty
- Daily automatic self-tests
- Full colour interactive display showing step-by step instructions for performing CPR and external defibrillation
- 8 year warranty
- Compact and lightweight – 1.4kg
- Long battery life with 125 shocks or 8 hours of continuous operation

Semi-Automatic

**\$2,995**

DCF-E2310 GST Free

Replacement Battery

**\$299**

DCF-2003



## Defibtech Lifeline Defibrillator

Defibtech's breakthrough design makes it possible for every budget to afford a life-saving AED

- The original award winning Defibtech AED
- Longest battery life on the market – 7 years
- Available in Semi-automatic & Fully-automatic models
- Simple to use – First responder to untrained bystanders
- High ratings for water and dust protection – IP54
- Manufactured in the USA
- Free software upgrades throughout 8 year warranty
- Daily automatic self-tests
- 8 year warranty
- Lightweight and durable



DCF-E110	Semi- Automatic	<b>\$2,595</b>
DCF-E130	Fully-Automatic	<b>\$2,795</b>
DCF-E210	Replacement Battery	<b>\$435</b>

GST Free





## **IPAD** Defibrillator Pads

BIOCU1007	Defib Adult Pads for IPAD CU-SP1	Pk/2	<b>\$65</b>
BIOA1102S	Paediatric Pads for IPAD CU-SP1	Each	<b>\$65</b>
BIOCF10A05	Defib pads for IPAD	Pk/2	<b>\$85</b>
BIOECUER1A	Defib pads for CU-ER1/CU-ER2	Each	<b>\$85</b>
BIOECUER1P	Defib pads paed for CU-ER1/ER2	Each	<b>\$245</b>

## **defibtech** Defibrillator Pads

DDP-100	Adult Pads	Each	<b>\$105</b>
DDP-200P	Paediatric Pads	Each	<b>\$205</b>

All GST Free




**Also for:**  
*Adult Pads for Lifeline VIEW*  
*Adult Pads for Lifeline PRO*

## Defibrillator Pads

MT9130	Defib pads for Powerheart AED adult	Pk/2	<b>\$109</b>
LADP1	Defib pads for HeartStart FR2 adult	Pk/2	<b>\$89</b>

LAM3870A	Defib pad HeartStart FR2 paediatric	Pk/2	<b>\$205</b>
LAM5071A	Defib pads for HeartStart first aid adult	Pk/2	<b>\$125</b>



# DEFIBRILLATORS SAVE LIVES.

Make sure you have one installed in your practice, workplace, school or club.

Visit [www.vitalmedicalonline.com.au](http://www.vitalmedicalonline.com.au) to place your order now!