

HEINE quality diagnostic instruments

Performance you can rely on in every situation



Experience The Difference

- 5 year guarantee
- State of the art optics and illumination
- Wall, desk or portable configurations





HEINE EN200® Wall Unit With Beta 400 Otoscope and Beta 200 Ophthalmoscope

The HEINE EN 200 can be customised to match your individual needs.

- EN200® Wall Unit
- BETA 400 LED F.O. Otoscope
- BETA 200 LED Ophthalmoscope
- 1 set of reusable tips (4 pcs.)



HEX9512250	Wall Transformer EN 200	\$769
HEB0011149	Allspec Tip dispenser/White	\$61
HEX9512005	Wallboard/EN200/BP/Tip Sender	\$438
HEA9512210	EN200/180 F.O. Otoscope + 180 F.O. Ophthalmoscope	\$1,575



\$2,019

See page 66
for Sphyg
accessories for
the EN200

HEA9512208

BETA 400® F.O. LED Otoscope & Ophthalmoscope Sets

**Set complete with: BETA 400 FO LED Otoscope,
BETA 200 LED Ophthalmoscope**

- 1 set of reusable tips
- 10 disposable tips
- Hard case
- Beta 4 USB rechargeable handle
- USB cord and plug-in power supply

\$1,842

HEB5328388

**Set complete with: BETA 400 FO LED Otoscope,
Beta 200 LED Ophthalmoscope**

- 1 set of reusable tips
- 10 Disposable tips
- hard case
- BETA 4 NT rechargeable handle
and NT 4 table charger

\$2,175

HEB5324420



LED NOW IN HEINE QUALITY.

HEINE

K 180® Diagnostic Sets

- K 180 Ophthalmoscope
- K 180 F.O. Otoscope
- 1 set of reusable tips
- Auto on-off recharge feature



**K 180® 2.5 V
Diagnostic Set** with
standard BETA battery
handle

\$895

HEA2791011

**K 180® 3.5 V
Diagnostic Set**
with BETA 4 NT
rechargeable handle
and NT4 table
Charger

\$1,610

HEA7923420





HEINE QUALITY
MADE IN GERMANY

HEINE mini 3000 series Maximum quality, minimum size



- Handy and ergonomic
- 5 year guarantee
- Quality optics and illumination

Experience The Difference



mini 3000[®] LED Diagnostic Set

LED HO

- mini3000 LED Ophthalmoscope
- mini3000 LED F.O. Otoscope
- 2 mini 3000 battery handles
- 1 set of reusable tips
- 5 disposable ear tips both 4mm and 2.5mm

\$875

HED8611021

Hard case



mini 3000[®] Diagnostic Set

- mini 3000 Ophthalmoscope
- mini 3000 Otoscope
- 2 mini 3000 battery handles
- 1 set of reusable tips
- 10 disposable ear tips

Also available with
mini 3000 F.O.
Otoscope

\$679

HED8731102

\$585

HED8732102

Hard case



10% GST will apply to all items unless they are marked GST Free



3.5V Ni-Cad Diagnostic Set

High quality standard rechargeable set with Ni-Cad power source

Features:

- 3.5V Ni-Cad Diagnostic Set (with Diagnostic Otoscope and Coaxial Ophthalmoscope)
- 4 reuseable specula



\$795

WL97206-VC

WelchAlllyn®

Pocket LED Otoscope/ Ophthalmoscope Soft Case Set

Features:

- LED illumination
- SureColour LED
- Lightweight
- Cool light
- 7 year warranty on on SureColor™ LED Lamp
- 3 year warranty



Soft case

*Also available in Blueberry and Plum

WL92871BLK	Pocket LED Set – Onyx	\$275	WelchAlllyn®
WL92871WHT	Pocket LED Set – Vanilla		

3.5V NiCad MacroView™ and Coaxial Diagnostic Sets

Features:

- MacroView Otoscope provides nearly a complete view of the tympanic membrane with 30% greater magnification
- Coaxial Ophthalmoscope for a quality view of the fundus
- 4 reuseable specula
- 3.5V NiCad power source for longer use and better performance



\$795

WL97206MVC

WelchAlllyn®

Pocket LED Otoscope/ Ophthalmoscope Hard Case Set

Features:

- LED illumination
- SureColour LED
- Lightweight
- Cool light
- 7 year warranty on on SureColor™ LED Lamp
- 5 year warranty



Hard case

WL92880BLK	Pocket Plus LED Set – Onyx with hard case	\$375	WelchAlllyn®
WL92880WHT	Pocket Plus LED Set – Vanilla with hard case		

*Also available in Blueberry and Plum

See the bigger picture with Welch Allyn Frontline Diagnostic Tools

Superior clarity - sharper view

WelchAllyn® 3.5 V MacroView™ Otoscope Get a larger view. A better view. A sharper view.



Delivers almost twice the field of view at higher magnification compared to standard otoscopes, making it easier to visualise the auditory canal and tympanic membrane. This instrument is a significant advancement in hand-held otoscopy, providing an improved view for enhanced learning and decision making.

WelchAllyn® 3.5 V PanOptic™ Ophthalmoscope Revolutionary technology for a difference you can see



Patient Eye Cup: Establishes and maintains optimal viewing distance, provides image stability, serves as a pivoting point for panning around the retina, helps with orientation and screens out ambient light. The eye cup is cleanable and reusable

Soft Grip Handle: Ergonomically designed for maximum comfort, balance and access to controls

Patented Glare Extinguishment System: Prevents interference from unwanted glare and reflections

Dynamic Focusing Wheel: Continuous, smooth action for more precise control and optimum view. Focusing range is -20 to +20 dioptres

Aperture Dial: Allows you to select from micro, small, and large spot sizes, slit aperture, red-free and cobalt blue filters

The fundamental challenge in ophthalmoscopy is to get a good view of the fundus in order to make a sufficient assessment. Let the Welch Allyn PanOptic Ophthalmoscope help you overcome this challenge.



Connex® Integrated Wall System

The Connex® Integrated Wall system provides everything vital for improved patient safety, speeding patient flow, and saving valuable time and space.

Features:

- Eliminates frustrating, time-consuming searches for equipment
- Immediately view patient data, then send it directly to an EMR/ HIS/EHR system
- Prevents tools from being lost, stolen or misplaced
- Enhances the visibility of information to your patients, improving communication & compliance
- Eliminates tangled, hanging cords
- Provides accurate data every time, speeding decision-making
- Open architecture design for custom configuration & expansion
- Includes: Blood Pressure, SpO2, SureTemp® Plus, Braun Pro6000, Coaxial Ophthalmoscope, MacroView Otoscope, KleenSpec Ear Specula dispenser & cuff storage.

\$6,885

WL84MTVEC



WelchAllyn

GS777 Integrated System

With 11720 Ophthalmoscope, 23810 MacroView, 76700 Locking Collars and LED Lamps

- Includes 11720 Coaxial Ophthalmoscope head with LED bulb
- Includes 23810 MacroView Otoscope head with LED bulb
- Includes 76700 locking collars
- OptiSense sensor automatically turns handles on and off when removed and replaced from the cradle
- 5 year warranty, (10 years on coiled cords, 7 years on LED bulbs)

\$1,450

WL77716-71



WelchAllyn

GS777 Integrated System

With 11820 PanOptic Ophthalmoscope, 23810 MacroView, 76700 Locking Collars and LED Lamps

- Includes 11820 PanOptic Ophthalmoscope head with LED bulb
- Includes 23810 MacroView Otoscope head with LED bulb
- Includes 76700 locking collars
- OptiSense sensor automatically turns handles on and off when removed and replaced from the cradle
- 5 year warranty, (10 years on coiled cords, 7 years on LED bulbs)

\$1,895

WL77716-82



WelchAllyn

3.5V Lithium Ion PanOptic™ MacroView Diagnostic Set

Revolutionary Technology For A Difference You Can See

Features:

- PanOptic Ophthalmoscope
- MacroView Otoscope
- Reusable Specula
- Hard Zipper Case
- LED Lamps
- Lithium Ion Handle with Charging Stand and Charger

\$1,395

WL97206MVP



WelchAllyn



K 180® Oscope 3.5V Head Only

- Light, high-impact polycarbonate.
- Very bright illumination and an unobstructed view of the ear canal and tympanum.
- Shock-resistant.
- XHL Xenon Halogen Technology.
- Extremely bright, white light.
- Acrylic viewing window with 3 x magnification
- Swivelling viewing window built into instrument.



\$249

HEB0211550



mini 3000® Oscope

- Viewing window with 3 x magnification produces distortion-free images of the examined object.
- Permanently attached swiveling viewing window allows for instrument use, and cannot be misplaced or lost.
- Bright and homogeneous illumination with excellent colour rendering.



Direct XHL Xenon Halogen illumination
2x AA Batteries (Not included)

\$185

HED0170210



K 100® Diagnostic Oscope

K 100® 2.5 V Diagnostic Oscope with standard BETA battery handle, 3 reusable tips and Soft pouch

- Direct Halogen Illumination; extremely bright, white light
- Extensive choice of accessories; very versatile e.g. for nasal examination
- Insufflation port; reliable test of tympanic mobility
- Ideal for minor procedures



\$340

HEB2451011



mini 3000® F.O. Oscope

- Viewing window with 3 x magnification produces distortion-free images of the examined object.
- Permanently attached swiveling viewing window
- Distal Fiber Optic (F.O.) illumination ensures an unobstructed view of the tympanum and ear canal.
- Pneumatic testing of tympanic mobility through insufflation port.



Fiber Optic (F.O) illumination
2x AA Batteries (Not included)

\$279

HED0170110



mini 3000® LED F.O. Oscope

LED HO

- Twice as bright as conventional XHL – Xenon Halogen instruments.
- Excellent colour rendering, even for red colours.
- Maintenance-Free no need to ever exchange the LED
- Exclusive Battery Performance Indicator
- Operation Time up to 10 hours.



\$365

HED0870110



**Over 17,000
Products**

www.vitalmedicalonline.com.au



Pocket LED Otoscope with Handle and Soft Case

Features:

- LED illumination
- SureColour LED
- Lightweight
- Cool light
- 3 year warranty for the unit
- 7 year warranty on on SureColor™ LED Lamp



WL22880BLK	Pocket Plus LED Otoscope – Onyx with handle & soft case	\$190
WL22880WHT	Pocket Plus LED Otoscope – Vanilla with handle & soft case	

*Also available in Blueberry and Plum



Pocket LED Otoscope and Handle

Features:

- LED illumination
- SureColour LED
- Lightweight
- Cool light
- 3 year warranty for the unit
- 7 year warranty on on SureColor™ LED Lamp



WL22870BLK	Pocket LED Otoscope – Onyx with handle	\$130
WL22870WHT	Pocket LED Otoscope – Vanilla with handle	

*Also available in Blueberry and Plum



K 180® 2.5v F.O. Otoscope with standard BETA battery Handle

- 1 set of reusable tips
- Soft pouch



\$435

HEB1811011



Beta 400® F.O. 3.5V Otoscope Head Only

All-metal housing and scratch-resistant glass optics which guarantee a long instrument life and high-quality patient assessments.

- 4.2 x Magnification
- One hand operation
- Exceptional depth of field
- XHL - Xenon Halogen Bulb



NEW

\$385

HEB0211400



Beta 400® LED Otoscope F.O.

- Exclusive 4.2x magnification.
- Fibre optic illumination.
- Light intensity can be individually adjusted
- The permanently attached optical system which can be flipped open/swivelled out of the way
- All metal housing and a scratch resistant glass window



5 YEARS GUARANTEE

\$509

HEB0811400



SHORT ON TIME?

Try our online ordering system. Placing your orders online is quick and easy!

www.vitalmedicalonline.com.au



mini 3000® 2.5v Ophthalmoscope

- Dustproof construction.
- A wide range of apertures and filters.
- Quality HEINE optics for high resolution views of the retina.
- Brightness can be reduced by using the fixation star with grey filter.
- 2x AA Batteries (Not included)



HED0171120	Black mini 3000® Ophthalmoscope	\$329
HED171120B	Blue mini 3000® Ophthalmoscope	\$329



K180 3.5v Ophthalmoscope - Head Only

- XHL Xenon Halogen bulbs for accurate colour rendering clear illumination.
- Ergonomical shape to fit the orbita perfectly.
- Maintenance-free dust-resistant housing.
- Constructed of high-impact polycarbonate.
- Precision spherical HEINE optics produce excellent quality images of the entire area under examination.
- 5 apertures.



\$399

HEC0230207



mini 3000® LED Ophthalmoscope

- Optimum brightness for ophthalmic examinations.
- Excellent colour rendering, even for red colours
- Maintenance free no need to ever replace the LED
- Exclusive Battery Performance Indicator
- Operation Time up to 10 hours.

LED HQ



\$435

HED0871120



K 180® Ophthalmoscope with standard BETA battery handle

- Standard aperture wheel No.1
- Soft pouch



\$555

HEC1821011



Beta 200 3.5v Ophthalmoscope - Head Only

Exclusive dustproof design protects the precision optical components and ensures maintenance-free operation. The BETA200 / 200 S Ophthalmoscopes are the only instruments with an aluminium frame system. Optical components are flush-mounted on a cast aluminum frame, eliminating misalignment and making the instrument shock-resistant.

- Eliminates corneal reflex
- Ideal for small pupils



\$580

HEC0230100



BETA 200 LED Ophthalmoscope

LED HQ

- Exclusive dustproof design .
- Aluminium frame system.
- Optical components are flush-mounted
- Eliminates corneal reflex
- Ideal for small pupils



\$705

HEC0830100





Pocket LED Ophthalmoscope with Handle

Features:

- LED illumination
- SureColour LED
- Lightweight
- Cool light
- 7 year warranty on on SureColor™ LED Lamp
- 3 year warranty



WL12870BLK	Pocket LED Ophthalmoscope – Onyx with handle	\$170
WL12870WHT	Pocket LED Ophthalmoscope – Vanilla with handle	

*Also available in Blueberry and Plum



Pocket PlusLED Ophthalmoscope with Handle and Soft Case

Features:

- LED illumination
- SureColour LED
- Lightweight
- Cool light
- 7 year warranty on on SureColor™ LED Lamp
- 3 year warranty



WL12880BLK	Pocket Plus LED Ophthalmoscope – Onyx with handle & soft case	\$220
WL12880WHT	Pocket Plus LED Ophthalmoscope – Vanilla with handle & soft case	

*Also available in Blueberry and Plum



Download the FREE Vital Medical Supplies APP



VISIT PAGE 10
FOR MORE
INFORMATION

FREE download available from:





WelchAllyn®

Disposable Specula

- Economical and convenient
- Made from non-toxic plastic and are available in a variety of sizes

KleenSpec Disposable Ear Specula for Diagnostic Oscopes



WL52432-U	2.75 mm	Tube/34	\$3.95
WL52434-U	4.25 mm	Tube/34	\$3.95

Kleenspec Specula Dispenser


\$62

WL52400PF



HEINE AllSpec® Disposable Tips

For BETA 400, 200, K 180, mini 3000, mini 3000 F.O. Oscopes



				+20
HEB0011128	Ø 2.5mm infants	Pkt/50	\$7.50	\$7.25
HEB0011127	Ø 4mm adults	Pkt/50		

AllSpec Tip-Dispenser

- A clean, new tip for each patient
- Holds up to 200 specs
- Instant selection of child or adult tip



\$61

HEB0011145

UniSpec Specula Dispenser

A clean, new speculum for each and every patient



\$67

HEB0011146

WelchAllyn®

Reusable Otoscope Specula

- Reusable specula are lightweight yet durable, and may be autoclaved, boiled, or sterilised in solution
- SofSpec™ reusable specula also have an elastomer tip specially contoured for an ultrapneumatic seal



Diagnostic/MacroView/Standard Reusable Specula

WL24400	Set of Four – 2.5, 3, 4 & 5 mm	\$38.00
WL24305	5.0 mm	\$9.50
WL24304	4.0 mm	
WL24302	2.5 mm	

Pneumatic/Operating/Consulting Standard Reusable Specula

WL22002	2.5mm	\$9.50
WL22003	3.0mm	
WL22004	4.0 mm	
WL22005	5.0mm	
WL22009	9.0mm	

UniSpec Disposable Ear Specula

For K100 and BETA 100 Operating otoscopes; (non F.O otoscopes) (Wide base ear tips)



HEB0011242	Ø 2.5mm infants	Pkt/50	\$9.25
HEB0011241	Ø 4mm adults	Pkt/50	

Desktop Dispenser – Pack of 250 pieces (grey)

Easy to place on desk; or move from room to room



				+10
HEB0011151	Ø 2.5mm Infants	Pkt/250	\$42.00	\$39.00
HEB0011150	Ø 4mm adults	Pkt/250		

Insufflator Bulb

For Digital MacroView™ Otoscope and MacroView™ Otoscope

\$30

WL23804



Sofseal Sealing Tips

WL24320	Small 80s	\$56.00
WL24330	Medium 80s	

WelchAllyn®

Insufflation Bulb

For Pneumatic Testing



HEB0011240	Insufflation bulb for BETA 200 and K 180 Oscopes	\$18.60
HEB0080102	Insufflation bulb for mini 3000 F.O. Oscope with connector	\$25.95



HEINE mini NT[®] charger set for mini 3000

With this set, existing mini 3000 handles can be updated to the 2.5 V rechargeable system.

- Innovative, compact design. Practical and space saving.
- Two independent charging compartments. Complete functionality for one or two handles at a time. Two handles can be charged at a time
- Fast charge feature. 4 hours to recharge a mini NiMH 2.5 V rechargeable battery
- Automatic charging. Independent of the state of charge of the rechargeable battery. Guarantees extended battery life
- Cannot be overcharged. Charger switches off automatically when the battery is fully-charged
- Charge status display. Light pulses from the charging compartment confirm the charging mode
- No memory effect. Full power at any time without loss of capacity.
- Auto On/Off feature
- Includes two batteries and two base caps



\$470

HEX0199485



HEINE NT4[®] Table Charger

The best in performance and Quality

- Innovative, compact design. Practical and space-saving.
- Two independent charging compartments.
- Complete functionality for one or two handles at a time.
- Automatic charging. Independent of the state of charge of the rechargeable battery.
- Guarantees extended battery life.
- Cannot be overcharged. Charger switches off automatically when the battery is fully charged.
- Charge status display. Light pulses from the charging compartment to confirm the charging mode.



\$485

HEX0299494



Universal Desk Charger with 2 Lithium Ion Handles

- Colored LEDs indicate power to the unit and that handles are charging
- Constant automatic safe charging and optimum lamp performance and life
- Accepts 3.5 V or 2.5 V NiCad, Lithium Ion and PocketScope™ handles
- Compact for convenience and mobility

WelchAllyn

\$1095

WL71736



HEINE Rechargeable Batteries



Used for 3.5V BETA rechargeable handles with NT300 or BETA TR chargers.

HEX0299382	NiMH rechargeable battery 3.5V	\$137
HEX0799383	Li-ion L rechargeable battery 3.5V Li-ion L battery only works in conjunction with BETA L bottom insert	\$155



Welch Allyn Rechargeable Batteries



WL72200	3.5 V NiCad Battery: Black for 71670 NiCad Handles and 71020-A Handles	\$105
WL72300	3.5 V NiCad Battery: Orange for 71000-C and AudioScope	\$105
WL71960	3.5 V Lithium-Ion Battery	\$299

WelchAllyn



BETA Battery Handle

- Standard 2.5V handle
- Uses 2 C-size batteries (not included)



Upgrading a HEINE handle from 2.5V to 3.5V

Heine 2.5V handles can be upgraded to a 3.5V rechargeable system anytime by purchasing the following accessories:

HEX0299494	BETA 4 NT4 Table Charger	\$485
HEX0299394	BETA 4 NT Bottom Insert	\$109
HEX0799383	Li-ion Rechargeable Battery	\$155
HEX0299395	BETA 4 USB Bottom Insert	\$145
HEX0099303	USB Cord	\$87
HEX0099305	E4 USB Medical Approved Plug-In Power Supply	\$49

\$137

HEX0199118



BETA Handle Belt Clip



\$52

HEX0099008



Replacement Globes



WL00200	\$44		WL04200	\$139		WL06500LED	\$144	
WL01300	\$16		WL04400	\$80		WL07600	\$35	
WL03000	\$67		WL04430	\$51		WL07800	\$64	
WL03100	\$41		WL04900	\$51		WL08500	\$60	
WL03300	\$67		WL06000	\$76		WL09500	\$780	
WL03400	\$67		WL06300	\$51		WL03100LED	\$139	
WL03800	\$51		WL06400	\$131		WL04900LED	\$145	
WL03900	\$75		WL06500	\$41		WL03800LED	\$145	

Replacement Globes - XHL® Xenon Halogen Bulbs



Otosscopes			Ophthalmoscopes			Dermatoscopes					
HEX0188037	Mini 2000 F.O	\$32		HEX0188084	K180 2.5 V	\$45		HEX0188109	Mini 3000	\$42	
HEX0188041	Mini Fibrilux	\$32		HEX0288070	Beta 200 3.5 V	\$45		HEX0288046	Delta 10 3.5 V	\$45	
HEX0188056	Mini 2000	\$32		HEX0188042	Mini 2000/Alpha	\$46		HEX0188034	Mini 2000	\$45	
HEX0188105	Mini 3000 F.O.	\$40		HEX0188069	Beta 200 2.5 V	\$46		HEX0188038	Delta 10 2.5 V	\$45	
HEX0188110	Mini 3000	\$40		HEX0188106	Mini 3000	\$48					
HEX0188077	K180 2.5 V	\$42		HEX0288086	K180 3.5 V	\$50					
HEX0288078	K180 3.5 V	\$50									



RETeval-DR™ Electroretinograph

Welch Allyn's RETeval-DR Electroretinograph is designed for healthcare providers who want to improve management of their patients with diabetes.

The RETeval-DR objectively screens for vision-threatening diabetic retinopathy with a 99.2% negative predictive value in just 3 minutes

- Non-mydriatic for patient comfort
- Quickly screen for vision-threatening diabetic retinopathy wherever patients with diabetes receive their care
- Time savings for the clinician and patient
- No advanced clinical training is required to perform the test or read the results. ReTeval-DR provides immediate test results and report
- High accuracy and low technical failure rate
- A lower cost option compared to dilated eye exams and teleretinal imaging services
- Handheld and battery powered



\$7,680

WLRETEVAL

WelchAllyn®

Spot™ Vision Screener

Welch Allyn's Spot™ Vision Screener is a handheld, portable device designed to help users quickly and accurately detect vision issues on patients from 6 months of age through adults.

- Fast binocular screening with a one-second capture time
- Ease of use – minimal user training is required
- A noninvasive 3-foot working distance requiring minimal patient cooperation
- Automated screening provides thorough, objective results in seconds
- Easy-to-understand results help determine follow-up care
- Successfully captures readings 97% of the time
- Reliable, automated screenings mean that you're more accurately referring patients who need to see eye care specialists for more thorough testing than if you were using another means for screening



\$8,595

WLVS100S

WelchAllyn®