



## Projet 101 Pulsatile Ear Irrigator

Designed by Doctors with innovative pulsatile irrigation

Utilising the most modern production methods the Guardian Projet 101 was designed to fulfil the most demanding surgeries and hospital ward requirements for a reliable, leak free and easy to use machine that is not expensive to run and does not require constant maintenance and subsequent down time.

- Extra long mains lead and handle hose so machine can be sited where it is most convenient
- Integral wall mount frees up valuable work space by wall mounting the Projet™ however easily removed for cleaning etc
- Sliding flow regulator conveniently mounted in handle. Water flow can be temporarily stopped by sliding regulator to off position or flow regulated by sliding control to desired position
- Ultra thin atraumatic tips help visualisation and direction of water flow
- Digital control of water pressure and flow is easy to operate and remains constant and reproducible

# \$379

PROJET101

**NEW**



PROJET-TIP Tips Box/100 **\$135.00**



## Ear Toilet Package

Microscope, Suction Unit, Sim Micro Cannula, Ear Gruber Set

- 1500mm Super long pantograph arm offers a wide operation space
- Four different working distances: 200mm, 250mm, 300mm, 400mm
- Smooth turret 3-step magnification pod
- Field diameter: 7.88mm-81mm
- Magnification: 2.4X-12X
- Green and blue filter selection
- 10W LED illumination
- 20,000 hours LED lifespan
- 50,000 Lux object illumination
- Power supply voltage: 90V/60HZ-250V/50HZ
- Includes 12 Months Parts Warranty

**Massive Medicare Rebate - 41647**

# \$5,995

ZOEAPACK

Category 3 - THERAPEUTIC PROCEDURES

41647	Group T	8-SURGICAL OPERATIONS
	Subgroup	8 -EAR, NOSE AND THROAT

EAR TOILET requiring use of operating microscope and microinspection of tympanic membrane with or without general anaesthesia.

**Fee:** \$109.00 **Benefit:** 75% = \$82.45 85% = \$93.45





## Ear Syringe NJ Phillips Precision Instrument

[www.njphillips.com.au](http://www.njphillips.com.au)

For use in the removal of waxy build-up in the ear

- Ergonomic design
- Extremely accurate
- Robust, reliable and durable
- Easily adjusted
- Calibrated from 0.5ml to 5ml
- Silicone tips

**\$275**

MCS1043



### ASSEMBLING METAL TIP UNIT:

1. Loosen needle nut until reasonably loose.
2. Slide needle tip base into the lower slot of the needle nut.
3. Screw the needle nut tightly back onto the main threaded mount.

### ASSEMBLING RUBBER TIP UNIT:

4. The rubber tip **MUST** always be sterilised prior to use.
5. Remove the metal tip and needle nut.
6. Push the rubber tip onto the main threaded mount.

Turn clockwise as you push to ensure it is firmly in place.

### Spares

MCNITIS020	Nozzle only		<b>\$39.95</b>
MCNIWX1187	Repair kit		<b>\$39.95</b>
MCS137	Rubber tips	Pkt/25	<b>\$49.00</b>

## Ear Syringe With O Ring

ParagonCare

- Ring fitted on plunger giving smooth movement
- Unit supplied with appropriate nozzles (2x)



### Spares

LRESORT	O-ring for nozzles	<b>\$6.75</b>
LRLESOR	O-ring for plunger	<b>\$7.95</b>
SAS018	Nozzle Set (Pair)	<b>\$10.50</b>
SAS001	Shield	<b>\$3.50</b>

**\$75**

LRLESWOR

## Ear Syringe

With Dual O Rings,  
2 Nozzles & Shield



Dual Tips

**\$52**

SY2010

SAS010R	Rubber O ring for tip	<b>\$3.75</b>
---------	-----------------------	---------------

## Ear Syringe

Metal Premium German Quality



**\$165**

HTESYR



## Propulse® Irrigator Accessories

KIT0099	Propulse Cleaning Tablets	Box/200	<b>\$45.00</b>
KIT0021	Propulse Noots Tank	Each	<b>\$39.00</b>
KIT0028	Propulse Ear Wax Removal Scoop	Box/50	<b>\$115.00</b>
KIT1012	Propulse Protect Cape	Box/36	<b>\$37.50</b>
KIT5002B	Propulse QrX™ Disposable Tips	Box/100	<b>\$195.00</b>

**mira3e**  
Health Group



You can place your order  
via fax on **1300 557 631**.