

*Stylish treatment
couches for all
practice applications*



ECLIPSE Electric Examination Couch

Features

- Electronic Hi Low
- Retractable 75mm castors
- Gas assisted head and foot sections
- High quality foam and upholstery
- Foot control
- Safe working load (SWL) 220Kg



2 Section
Model: EE220



3 Section
Model: EE320

ACERO Electric Examination Couch

Features

- Electric operated variable height
- Electric operated head section
- Electric operated foot section
- Padded, upholstered top, with face cut out
- 100mm castors; fully lockable
- Safe working load (SWL) 200kg



2 Section
Model: AX46002SECT



3 Section
Model: AX46103SECT






Acero All Electric Examination Couch

ParagonCare

Acero examination couches offer a very cost effective solution to any consultation facility. The all electric operation not only provides ease of use, but is ergonomic and OH&S friendly. The couches are fitted with 100mm castors, fully braked which enables transportation where necessary. Twin support frame design enables excellent stability and weight distribution. Unique base design gives clear access for patient lifters and over bed tables.

2 Section

\$1,795



AX4600


- Electric operated variable height
- Electric operated head section
- Padded, upholstered top, with face cut out
- 100mm castors; fully lockable
- Safe working load: 200kg
- Height Range: 490 – 900mm
- Head Section: 595mm, Positive tilt to 85°
- Foot Section: 1300mm
- Size: 1900mmL x 700mmW
- 3 year warranty

Optional:

- Foot Control
- Side Rails



3 Section



\$1,895

AX4610

- Electric operated variable height
- Electric operated head section
- Electric operated foot section
- Padded, upholstered top, with face cut out
- 100mm castors; fully lockable
- Safe working load 200kg
- Height Range: 490 – 900mm
- Head Section: 595mm, Positive tilt to 85°
- Foot Section: 720mm
- Mid Section: 530mm
- Size: 1900mmL x 700mmW
- 3 year warranty

Optional:

- Foot Control
- Side Rails



AX5570	Pair of Side Rails	\$309.00
AX3021	Foot control for Acero 2 Section	\$165.00
AX3017	Foot Control for Acero 3 Section	\$229.00

Freight charges will apply



Eclipse Electric Examination Couch - 2 Section

ParagonCare



- Electronic Hi Low
 - Retractable 75mm castors
 - Gas assisted head and foot sections
 - Breath holes on back rest and foot rest (where applicable)
 - High quality foam and upholstery (60mm deep)
 - Foot control (standard)
 - Safe Working Load (S.W.L) 220Kg
- Overall Size: 1950 x 700mm
 - Back Rest: 700 x 700mm
 - Mid Section 1250 x 700mm
 - Height Range: 460 - 935mm

\$1,195

AX220

Freight charges will apply

Eclipse Electric Examination Couch - 3 Section

ParagonCare




- Electronic Hi Low
 - Retractable 75mm castors
 - Gas assisted head and foot sections
 - Breath holes on back rest and foot rest (where applicable)
 - High quality foam and upholstery (60mm deep)
 - Foot control (standard)
 - Safe Working Load (S.W.L) 220Kg
- Overall Size: 1950 x 700mm
 - Back Rest: 800 x 700mm
 - Mid Section 480 x 700mm
 - Foot Rest: 650 x 700mm
 - Height Range: 460 - 935mm

\$1,250

AX320

Freight charges will apply



New Facility or upgrading your practice?

We can help you make the right decisions and bundle a special equipment package for you at a great price!

Call 1300 557 651



Altona 3 Section Couch

3 Section Variable Height Examination Couch

The Medical 3 Section has an adjustable back rest, with the added versatility of a gas strut operated foot section. This allows more versatility to the clinician. Standard width is 70cm, and this model comes with no recess.

Specifications	
Colour	Balmain Blue
Min/Max Height	46cm/91cm
Cushion Width	70cm
Head/Mid/Leg Section	68.5cm/56.5cm/68cm
Total Length	195.5cm
Head Tilt Angle/Leg Tilt Angle	-35 to 87°
Lift Mechanism	Electric
Safe Working Load	225kg
5 Year Frame warranty	
3 Year Electronic warranty	

\$1,395

AC-001
Freight charges will apply


AC-025 Paper Towel Holder **\$59.00**



Preston 3 Section Couch

3 Section Variable Height Examination Couch

- Equipped with one motorised actuator for high/low adjustment
- Individual "levelling foot" on all four corners
- Four retractable heavy duty castors
- Twin gas struts for solid head section support
- Classic dual lift construction through all ranges of height adjustment
- 62cm wide
- Blue only
- 225 kg SWL



\$1,395

ACPRESTON

Freight charges will apply




Regal Bariatric 3 Section Couch

3 Section Variable Height Examination Couch

An excellent solution for supporting bariatric patients in both sitting and lying positions, a must for all Hospital Departments, Health Centres and Care Centres.

- Full electric operation
- Foot switch controls with up and down functions
- Hand switch controls with back and leg functions
- 320 kg SWL



\$2,195

AC-004

Freight charges will apply





Seers Bariatric Couches

ParagonCare

These products meet the increasing requirements for providing examination couches with the capability of supporting extremely heavy patients. All these models have a safe working load and lifting capacity, from minimum height of 325kg.



Bariatric 3 Section - 850mm wide

\$4,095

AX3573

Bariatric Model Specifications

- 325kgs SWL (from minimum height)
- Adjustable angle electric backrest from horizontal to +80° positive on all models
- Electric height elevation from 43-97cm and 45-98cm on large wheel model
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on all models
- 3 year warranty



Bariatric 2 Section - 1050mm wide

\$4,895

AX4642



Freight charges will apply


Seers Couch Range

ParagonCare

Gynaecology Couch

Designed for general gynaecology procedures featuring a fully electric height adjustment. This unit features a crescent shaped seat design and are supplied with fully adjustable leg supports, paper towel holder and tray.

- 240kg SWL
- Back rest adjustable from horizontal to 80 degrees
- Crescent shaped seat design
- Fully adjustable leg supports
- Side support loops
- Tray & paper towel holder
- Height range 50-102cm
- Compact frame design with retractable wheels
- Adjustable foot for stability on uneven surfaces
- 3 year warranty
- Optional Seat Extension: AX6043



\$3,495


AX8563



Multi Couch

Featuring an adjustable angle back rest and with patient arm rests and dual adjustable leg rests. This unit is suitable for a wide range of seated examination procedures such as podiatry, leg ulcer clinics and orthopaedic casting rooms. Ideal for GP surgeries and Health centres.

- 240kg SWL
- Back rest adjustable from horizontal to 80 degrees - electric operation
- Adjustable angle leg rests are adjustable from horizontal to -90 degrees and incorporate anti-trap design gas struts
- The fold down armrests have angular adjustment
- Paper towel holder
- Height range 50-105cm - electric operation
- Compact frame design with retractable wheels
- Adjustable foot for stability on uneven surfaces
- 3 year warranty



\$3,195

AX0565



Freight charges will apply



Gynaecology Table

With Leg Extensions (Stirrups) and Foot Pads

- Operation by dual foot switch.
- Electric Hi Lo with the other motor for adjustable seat angle
- Sleek design with rounded frame legs for safety
- Minimum height to maximum height in 30 seconds
- Leg plug for examination
- Gas struts control the back section
- Adjustable angle backrest from horizontal to +80 degrees
- Adjustable angle seat from horizontal to +30 degrees
- Excellent height range of 60cm – 100cm
- Standard width of 52cm
- Fitted with stainless debris tray
- 7.5cm castors for safety and mobility
- Retractable wheel system for ease of mobility
- T Motion motor
- 250 kg SWL

\$2,795

PMGCNB

Freight charges will apply



Podiatry Multipurpose Chair

- Electric Hi Lo operation by foot switches
- Sleek design with rounded frame legs for safety
- Minimum height 55cm to maximum height 90cm in 30 seconds
- Dual motors: Electric Hi Lo. 2. Electric mid-section
- Adjustable backrest featuring double gas struts, angle adjustable from 0o to +60o
- Adjustable seat section featuring one motor, angle adjustable from 0o to +20o.
- Adjustable foot section featuring single gas strut
- Double removable arm supports
- Able to lift up to 250kgs working weight from its lowest height.
- Retractable castors 7.5cm for excellent mobility and safety.
- Table goes horizontal for people recovery in emergency situations.

\$2,495

PMPCNB

Freight charges will apply



Electric Contoured Massage Couch

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Minimum height 50cm to maximum height 90cm in 30 seconds
- Adjustable head section with manual operation, angle from +45° to -45°
- Arm rest and cradle
- Adjustable foot for uneven surfaces
- 5cm castors for safety and mobility
- Retractable locking castors
- T Motion motor
- Padding thickness 5cm
- Foam Density: 36g/m2
- 250 kg SWL

\$1,750

PMCMNB

Freight charges will apply





Examination Couches

ParagonCare

Standard

- 1810L x 600W x 820H cm
- Adjustable head section
- Detachable legs
- Powder coated frame
- Heavy-duty vinyl upholstery
- Comes flat packed
- 125 kg SWL

\$375

AX357G



Heavy Duty

- Designed for greater weight capacity 250kg SWL
- Wider patient surface
- Comfortable thick padded upholstery
- 1810L x 750W x 820H mm
- Dual strut back rest system
- Powder coated fully welded metal frame
- Dark grey upholstered top
- Optional modular drawer system available (see below)

\$659

AXEXAM0385



Optional modular drawer system

AXEXAM0382



Freight charges will apply

Free Standing Couch

No face hole

- Fixed height 82cm
- New easy lift and lock mechanism for easy assembly
- Lift up back section for assessment
- Comes flat packed
- Total Length: 195cm
- Width: 71cm
- Height: 82cm
- Loading weight: 280Kgs

\$385

PMFSGR



Freight charges will apply



Brooksbuilt Cabinets

All cabinets are made from Australian produced material. All the functioning parts are made in the traditional cabinet making method to ensure continued reliability for the life of the couch. A variety of colours are available.



Dimensions: 178L x 56W x 85H cm






Photo: Teak with dark brown top

www.brooksbuiltcabinets.com.au

Cupboard Couch		
18000225	Plain	\$1,495.00
18000222	Wood	\$1,525.00
Half Multibank Couch		
18000445	Plain	\$1,720.00
18000442	Wood	\$1,750.00
Multibank Couch		
18000665	Plain	\$1,875.00
18000662	Wood	\$1,895.00

Upholstery Colours	Black, Light Grey, Dark Grey, Dark Brown, Blue, Burgundy, Green, Beige, Melon
Cabinet Colours Plain	White, Light Grey, Dark Grey
Cabinet Colours Wood Grain	Ash, Oak, Beech, Teak
Pricing depends on the wood colour you choose (plain or wood grain)	

Couches made to order; delivery approx 6-8 weeks

Freight charges will apply

Blood Collection Chair

ParagonCare

This comfortable chair is widely used in pathology centres and may be used in a variety of other applications such as minor procedures and general examinations.

- 80L x 75W x 105H cm
- Seat width 46 cm
- Powder coated frame
- Upholstered arms and seat
- Non-slip step surface



AX468	Grey	\$575
AX468B	Black	

Freight charges will apply

Blood Collection Chair

With Foldable Foot Rest

- Comfortable, soft and supportive high density foam
- Seat width 55cm
- Width including arm rests 91.5cm
- Arm rests width 18.5cm
- Length 65cm
- Length including foot rests 84cm
- Foldable foot rest
- Weight capacity 175kg



\$525

PMBCC

Freight charges will apply

Portable Massage Table

- Aluminium frame and legs
- Manufactured without electric welding to join the aluminum
- Auto lock cabling system
- Rounded corners and deluxe no voice handle
- Luxurious PU vinyl
- Multi density 5cm foam system



\$325

PMPMTNB

Freight charges will apply



You can place your order via fax on
1300 557 631.

HEALTH. COVERED.

ParagonCare

paragoncare.com.au

Your Partner in Healthcare

With an extensive range that combines Australian design, local manufacturing and market leading global brands, Paragon Care have your healthcare solutions covered.



Paragon Care continues to build it's strong representation within the following healthcare markets:

- Hospital & Acute Care
- Aged Care
- Primary Care
- Ophthalmology / Optometry
- Materials Handling
- e-Health
- Service & technology

Paragon Care product range includes:

- Diagnostic Equipment
- Doctors Bags
- Surgical Instruments
- Carts & Trolleys
- Examination Couches
- Privacy Screen Systems
- Stools, Steps, & Skips
- IV Stands & Oxygen Bottle Trolleys
- Wheelchairs
- Stainless Steel Ware
- General Medical Products



ParagonCare

Spectra Mobile Storage Carts

The Spectra range is a cart system that combines mobility and drawer flexibility into a lightweight colour coded package which has applications throughout all healthcare environments.

The unique drawer system is based on removable 60 x 40 tray systems that offer different depths and the option of integrated drawer dividers. This allows flexibility with the carts applications, particularly in the ability to restock the carts with transfer tray systems.

The carts shelves are one piece moulded polymer which are easy to clean and eliminate infection control issues due to the rounded internal mould design, whilst integrating a raised lip 3 side to prevent items falling off.

The Spectra range offers unlimited drawer configurations and a range of colourful drawer fronts.

Multi Drawer Trolley

700 x 500 x 1000 mm
AX-M0605 - 1 x 50mm deep tray
AX-M0601 - 2 x 100mm deep trays
AX-M0602 - 2 x 200mm deep trays

Interlocking dividers optional
Work surface area 620 x 460mm
125mm castors

2 Drawer Trolley

700 x 500 x 1000mm
AX-M0601 - 2 x 100mm deep trays
Interlocking dividers optional
Work surface area 620 x 460mm
125mm castors

1 Drawer Trolley

700 x 500 x 1000mm
AX-M0601 - 1 x 100mm deep tray
Interlocking dividers optional
Work surface area 620 x 460mm
125mm castors

2 Shelf Trolley

700 x 500 x 1000mm
Work surface area 620 x 460mm
125mm castors



AXH0030R	Red	\$1,995
AXH0030B	Blue	

AXH0020R	Red	\$1,150
AXH0020B	Blue	

AXH0010R	Red	\$795
AXH0010B	Blue	

AXH0001	\$519
---------	--------------

Freight charges will apply

Spectra Mobile Storage Cart Trays

*Purchased separately. Call us to order



50mm deep tray*

AX0605
Interlocking divider 600mm
Interlocking divider 400mm



100mm deep tray*

AX0601
Interlocking divider 600mm
Interlocking divider 400mm



200mm deep tray*

AX0602
Interlocking divider 600mm
Interlocking divider 400mm



Treatment Trolleys All models supplied with a basket, waste bin and a push handle

Freight charges will apply

1 Drawer + Accessories



650 x 480 x 900mm

\$875

PMTTR01

2 Drawer + Accessories



650 x 480 x 900mm

\$960

PMTTR02

3 Drawer + Accessories



650 x 480 x 900mm

\$1040

PMTTR03





Instrument Trolleys

Paragon Care instrument trolley design methods are class leaders in providing high quality carts for the surgical, clinical and consulting environments.


Stainless Steel Series



- Stainless steel construction
- Safety folds to all edges
- Reinforced shelves
- 75mm castors
- Stainless steel drawers
- Full extension precision bearings drawer runners

ParagonCare


**Guard Rail
No Drawer**



AX033	490 x 490 x 900mm	\$445.00
AX036	600 x 490 x 900mm	\$495.00
AX039	750 x 490 x 900mm	\$535.00
AX042	900 x 490 x 900mm	\$575.00

Freight charges will apply


NOVA Series



- Powdered coated frame
- Stainless steel shelves and drawer fronts
- 50mm castors

ParagonCare

**Guard Rail
No Drawer**




490 x 490 x 900mm

\$345

AX321

Freight charges will apply


Powder Coated Series



- Powder coated frame and shelves
- Safety folds to all edges
- Reinforced shelves
- 75mm castors

ParagonCare

**Guard Rail
No Drawer**



490 x 490 x 900mm

\$369

AX342


Freight charges will apply



ParagonCare

ParagonCare

1 Drawer




Freight charges will apply

AX054	490 x 490 x 900mm	\$575.00
AX057	600 x 490 x 900mm	\$659.00
AX060	750 x 490 x 900mm	\$750.00
AX063	900 x 490 x 900mm	\$815.00

ParagonCare

2 Drawer



Freight charges will apply

AX066	490 x 490 x 900mm	\$825.00
AX069	600 x 490 x 900mm	\$895.00
AX072	750 x 490 x 900mm	\$1,015.00
AX075	900 x 490 x 900mm	\$1,085.00

ParagonCare

3 Drawer



490 x 490 x 900mm
with 125mm castors

\$1,150

AX076

Freight charges will apply


ParagonCare

1 Drawer

490 x 490 x 900mm

\$425

AX345



Freight charges will apply


ParagonCare

2 Drawer

490 x 490 x 900mm

\$485

AX348



Freight charges will apply


ParagonCare

3 Drawer

490 x 490 x 900mm

\$545

AX349



Freight charges will apply


ParagonCare

1 Drawer

490 x 490 x 900mm

\$409

AX324



Freight charges will apply


ParagonCare

1 Drawer

750 x 490 x 900mm

\$475

AX333



Freight charges will apply


ParagonCare

2 Drawer

490 x 490 x 900mm

\$529

AXH0327



Freight charges will apply



Multi Drawer Trolley

ParagonCare

- Cabinet Size: 580 x 480 x 980mm
- Overall: 650 x 530 x 985mm
- Soft closing drawers: 2 x 80mm, 2 x 160mm, 1 x 240mm
- Sliding extension platform
- ABS Top surface
- Corner protectors
- Central locking system
- Castors with brakes

NEW



\$1,450

AXCAM201

Freight charges will apply

Anaesthesia Cart

ParagonCare

- Cabinet Size: 580 x 480 x 980mm
- Overall: 650 x 530 x 985mm
- Soft closing drawers: 2 x 80mm, 2 x 160mm, 1 x 240mm
- Multi function top shelf with eight tilt bins
- Sliding extension platform
- ABS Top surface
- Corner protectors
- Central locking system
- Castors with brakes

NEW



\$1,525

AXCAA005

Freight charges will apply

5 Drawer Emergency Trolley

ParagonCare

- Construction of sturdy stainless steel frame
- Durable stainless steel top shelf rails on 3 sides
- Colour coordinated powder coated panels & drawers
- 3 x drawers 125mm deep
- 2 x drawers 250mm deep
- Precision bearing drawer runners, full extension
- Heavy duty 125mm castors
- Corner buffers 125mm diameter
- 2 year warranty
- Colour options: flame red or blue

Accessories

AX0120	Drawer insert container	\$119.00
AX0126	Drop down work surface	\$155.00
AX0142	Waste Bin	\$115.00
AX0123	Vertical Oxy bottle holder and clamp	\$135.00
AX0129	IV pole and clamp	\$135.00
AX0135	Power board 4 outlets	\$105.00
AX0143	Push Handle	\$69.00



\$1,495

Red: AX119R
Blue: AX119B

Freight charges will apply

Emergency Trolley

- Defibrillator shelf
- Extension table
- 304 Stainless steel IV pole
- Cardiac CPR board
- Waste bin
- External power holder & plug board
- Oxygen tank holder
- 304 Stainless steel wire basket
- Central or seals lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawers
- Dimensions: 650 x 480 x 980mm
- Drawer Sizes: 3 x 110mm, 1 x 135mm, 1 x 195mm
- G.W 55kg
- Central locking system through side



ALL ACCESSORIES INCLUDED

Freight charges will apply

\$1,895

PMEMTR



5 Drawer Emergency Trolley - All accessories included

ParagonCare

- Dimensions: 650 x 480 x 980
- Drawer Sizes: 2 x 110mm, 2 x 135mm, 1 x 195mm
- Soft closing drawers
- Lever action central locking
- Castors with brakes
- Defibrillator shelf
- Extension table
- Oxygen tank holder
- Waste bin
- Cardiac CPR board
- Non slip mat on top surface
- Push handle
- Wire basket



\$1,825

AXCAE005

Freight charges will apply



Chrome Wire Shelving

ParagonCare

Flexible shelving system with the ability to add or reduce the amount of shelves required.

- All shelves can be adjustable and positioned to suit storage requirements
- High quality chrome plated finish
- Adjustable foot for uneven surfaces
- Flat packed and easy to assemble
- Standard size 1524x533x1750mm
- Many other sizes available



\$545

AX940

Freight charges will apply

Trolley with Wire Baskets

ParagonCare

A versatile solution for mobile storage. This unit comes standard with a solid stainless steel top, baskets are completely supported with high quality full extension drawer runners.

- Size 610x460x1020mm high
- 4x100mm deep baskets
- 1x200mm baskets
- Different basket configurations are available



\$1,095

AX2701

Freight charges will apply

Chrome Wire Trolley

ParagonCare

- Size 460 x 760 x 1020mm high
- Chrome plated shelves
- Shelves adjustable if required
- 125mm castors and corner buffers
- Built in push handles

\$379

AX2350

Freight charges will apply



Wall Panel Baskets


ParagonCare

This basket is designed to replace the plastic tray system to eliminate dust collecting on the base of the unit.

- Close mesh design allows storage of small items
- Light weight robust design
- Replaces plastic tray system

AXWB134	105 x 125 x 85mm	\$14.00
AXWB132	105 x 190 x 85mm	\$17.00
AXWB122	415 x 350 x 200mm	\$42.00

Metal Wall Panels		
AXWP460M	460 x 460mm	\$67.00
AXWP900M	460 x 900mm	\$140.00



Freight charges will apply

Tilt Bins

ParagonCare

- Clear removable bins keep your supplies visible.
- The bins tilt out effortlessly
- Mounting hardware and labels included
- Mix and match sizes for your own requirements



AXH20303	3 Bin Unit	60W x 24H x 20D cm	\$145.00
AX20403	4 Bin Unit	60W x 21H x 17D cm	\$98.00
AX20503	5 Bin Unit	60W x 16.5H x 13.5D cm	\$95.00
AX20603	6 Bin Unit	60W x 11.5H x 9.5D cm	\$54.00



Single Linen Skip

ParagonCare

- Frame designed to support bag
- Easy removal of bag
- Sturdy frame construction
- Heavy duty 75mm castors



54W x 38W x 85H cm

\$245

AX286

Freight charges will apply

Single Linen Skip with Foot Operated Lid

ParagonCare

- Frame designed to support bag
- Easy removal of bag
- Sturdy frame construction
- Heavy duty 75mm castors
- Foot operated polymer lid



54W x 38W x 85H cm

\$309

AX291

Freight charges will apply

Polymer All Purpose Carts

ParagonCare

The polymer cart range is ergonomically designed and has numerous applications. Constructed of reinforced nylon fibres, the carts have a 250kg load capacity and will withstand the harshest working conditions. They are also an excellent cart for conveying lightweight items with little effort.



Freight charges will apply

AXH0801	2 tier small	85L x 48W x 95H cm	\$229.00
AXH0802	2 tier large	10L x 52W x 95H cm	\$379.00
AXH0803	3 tier small	85L x 48W x 102H cm	\$295.00
AXH0804	3 tier large	110L x 52W x 102H cm	\$495.00

M40 Computer Cart

ParagonCare

The M40 Computer Cart is an entry level cart offering a comprehensive set of nursing features – 16 inch height adjustment, internal storage for small form CPUs and a keyboard light – making them ideal for point of care applications.

\$1,325

AXM40



Download the **FREE** Vital Medical Supplies **APP**

VISIT PAGE 10 FOR MORE INFORMATION

NEW

FREE download available from:





Step Up Stools

Powder Coated

- Powder coated frame
- Non-slip step surface

Freight charges will apply



ParagonCare

Single		
AX375G	42L x 30W x 20H cm	\$95.00
Single with Handle		
AXH0378	With Handle 42L x 30W x 20H cm	\$205.00
Double		
AX393G	54L x 43W cm Step Height 20cm and 40cm	\$165.00
Double with Hand Rail		
AXH0402L (Left)	With Hand Rail 50L x 45W cm 120H cm Step Height 20cm and 40cm	\$295.00
AXH0402R (Right)		

Hand rail option left or right

Step Up Stools

Stainless Steel

- Stainless steel frame
- Aluminium tread plate step surfaces

Single			
AX387	Model 1060 Chrome	40L x 30W x 20H cm	\$175.00
AX390	Model 1061 Powdercoated	70L x 30W x 20H cm	\$205.00
Double			
AX405	Model 1062 Chrome	57L x 40W cm Step Height 19cm and 37cm	\$219.00



Freight charges will apply

Brooksbuilt Timber Step-Up Stools

- Timber colours: Ash, Oak, Beech or Teak
- Stool tops available in Black or Grey

Single		
18000950	45L x 30W x 23H cm	\$119.00
Double		
18000970	50L x 54W x 30H cm	\$189.00



Sandel® Ergo-Step™ Stool

- Lightweight (3.2kg) made of moulded plastic with a cushioned top for added comfort
- Can be stacked for added height (max: 2 high)
- Can be interlocked (side by side) to create work platform
- Anti-Slip surface



Sandel® Ergo-Step™ Stool			Each
AN1170T	Top* or	45.7cm(L) x 35.6cm (W) x 12.7cm (H)	\$99.00
AN1170B	Bottom*		\$138.00

*Top step stools (AN1170-T) should ONLY be stacked on bottom step stools (AN1170-B). Tops are not designed to be stacked on one another.

Sandel® Ergoplus Anti-fatigue Mat

SANDEL®

Helps reduce the stress and strain on muscles and joints caused from standing during long procedures.

- Premium padding for additional comfort
- Fluid resistant
- Stain resistant
- Anti-Slip surface



AN2317-P	Small 33cm x 43cm	\$193.00
AN2332-P	Medium 50.8cm x 81.3cm	\$495.00
AN2348-P	Large 50.8cm x 122cm	\$739.00

Ansell



Theatre Stool

ParagonCare

- Height Range 530-730mm
- Black vinyl upholstered seat
- Five star base
- 50mm twin castors
- Hand operated height adjustment



\$175

AX262

Freight charges will apply

Theatre Stool with backrest

ParagonCare

- Black vinyl upholstered seat and backrest
- Backrest with lumbar adjustment
- Five star base
- 50mm twin castors
- Hand Operated height adjustment



\$409

AX265

Freight charges will apply

Surgeon Stool Foot Operated

ParagonCare

- Heavy duty weighted base for stability
- Enclosed base cover for infection control
- Easy to use and hands free. Height adjustment via foot control
- Height range 530 – 710mm
- Black vinyl upholstered seat
- 75mm heavy duty castors



\$1,095

AX267

Freight charges will apply

Saddle Stool

ParagonCare

- Hand operated height adjustment
- Height adjustable
- 530mm-730mm
- Five castor polished base
- 50mm castors



\$295

AXH0271

Freight charges will apply

Round Black Stool

- Adjustable height: 49cm-63cm
- Cushion diameter: 35cm
- Foam thickness: 6cm



\$99

PMRSBL

Freight charges will apply



Gas Lift Saddle Stool

- Adjustable height: 49cm-63cm
- Cushion diameter: 35cm
- Foam thickness: 6cm



\$149

PMSSBL

Freight charges will apply





RTS Swivelling Telescopic Arm

ParagonCare

You won't trip over anything hanging on the wall. This is the best reason for the use of telescopic arms, the screening system which is only there when it is needed. Following use of the screen, simply push together and turn toward the wall. Your room is then restored to its usual appearance.

The RTS Telescopic Arm comes in 3 different lengths and 3 great colours. Please indicate your preferred colour when ordering.

RTS 1.2 Telescopic arm		
AXHRTS12	55-120cm	\$375.00
Curtains		
AXHTCS123W	White 140W x 175L cm	\$139.00
AXHTCS123B	Blue 140W x 175L cm	

RTS 2.1 Telescopic arm		
AXHRTS21	85-210 cm	\$435.00
Curtains		
AXHTCS213W	White 230W x 175L cm	\$155.00
AXHTCS213B	Blue 230W x 175L cm	

RTS 1.5 Telescopic arm		
AXHRTS15	65-150 cm	\$425.00
Curtains		
AXHTCS153W	White 170W x 175L cm	\$149.00
AXHTCS153B	Blue 170W x 175L cm	



Freight charges will apply

Ropimex is a versatile system. A number of combinations can be made. Above we have listed the most common combinations. Curtains are also available in 120cm and 140cm length. More colours are available but prices could differ. Delivery will be approx 4 weeks for items listed above. Different colours etc. will take 6-8 weeks. If you wish to have a different configuration, please do not hesitate to ask.



RLP Lightweight folding screen

ParagonCare

Ropimex® lightweight folding screens convince by their ease of handling. With the 360° hinges, every panel may be moved into the desired position

AXHRLP2	2 Section	101W x 165H	\$570.00
AXHRLP3	3 Section	156W x 165H	\$750.00
AXHRLP4	4 Section	204W x 165H	\$985.00
AXHRLP5	5 Section	259W x 165H	\$1,150.00

Freight charges will apply

Screen Colour

Please indicate your preferred colour when ordering.

Screen Material

Staph Chek® fabric. The material is antibacterial, flame resistant, self-extinguishing, antistatic, dirt-resistant and long lasting. It does not need to be washed. Marks are simply wiped off with sponge and cleaning agents.



Privacy Screens

ParagonCare

Powder coated frame, 50mm castors and includes curtains.

AX420	2 Fold	Panel Size 90W x 1401H cm	Unfolded 180W x 170H cm	\$425.00
AX411	3 Fold	Panel Size 90W x 1401H cm	Unfolded 270W x 170H cm	\$585.00

Freight charges will apply

Cubicle Curtains & Tracks

- Supply only
- Easy to install
- Fittings included
- Fire retardant fabric certified to meet Australian standards
- 100% polyester

Curtain Tracks			
LNH600SJ	J Satin Silver	2400x1500mm	\$320.00
LNH600SL	L Satin Silver	2400x1500mm	\$320.00
LNH600SS	Straight Satin Silver	3800mm	\$275.00
Curtains			
LNGMS8S	Curtain For 4M Track Snow (Off White)		\$395.00
LNGMBB6	Curtain For 4M Track Bluestone		\$395.00
LNGML29	Curtain For 4M Track Lichen		\$395.00



Bluestone




Snow



Lichen

L&N Interiors



PLAN VIEW



Opening an account with us is easy.

Call **1300 557 651** and speak to one of our Customer Service team members. Alternatively you can open an account online at www.vitalmed.com.au/new-accounts



Mayo Instrument Table

- 304 Stainless Steel Construction
- Telescopic height from 78 - 120cm
- 75mm smooth castors with locking front wheels
- Removable tray 60 x 40



NEW

\$395

PMMIT

Freight charges will apply

ParagonCare

Mayo Instrument Table

- 4 wheel design
- Stainless steel construction
- Height range 82-112cm
- Tray size 56 x 41 cm
- 50mm castors



**EXCELLENT
QUALITY!**

\$589

AX256

Freight charges will apply

Mayo Instrument Table

ParagonCare

- Mayo Table - Three Leg
- Stainless steel construction
- Height range 82-112cm
- Tray size 56 x 41cm
- Three leg base
- 75mm castors



\$615

AX254

**EXCELLENT
QUALITY!**

Freight charges will apply

Overbed Table

ParagonCare

- Top Size: 790 x 360mm
- Height Range: 73 -1050mm
- Height adjustment
- Spring assisted locking knob
- Castors: 50mm Twin
- Malamine top with PVC "T" edging
- Powder coated frame - Ivory



\$225

AX9700

Freight charges will apply



IV Stand ParagonCare Stainless Steel

- Telescoping
- Height range 135–205cm
- Stainless steel construction
- 50mm castors
- Stable five-star weighted base

\$249

AX217

Also available
with black polymer base

\$249

AXH0216

Freight charges will apply



Smart IV ParagonCare Stand

- One hand adjustment stand which enables quick and easy raise and lower of the pole
- Height range 130–210 cm
- Stable five star weighted base
- Ideal for use with IV pumps and syringe pumps
- 80mm castors
- Telescoping
- 4 hook IV holder
- Aqua seal

\$225

AX2380

Freight charges will apply



IV Stand ParagonCare Polymer Base

- Lightweight – 3kg
- Height range 135–215cm
- Height adjustment with locking screw
- Base diameter 635mm
- Bottle holder 4 hook
- 50mm Castors (2 lockable)
- Colour blue
- Easy to assemble

\$135

AXH2386

Freight charges will apply



IV Pole ParagonCare Foldable

Easy and quick to
use in emergency
and homecare

\$249

AXH84132

Freight charges will apply



Mobile ParagonCare Base Adjustable Height

- Single hand operation
- Easy to assemble
- Cable clip
- White metal powder coated
- Base diameter – 680mm
- 50mm Castors (2 lockable)
- Height adjustable – 1070mm-1470mm

\$215


AXH2384

Also available:
Non-adjustable
mobile base
(1070mm)

\$165

AXH2382

Freight charges will apply



Mobile Base

- Stable 5-star solid base
- Black 50mm twin wheel castors
- Chrome plated stem with key wheel
- Height – 850mm stem
- Diameter – 12.5mm

\$149

VIDCD

Freight charges will apply





The Play Panel

Features:

- Play panels are designed for durable use so they are safe, durable and easy to maintain
- Conform to Australian Standards (AS/NZ ISO:8124)
- Wall mounted so they don't take up valuable floor space
- Ensure there is no ongoing cost of replacing broken or stolen toys
- Can be wiped regularly with antibacterial wipes to help prevent cross-infection
- Leave no messy toys for staff to continually tidy
- Do not create a tripping or choking hazard
- Are the preferred play option for many Australian hospitals including Royal Children's Melbourne, Mater Brisbane, Sydney Children's, Royal North Shore, Children's Hospital at Westmead and over 500 medical practices.

Wall Mounted		
CHWUR	Red	\$425.00
CHWUB	Blue	
CFS100Y	Yellow	
CHWUP	Purple	

Please specify colours when you order



The Play Centre


- Free-standing
- Under 1m² floor space required
- All attachments are included
- Available colours: red, blue, purple or yellow

CHFS3	3 Play Panels	\$1,350.00
-------	---------------	-------------------

Please specify colours of all 3 panels at time of ordering



A safe, tidy, hygienic
play solution
for waiting rooms



Panels come in your
choice of four vibrant
colours

Freight charges will apply