

Whiter, Brighter AND Greener



Welch Allyn Green Series™ LED Medical Exam Lights

- Lower cost of ownership due to elimination of bulb replacement and lower energy consumption
- Helps you see true colour tissue rendition, and provides a uniform light spot
- Reduced risk of cross-contamination for your patients with touchless on/off and optional disposable sheaths
- 5-year limited warranty†

GS 300

44456 GS 300 General Exam Light with Mobile Stand

44416 GS 300 General Exam Light with Table/Wall Mount



GS 600

44606 GS 600 Minor Procedure Light with Mobile Stand

44616 GS 600 Minor Procedure Light with Table/Wall Mount

GS 900

44906 GS 900 Procedure Light with Mobile Stand

44900-C GS 900 Procedure Light with Ceiling Mount

44900-W GS 900 Procedure Light with Wall Mount



EXAM LIGHT

48816 GS Exam Light IV with Table/Wall Mount






Green Series™ 300 General Exam Light

Utilising a single white LED, the GS 300 Exam Light provides a whiter, brighter light than halogen and consumes less energy.

Features

- Perfect combination of intense light output and cool operation in a compact design
- Affordable price and rugged performance
- Touchless On/Off
- Multiple mounting options for floor space maximisation
- Lower cost of ownership and energy consumption due to elimination of bulb replacement
- Easy assembly
- 5 year limited warranty to include the LED light source and power supply only



WL44416	GS 300 Exam Light with wall/table mount	\$645.00
WL44456	GS 300 Exam Light with mobile stand	\$750.00


WelchAllyn

Green Series™ 600 Minor Procedure Light

Featuring three LEDs and a wider head, the GS 600 is designed specifically for demanding viewing conditions.

Features

- Three white LEDs for greater intensity
- Touchless power controls and adjustable light output
- Disposable sheath available to help reduce the risk of cross-contamination
- Rugged, ergonomic design allows for improved performance and long-lasting durability
- Low cost of ownership due to elimination of bulb replacement and lower energy consumption
- Easy to assemble and clean
- 5 year limited warranty to include the LED light source and power supply only



WL44616	GS 600 Procedure Light with wall/table mount	\$1,149.00
WL44606	GS 600 Procedure Light with mobile stand	\$1,295.00

WelchAllyn

Green Series™ 900 Procedure Light

This procedure light, featuring six LEDs, combines outstanding durability, extraordinary maneuverability, unparalleled spot quality and multiple mounting options.

Features

- Six white LEDs for greater intensity
- Touchless On/Off with dimming option
- Multiple mounting options for floor space maximisation
- Disposable sheath available to help reduce the risk of cross-contamination
- Rugged, ergonomic design allows for improved performance and long-lasting durability
- Lower cost of ownership due to elimination of bulb replacement and lower energy consumption
- 5 year limited warranty to include the LED light source and power supply only



WL44900-C	GS 900 Procedure Light with ceiling mount	\$2,995.00
WL44900-W	GS 900 Procedure Light with wall mount	\$2,995.00
WL44906	GS 900 Procedure Light with mobile stand	\$2,995.00

WelchAllyn

Green Series™ Exam Light IV

Perfect for the OB/GYN office or for specialties that require very focused light in specific areas. The GS Exam Light IV provides you with 3x the typical light output and features a highly flexible light pipe that is twice as long as conventional pipes. With a compact, unobtrusive head, it ensures pinpoint accuracy and provides true tissue colour rendition during an exam.

Features

- One LED - 3x the typical light output; pinpoint accuracy; true tissue colour rendition
- Highly flexible light pipe for easy maneuverability
- Lightweight pipe is twice as long as conventional pipes
- Touchless On/Off
- Multiple mounting options for floor space maximisation
- Disposable sheath available to help reduce the risk of cross-contamination
- Lower cost of ownership due to elimination of bulb replacement and lower energy consumption
- Easy to assemble
- 5 year limited warranty to include the LED light source and power supply only



\$1,175

WL48816

WelchAllyn

HEINE QUALITY
MADE IN GERMANY

HEINE EL10 Examination Light

A real highlight in all respects with
LED in Heine Quality Illumination

Experience The Difference

LED NOW

LED NOW IN HEINE QUALITY.

- Continuous adjustment
- One-hand operation
- Compact design





HEINE EL10 LED® Examination Light

Only a few lights are as bright, none are as homogeneous as this. Illumination without compromises for every application. Stationary and yet flexible. Position, brightness, and spot size can be adjusted with one hand.

LED –economical with high colour fidelity. The new examination light combines many characteristics that facilitate everyday working life.

- 5 year guarantee
- Bright and homogenous light
- Adjustable spot size
- Brightness control
- Easy to position
- One handed operation
- Compact design
- Hermetically sealed housing
- Various mounting options available with Stand, Wall and Clamp mounts



HEJ0827001	EL10 LED Examination light with Wall Mount	\$2,215
HEJ0827002	EL10 LED Examination light with Clamp Mount for Table Top Mounting	\$2,365
HEJ0827003	EL10 LED Examination light with Wheeled Stand & Metal Base	\$2,575



LED HQ
LED NOW IN HEINE QUALITY.



Ri-magic® Examination LED Light

The flexible examination light ri-magic® LED illuminates the examination and operation area in the best possible way.

- The most modern LED technology, similar to daylight and high-performance; ideal for examinations and smaller procedures
- Lifetime of the lamp up to 50,000 hours, replacements no longer necessary
- Significantly less heating and energy loss
- The focusing device at the front of the light head allows exact adjustment of the size of the illuminated examination area
- Highly flexible, long special arm can be bent in all directions, the light remains fixed in the chosen position
- Simple to use handle with rheostat for infinitely variable regulation of light intensity
- Wide-range power cord for 230V/120V



RIE1010235	Wall mount	\$875.00
RIE1010236	Mobile mount	\$995.00
RIE11343	Mobile stand and adaptor	\$325.00



HEINE EL3 LED Examination Light

Simply bright LED Illumination

The EL 3 LED provides bright (> 42,000 lux / > 350 lumen), shadow-free illumination precisely at the point of interest. at a working distance of 30 cm.

- Shadow-Free Illumination. The compact illumination head (Ø approx. 60 mm) allows convenient, precise positioning, particularly in difficult application situations. Optimal positioning eliminates unwanted shadows in the illuminated field.
- Maintenance-free and easy to clean.
- The light weight head and stable goose neck allow easy and intuitive adjustment and keep the light exactly in the required position.
- LED in HQ. The new standard in LED illumination defining optimal light intensity, homogeneity and colour rendering for each individual examination situation – Colour temperature: 5,000 K, Colour Rendering Index (CRI) > 90 on a scale of 100. This ensures true tissue colour for the most accurate diagnosis.
- 134cm Goose Neck allowing for head to toe examination
- Hermetically Sealed Housing
- No dust traps
- Easy to clean



HEJ0827011	EL 3 LED Examination Light with wall mount	\$1,029.00
HEJ0827013	EL 3 LED Examination Light with clamp for table-top mounting	\$1,175.00
HEJ0827014	EL 3 LED Examination Light with wheeled stand, metal-base	\$1,385.00

LED HQ
LED NOW IN HEINE QUALITY.



Derungs Visiano LED Examination Light

The Visiano LED exam light made by Derungs offers the best in LED medical lighting technology, combining low heat output with adjustable colour and brightness. This provides medical staff complete control when carrying out minor procedures and general examinations.

Outstanding features of these LED Medical Lights include:


- Diamond Optic Technology producing 60,000 lux at 0.5m
- 16 premium LEDs with life-span > 50,000 hours
- 4 dimming levels
- Changeable colour temps. (4500K cool/3500K warm) for bringing out the most subtle colour tones and details
- Total freedom of movement with 360 degrees rotatable light head and integrated handles
- Large light field diameter of Ø 21cm. 20% more than conventional luminaires
- Shadow removal with wider light head coverage
- Sanitized® material with silver ions gives protection against spread of bacteria
- Fully enclosed light head, rated according to IP20. Cleanable with disinfectant spray
- Highly adjustable and robust spring-balanced arm, tested to 20,000 movement cycles
- The light colour temperature of these lamps can be changed by a simple touch of the control panel on the lamp head. The user can select the optimum colour for the examination depending on the skin area and type of tissue or wound. The two colour temperatures are selectable on the touch panel



Ceiling Mount



Wall Mount



Mobile Mount

MGL514	With Ceiling Mount	\$2,995.00
MGL515M	With Mobile Mount	\$1,995.00
MGL515W	With Wall Mount	\$1,975.00



Colour Temp 3500K



Colour Temp 4500K





Derungs Visiano 10-1 LED Examination Light

For everyday use, the Visiano 10-1 will win you over with its excellent illumination intensity and large light field. Having a sturdy flexible arm ensures positional accuracy and good durability. The sealed construction protects against contamination.

- High illumination strength 50'000 lux / 0.5 m.
- Illuminated field 180 mm diameter
- Colour temperature 4400 K.
- 280° rotation of the luminaire head.
- Slimline design with ergonomic handle.
- Long arm reach.
- Maintenance-free.
- 2 models- wall mount or mobile complete with stand



MGL501W	With Wall Bracket	\$1,250.00
MGL501S	With Mobile Stand	\$1,325.00



Visiano 10 Exam Light

- Long gooseneck arm reach
- Slimline design
- Wall / table mount

Opticlux - UV Magnifier

- Ceiling / wall mount option
- White / UV light
- Large 15cm lens



OPTICLUX COMING SOON!



Halux LED Examination Light

The Halux LED is an ergonomically designed, easy to use light which gives 40,000 lux at 0.5 metres. Selectable light colour is ideal for different procedures with dimmable light output.

Features:

- 40,000 Lux at 0.5 metres
- 180mm Ø beam width
- 7 Powerful LEDs
- 3 colour levels – 3300K, 3800K and 4400K
- 5 dimming levels from 100% down to 20%
- Easy to clean smooth edges
- 950mm reach, Spring-balanced arm system for precise light positioning
- Dimming/colour change switching on lamp head
- Option with wall mount or mobile adapter bracket to suit most mobile stands



Mobile Mount



Wall Mount

MGL509M	With Mobile Mount*	\$1,495.00
MGL509W	With Wall Mount	\$1,475.00

*Mobile Base not included



Halogen Examination Light

Cost-effective solution for surgery examination lighting

- 35 watt globe, 12V (approx 4000hrs)
- 70 cm long gooseneck
- Easy grip handle with on/off switch
- Additional glass filter reduces heat

\$295*

MGL308

Brackets to be purchased separately



LED Examination Light

- 700mm Gooseneck
- Sure grip handle with on/off switch
- Cool white light, 5500°K colour,
- Concise beam angle of 30°
- 3 x 1 watt high power LED giving 16,000 lux at 30cm
- Low energy consumption, only 5 watt
- Low heat output
- Wall or mobile mounted
- 40,000 hour bulb life
- 2 year warranty

\$495

MGL308LED





LUXO Quality Examination Lighting



GET VITAL APP
Scan Here

LHH G2 LED

Light source: LED
CRI: 95 (R9:98)
CCT: 4000K
Lifetime: 50,000 hours
Lux level @ 0.5m: 120,000 lx
Lux level @ 1.0m: 35,000 lx
LFD@ 1.0m: 300mm d10
Control: On/off switch
Seamless dimming
Night light:
Starts at 1% light

LHH G2-TW

Light source: LED
CRI: 95 (R998 4000K)
CCT: 3500K,4000K,4500K
Lifetime: 2:50,000 hours
Lux level @ 0.5m: 65,000 lx
Lux level @ 1.0m: 25,000 lx
LFD@ 1.0m: 300mm d10
Control: On/off switch, dimming
Control: Tuneable white cycle

The best in its class.

LHH G2 LED Examination Light

Its highly flexible self-balancing arm in conjunction with excellent light output gives the LHH G2 diagnostic luminaires vast popularity among medical professionals. The frictionless 3d movement and ease of use make the LHH G2 the most flexible examination luminaire on the market... and the best in its class.

Also available with **Tuneable white** illumination, offering a high output LED module that can be adjusted to one of three colour temperatures.

Features:

- SteriTouch® antimicrobial white powder-coated aluminium and steel construction.
- 3D positioning with 3-axis movement and perfect balance.
- Protective covers on arms to prevent crush-points.
- Smooth membrane switch for ease of use and cleaning.
- Fully dimmable with night light function.
- High output LEDs with refractor lenses for increased light output.

LHG026530	Wall/Ceiling mount	\$1,275.00
LHG026531	Table/Trolley mount	\$1,275.00

LHH G2 Tunable White Examination Light

LHG026542	Wall/Ceiling mount	\$1,395.00
LHG026543	Table/Trolley mount	\$1,395.00





Green Series Mobile Stand Examination Light IV

\$349

WL48950

Examination Light Tall

\$349

WL48960



WelchAlllyn

Freight charges will apply

Mobile Base for LHH Examination Light

Freight charges will apply

\$345

LLTS



LUXO

Mobile Base

- Stable 5-star solid base
- Black 50mm twin wheel castors
- Chrome plated stem with key wheel
- Height 850mm stem
- Diameter 12.5mm

Freight charges will apply

\$149

VIDCD



Extended Wall Mount for GS300/600 Examination Light



\$255

WL44215

WelchAlllyn

Ceiling Mount Bracket for LHH 60cm



\$355

LU2529

LUXO

Wall Bracket and Base



\$75

LUX899

LUXO

Maggylamp Wall Bracket



\$65

MGL307

MAGGYLAMP
medical

Maggylamp Mobile Base Bracket



\$52

MGL309

MAGGYLAMP
medical

ri-magic Mobile Stand + Adaptor

Riester

\$385

RIE11343



Disposable Sheaths for Green Series™ 600/900

\$205

WL52630
Box 250



WelchAlllyn

Disposable Sheaths for Green Series™ Exam Light IV

\$39

WL52640
Box 25



WelchAlllyn

35 Watt Halogen Bulb for Maggylamp

\$24

MGL306



MAGGYLAMP
medical



Tevisio LED Magnifier Light

A visual field designed for optimum performance, the TEVISIO magnifier luminaire combines high efficiency with ergonomic design.

Features:

- Unique lower arm assembly allows head to swing below bench/shelf height
- 48 premium SMD LEDs
- Continuous dimming to 10%
- 2 x magnification
- 6,000 lux output using only 14 watts
- Visualiser function for 3D shadow effects
- Large arm radius with innovative flexible head joint
- Balanced arm, easily adjustable
- Maintenance-free design with LED life up to 50,000+ hours
- Wide-angle lens 153mm (diameter) for distortion-free vision
- 784mm reach
- 2 year warranty



\$995

MGL160

MGL162	Add-on Lens	\$295.00
--------	-------------	-----------------

Derungs
MEDICAL LIGHTING

WAVE LED Magnifier Light

Shadow-free magnification with light from both sides, as well as three-dimensional magnification with light coming from either the left or right. The 3D feature is especially welcome when working with circuit boards and similar delicate objects.

Features:

- Steel arm, shade in aluminium
- Standard with 3.5 dioptre (1.9x magnification) crown glass lens, 175x108mm
- Secondary 4D, 6D or 10D lens can be attached for 2.9x, 3.4x or 4.3x magnification respectively
- 784mm reach
- 2 year warranty
- Supplied with 2 x integrated 6W LED modules, individually switchable
- Step dimming (0-50-100%)
- Automatic shut off after 4 or 9 hours



\$575

WAL025960

LUXO

MAGGY-VUE Magnifying Lamp

- 4 Diopter circular lens, Ø 120mm
- 2 x magnification
- Shadow-free diffused lighting
- 22 Watt fluorescent lamp (included)
- Spring balanced joints
- Optional clip-on lens for increased magnification (3x)



QUALITY
MADE IN
GERMANY

MAGGYLAMP
medical

MGL150WM	Wall Model	\$425.00
MGL150MM	Mobile Model	\$425.00
MGL151	Add-on Lens	\$155.00

LAMP EXAMINER LFM30-101

- Magnification and cool lamp
- Round 4 diopter lens
- Internal spring-balance arm with parallel movement
- Magnification ratio: x 100% actual size



\$225

LULKMAGW

LUXO

Hand-Held LED Maggylamp With Wall Bracket

Designed for the medical profession with the added features of a wall bracket and extendable curly cord. The wall bracket gives the user more bench/table space. The cord stretches up to 2.5 metres from the mounted wall bracket.

Ideal for viewing:

- Skin cancers
- Lesions
- Fungi
- Warts
- Rashes/hives
- Stitches
- Ingrown hairs



2x
Magnification

RB216	MAGGYLAMP hand-held	\$595.00
RB218	MAGGYLAMP with wall bracket	\$650.00

MAGGYLAMP
medical



HEINE mini-c® Clip Lamp



Compact, long-life pocket diagnostic light

Sturdy, long-life pocket diagnostic light in a compact, modern design with chrome head. Bright, concentrated light thanks to novel XHL Xenon Halogen lens-lamp. Mini-c clip lamp only fits a mini-c handle. Not compatible with mini 3000 instruments.

- Compact, oval, attractive light
- Fits ergonomically in the hand
- High-quality handle: Shockproof, sturdy, non-slip
- Switches off automatically when put back in pocket
- Replaceable batteries. Size AAA

\$94

HED0173109



Ear Light for mini-c Clip Lamp

6 pieces in blister pack

HED0073105

Pkt/6

\$114.00




Diagnostic Penlight

- Modern satin finish aluminium body – delux model
- Pocket clip doubles as on/off switch
- Useful 5cm rule and dilation guide
- 2 x AAA batteries included

\$24.⁵⁰

SAS233



Promar Torches

ParagonCare

- Sealed unit with 2xAAA size batteries included - non rechargeable, non replaceable batteries
- Concentric beam bulb with pocket clip/switch
- Colour - White



	Each	+6
LRDDT	\$3.95	\$3.50

WOODS LIGHT Ultra Violet

For the detection of various skin disorders

- Convenient hand size
- Rugged plastic construction
- Plastic lens over tube to protect patients' eyes
- Includes torch



\$25

LRUVL

Tube for Woods Light

LRUVT

\$9.50

ParagonCare

Handheld LED Illuminated Magnifier

\$29

MGL485

- 85 x 50mm lens
- Crystal clear acrylic lens with
- 2.0 X mag and 4.0 X mag bubble
- Complete with batteries and pouch
- Ergonomic shaped handle fits comfortably in the hand



MAG-LITE®

High Quality Examination Torches

Mini Maglite

AAA

\$29.⁹⁵

SH85200

The Mini Maglite® 2-Cell AAA flashlight is a smaller and lighter version for ease of use and storage.



- High-intensity light beam
- 1/2 turn, twist focus, spot-to-flood
- Spare lamp safely secured in the tailcap
- High-strength aluminium alloy case
- O-ring sealed for water resistance

AA

\$34.⁰⁰

SH85700

MagLite 3 Cell

A very popular D cell Mag-Lite. Featuring a long handle, extending its length to 31cm.



\$49

SH86230

MAGGYLAMP UV Woods 4x4 WATT

- Fitted with 4 x 4 Watt UV light
- UV Wavelength Lens: black light blue, 365nm biconvex glass, rectangular 100 x 50mm
- Magnification: 2.25 x 5 dioptre

\$549

RB215

MAGGYLAMP
medical

NEW





LoupeLight 2 with mPack mini for i-View loupe mount (HR / HRP)

Compact LED HQ Headlight – now even brighter

The combination of light and loupe delivers the ultimate in performance as high quality illumination with LED HQ maximises the image quality of the loupe. For a perfect view and a precise diagnosis, the LoupeLight 2 is now 25% brighter than its predecessor model, with a homogeneous edge-to-edge illumination of typ. 55,000 lux.

For intensive use, the LoupeLight 2 now features a replaceable connecting cord which connects the device to the mPack mini – the new and compact power source for the LoupeLight 2.

- Compact and Lightweight
- Bright, homogeneous and coaxial illumination (typ. 55,000 lux)



HEINE mPack Mini

The new, high performance power source HEINE mPack mini for LoupeLight 2 and MicroLight 2 offers freedom from cables that tie the users to desk, wall or static power sources. With a weight of only 95 g, this mobile power source is now approx. 70 % smaller and lighter than its predecessor.



HEC0832265	LoupeLight 2 F HR/HRP + mPack Mini	\$2,275.00
HEC0032241	Yellow filter F.LED for LoupeLight 2	\$132.00
HEX0099007	mPack Mini Belt Clip	\$55.00
HEC0032240	Connecting Cord SC1 (1,5m/Ø 2,4mm)	\$197.00
HEC0032243	Connecting Cord SC2 1.5m (1,5m/Ø 3,2mm)	\$197.00

LED HQ
LED NOW IN HEINE QUALITY.



HEINE LoupeLight 2 Sets with HR Binocular Loupes

For a perfect view and a precise diagnosis, the LoupeLight 2 is now 25% brighter than its predecessor model, with a homogeneous edge-to-edge illumination of typ. 55,000 lux.

For intensive use, the LoupeLight 2 now features a replaceable connecting cord which connects the device to the mPack mini – the new and compact power source for the LoupeLight 2.

- Compact and lightweight
- Bright, homogeneous and coaxial illumination (typ. 55,000 lux)
- Optional polarisation filter system

Includes:

- LoupeLight2
- mPack mini
- mPack mini belt clip
- Yellow Filter
- E4 USB Plug-In Transformer
- 2 x Protective Lenses
- 2 x Sterilisable Swivel Levers
- Retaining Cord
- Cleaning Fluid
- Cleaning Cloth
- Carry Case



HEC0832450	HR 2.5x/340mm, LL2, I-View A.S-Frame	\$4,655.00
HEC0832451	HR 2.5x/420mm, LL2, I-View A.S-Frame	
HEC0832452	HR 2.5x/520mm, LL2, I-View A.S-Frame	

LED HQ
LED NOW IN HEINE QUALITY.





HEINE S-Frame + HR® 2.5 x Binocular Loupe with i-View

HEINE S-Frame®

- Innovative spectacle frame for HR Binocular Loupe
- Appealing and functional design
- Can be customised to fit any shape of head and face

HEINE HR® Binocular Loupe 2.5X

- The brightest image – Heine’s exclusive coating technology reduces reflection to 0.3%, providing the brightest image of any comparable loupe
- Crispest image with colour correction
- Super-wide, super-deep field of view
- Dust and water resistant



HR Binocular Loupe Set with i-View Optics, S-Frame with i-View loupe mount for HR 2.5 x loupes, detachable protective lenses, 2 sterilisable swivel levers, retaining cord, cleaning fluid, carrying case.



HEC0032355	340mm working distance	\$2,190.00
HEC0032356	420mm working distance	
HEC0032567	520mm working distance	

Optional Polarisation Filter



HEC0032535	Polarisation filter P2 for LED Loupelight and HR Loupes	\$180.00
HEC0032536	Spare polarisation filter P2 for HR Loupes (2 pairs)	\$90.00



This Catalogue has been printed on paper sourced from Sustainable Forests under the PEFC accreditation.



Here's a bright idea!
Why not visit
www.vitalmedicalonline.com.au
to place your lighting
order today.



HEINE ML4 LED Headlight

Ultra Bright Coaxial Illumination – Now with LED ^{HQ}

- **LED in HQ** - the new standard in LED illumination defining optimal light intensity, homogeneity and colour rendering for each individual examination situation.
- **Adjustable illumination spot size** - 30 mm to 80 mm spot size range (at 420 mm working distance) to fit any examination situation.
- Bright and homogenous illumination. Absolutely bright light spot that is uniform from edge-to-edge for the perfect illumination in all examination situations.
- **Stepless light intensity control** - Optimal brightness setting prevents reflexes.
- **Flexibility**. Flexible power source options: 100% mobility with a choice of cable-free mPack UNPLUGGED head-worn battery or mPack belt-worn battery pack.
- **Coaxial illumination** - Coaxial design ensures a completely shadow-free image and allows for excellent illumination of difficult to see areas.
- **Comfortable and secure fit** - The ML4 Professional headlight has multiple adjustment points and soft padding for comfort and stability – even during long examinations.
- **Optional Filters** - Polarisation filter P 2 improves contrast and yellow filter reduces blue light.

HEINE QUALITY
MADE IN GERMANY



\$1,915

HEJ0831410

HEJ0831413	ML4 Headlight with mPack	\$2,795.00
HEJ0831415	ML4 Headlight with mPack unplugged	\$2,950.00
HEX0799666	mPack Unplugged + Plug Transformer	\$1,095.00

HEINE mPack UNPLUGGED

Headband-mounted rechargeable battery

- **Wireless** - 100% mobility for absolute freedom of movement.
- **Rechargeable battery integrated on the headband** - No inconvenient cable connections.
- **Lightweight with balanced weight distribution** - High level of comfort.
- **Charge Status Display** - Indicates remaining operating time.
- **Intelligent** - Can be used as mains power supply even in "charge" mode.
- **Fast Charge** - Charges to full capacity in only 2 hours. Operating time 3.5 hours of continuous ON time at full power.
- **Top Safety Standards** - Integrated charge and performance checks during operation. Lithium polymer technology. UL compliant.



HEINE mPack

- **Portable power pack** - Does not tie the user to desk, wall or static power supplies.
- **Li-ion Technology** - Triple capacity compared with current NiCd systems. No memory effect, autonomy for 8.5 hours of continuous ON time at full power.
- **Fast Charge** - Charges to full capacity in 2 hours.
- **Charge Status Display** - Always lets you know how much operating charge you have left.
- **Stepless light intensity control** - Optimum brightness setting prevents reflexes.
- **Flexible Charging Options** - EN 50 Desk and wall charger, or plug-in transformer.
- **Intelligent** - Can be used as mains power supply even while in "charge" mode.
- **100 – 240 V Compatible Transformer** - incl. 4 plugs for world wide compatibility.



Mag-eyes The hands-free, head mounted magnifier

Mag-Eyes with 2 Lenses
(1.6/2.0 X Magnification)

\$79

MGL400

Mag-Eyes with 2 Lenses
(2.25/2.75 X Magnification)

\$95

MGL410



Working distances (Lens to Subject)	
1.6 X (30-35cm)	2.0 X (25-30cm)
2.25 X (18-25cm)	2.75 X (10-18cm)

Guthrie Loupe

ParagonCare

- Adjustable head band padded for user comfort
- Plastic, Precision Lens
- Retractable visor
- 2.7x Magnification



\$115

LRM2001