



Deluxe Oxygen Therapy Bag

This bag can carry a D or C size oxygen cylinder with the regulator installed and ready for action. Carry handles and fold away backpack straps

- High Vis Reflective straps and tie down loops
- Incredibly durable 1200 Denier coated TPE, water and stain resistant
- Removable waterproof skid sheet
- Removable top bag for general First Aid equipment



\$295

AD12016

Silicone Resuscitator

Made from transparent durable silicone

Kit includes:

- 1500ml resuscitator
- 2500ml oxygen reservoir
- 2 silicone masks – small adult and large adult
- 3 Guedel Airways
- See-through carry case

\$175

VIFN045
GST Free



AUTOCLAVABLE TO 134°

Multigate Ark Resuscitation

The Airway Resuscitation Kit helps you keep focus.

There is no need to find all the necessary devices for airway resuscitation. They are already in a single case.

Contents include:

- Two Yeescope Laryngoscopes
- Endotracheal Tubes (cuffed)
- Laryngeal Mask Airways
- Intubating stylet

\$245

AN010004



Disposable Resuscitators

Made From Transparent Durable Vinyl

Each resuscitator includes: bag, mask, patient & intake valves, reservoir bag and 2 metres oxygen tubing.

Patient valve incorporates vomit deflector/peep valve adaptor.



FN040	Infant	\$29.00
FN041	Child	\$29.00
FN042	Adult	\$29.00

Silicon Masks



FN030	Small Infant	\$10.00
FN031	Large Infant	\$13.00
FN032	Small Child	\$16.00
FN033	Large Child	\$22.00
FN034	Small Adult	\$25.00
FN035	Large Adult	\$27.00

AUTOCLAVABLE TO 134°

Oxygen Resuscitation Kit

The Oxygen Resuscitation Kit is for rapid response administration of oxygen. A variety of sized internal pouches gives many options for accessories and provides for the supply of oxygen therapy.

Contents include:

- Multiflow regulator 0.25-25 LPM for C cylinder
- Adult resuscitator with reservoir and masks
- Sizes 00-4 Airways
- Cylinder key wheel
- Oxygen tubing
- 3 x adult therapy masks
- 3 x child therapy masks
- CPR Pocket mask



Package does NOT include cylinder. Please contact BOC gases. Ph: 131262 for cylinder

FN001	Contents + Carry Bag	\$625.00
-------	----------------------	-----------------

GST Free

Oxygen Masks

		Each	50+
BIA1624EE	Oxygen Mask Adult	\$1.50	\$1.25
BIA1625EE	Oxygen Mask Adult + 2m tube	\$2.75	\$2.30
BIA1626EE	Oxygen Mask Child	\$1.50	\$1.25
BIA1627EE	Oxygen Mask Child + 2m tube	\$2.75	\$2.30
BIA1623EE	Oxygen Tubing 2m	\$1.75	\$1.50



Nasal Oxygen Cannulas

		Each	50+
IN308MM	Nasal Oxygen Cannula Curved Flared Tip 3.0m	\$3.20	\$2.75



AccuREG Multi-Flow Piston O2 Regulator

The AccuReg Multi-Flow Piston O2 Regulator is constructed from a solid piece of brass, creating a uni-body style regulator that provides added strength and safety. The quality crafted internal brass components provide a safe and effective regulator. It comes with 2 Self Sealing high pressure outlets. This device is designed for operation with compressed medical oxygen. The fixed pressure regulator reduces the high pressure of the oxygen cylinder to a pressure compatible with the adjustable flow portion of the regulator. The flow meter is intended for use where masks, catheter and cannula administration of gas is required. The numbers on the flow meter dial indicate litres per minute (lpm) flow rate.

Flow settings: 2, 3, 4, 5, 6, 8, 9, 10, 15, 25 lpm



\$375

FN007
GST Free

CGA 870 PIN index fitting for all cylinders

Oxygen Therapy Regulator 0-25 LPM

The Piston Oxygen Therapy Regulator is a reliable oxygen therapy delivery device. It offers lightweight and compact design allowing it to be used in most therapy applications.

Features:

- All brass internal
- Barb outlet
- 11 therapy flow settings
- Weighs approx. 340g

\$159

FN008



Flowmeters

Features:

- Precision engineered
- Pressure compensated design accurately measures flow
- Easy to read
- Scratch resistant
- Virtually unbreakable transparent cover tube
- Colour coded handwheel



		Each
E86615AIR	Air, 0-15LPM, Ball Type, Model 866	\$145.00
E86615LOX	Oxygen, 0-15LPM, Ball Type, Model 866	

Gas Bottle Holder

ParagonCare

- 13L x 13W x 35H cm
- Hooks over trolley
- To suit C size cylinder

\$115

AX534



Oxy Bottle Trolley

ParagonCare

Powder coated metal frame with 150mm wheels

C Size

\$225

AX537



D Size

\$259

AXH0546



Regulators

Features:

The Model 801 Single Stage Preset Regulator possesses the sensitivity and flow capacity of much larger regulators in a compact lightweight design. Preset at 400kPa, the Model 801 is engineered to provide equipment downstream with constant pressure.

Features:

- Proven seat design
- Fully chromed
- Large, easy to read gauges Extensive monel filter
- Complies to AS 2896 & AS 4267



		Each
E801Z258	Nitrous Oxide, Yoke Inlet, Model 801	\$295.00
E801Z259	Air, Yoke Inlet, Model 801	\$385.00
E801Z260	Oxygen, Yoke Inlet, Model 801	\$310.00



HEINE® Classic+ Fiber Optic (F.O.) Laryngoscope Sets

Set complete with:
Paed 1, Mac 2, Mac 3 blades, spare XHL Xenon bulb, zipper case

		Each
HEF1191086	With Standard F.O. Battery Handle	\$1,375.00
HEF1923865	With Standard F.O. 4NT Rechargeable Handle (Li-ion) and NT4 Table Charger	\$2,055.00



As above, but with LED illumination, without spare bulb

		Each
HEF1918860	With Standard F.O. LED Battery Handle	\$1,515.00
HEF1994865	With Standard F.O. 4 LED NT Rechargeable Handle (Li-ion) and NT4 Table Charger	\$2,175.00
HEF1918820	Easy Clean Laryngoscope Set LED Battery Handle 3 blades	\$1,485.00



HEINE® XP Emergency Set with XP Disposable Blades

Specially-developed for emergency-use featuring the HEINE XP® Disposable Laryngoscope Blades. These sets are equipped with XP blades of the same anatomically-correct shape as the Classic+ blade, the world's best-selling F.O. blade with integrated fiber optic bundle. The HEINE XP Emergency Laryngoscope Set is complete with the 5 most frequently used blade types and sizes for infants and adults. Instantly-ready for use with F.O. SP Handle and replaceable batteries.

		Each
HEF5710815	with blades Miller 0, Miller 1, Mac 2, Mac 3, Mac 4, F.O. SP Handle 2.5 V, two alkaline batteries, in case	\$242.00



HEINE® XP Emergency Set Accessories

		Each
HEF0122815	F.O. SP Laryngoscope Handle 2.5 V, with bulb	\$137.00
HEX0188035	XHL Xenon Halogen spare bulb for F.O. SP Handle	\$36.00
HEF0022740	Refill blades	\$42.00



Opening an account with us is easy.

Call **1300 557 651** and speak to one of our Customer Service team members. Alternatively you can open an account online at www.vitalmed.com.au/new-accounts



Conventional Disposable Laryngoscope Set LED

Conventional disposable laryngoscope kit solves contamination problems, eliminates the cost and time spent cleaning laryngoscope blades and returning them for use.

- Complete set of most popular laryngoscope blades
- Blades made of sturdy surgical stainless steel
- Kits supplied with chrome plated handle
- Convenient, ready-to-use
- Most economical - sterile
- Vinyl zippered pouch

Kit includes:

- Medium handle, MacIntosh blades size: 2, 3, 4
- Miller blades size: 0, 1, 2, 3

\$140

EC55133



Single Use Fibre Optic Laryngoscope Set

Disposable plastic laryngoscope set is manufactured with high quality polycarbonate. It consists of a AA handle, three Macintosh blades, 3 Miller blades and a plastic storage case.

A strong single use item designed for field use.

\$115

AD1151



Greenline®/D™ F.O Disposable Laryngoscope Set Xenon

Disposable sterile, fibre optic, Stainless Steel blade, with chrome plated handle in zippered pouch

Kit Includes:

- Medium handle,
- MacIntosh blades size: 2, 3, 4
- Miller blades size: 0, 1, 2, 3

\$189

EC55333



Endotracheal Tubes

	Size	Each	20+
BO183050	5.0	\$7.95	\$7.50
BO183060	6.0		
BO183065	6.5		
BO183070	7.0		
BO183075	7.5		
BO183080	8.0		
BO183085	8.5		
BO183090	9.0		

	Size	Each	20+
BO111020	2.0	\$11.50	\$10.25
BO111025	2.5		
BO111030	3.0		
BO111035	3.5		
BO111040	4.0		
BT111045	4.5		
BT111050	5.0		
BT111060	6.0		

Cuffed



Uncuffed



Guedel Airways

	Size	Each	50+
MA7838	Sz 00 Orange 40mm	\$0.65	\$0.50
MA7636	Sz 0 Blue 50mm		
MA7635	Sz 1 Purple 60mm		
MA7634	Sz 2 White 70mm		
MA7633	Sz 3 Green 80mm		
MA7632	Sz 4 Yellow 90mm		
MA7631	Sz 5 Red 100mm		



All GST Free

Resuscitation Face Shield Key Ring

You are always prepared should an emergency situation arise

This nylon key ring case contains one transparent plastic mask with one-way valve.

\$5.95

FN006



Pocket Mask With Oxy Port & Tubing

Pocket mask with Oxy Port and Tubing for emergency CPR.

\$9.50

MDI001
GST Free





Econ-o-mist Forte

Nebuliser Therapy System

- Recommended for occasional use
- Impact resistant case
- Compact size

Includes:

- Tubing, jet nebuliser bowl and masks for child and adult
- Thermal fuse for protection



\$109

AE650-040
GST Free



Ventalair MAX

Allersearch Nebuliser

- Reliable nebuliser therapy system for regular use
- All accessories included
- 4 year warranty



\$199

AE650-044
GST Free



Compact Space Chamber Plus

- New Cross Valve Technology™
- May assist in more medication being delivered into the lungs
- As used in Australian Hospitals
- Designed and developed in Australia
- Single Patient Use
- 1 Year Warranty



\$13

AE490-502
GST Free



Nomad NexGen Travel Nebuliser System

Includes:

- Multi-volt power supply unit
 - 12 volt DC cable for car lighter adaptor
 - Accessories kit with: jet nebuliser bowl, adult mask, child mask, mouthpiece, air tube and elastics
 - Carry Bag and Spares Kit
 - Warranty 4 years (does not include consumables)
- Battery purchased as a separate item



\$225

AE650-043
GST Free

The NOMAD NexGen nebuliser system has the option of being operated from its own rechargeable battery.

Nomad NexGen Battery

For Nomad Travel Nebuliser



\$175

AE650-059N
GST Free



Nebuliser Masks Bowls & Tubing All GST Free

		Each	50+
BIA162900	Nebuliser Mask Adult	\$2.95	\$2.50
BIA163200	Nebuliser Mask Child	\$2.95	\$2.50
BIA1623EE	Oxygen tubing 2m	\$1.75	\$1.50
BIA1636EE	Nebuliser Bowl Aerflo	\$3.50	\$3.25



Breath A Tech Spacer

\$23.95

ZZ2050



Nebuliser Kits

Mask, Tubing and Nebuliser bowl

All GST Free		Each	50+
BES31612	Adult	\$2.95	\$2.50
BES31622	Child	\$2.95	\$2.50



Download the **FREE** Vital Medical Supplies **APP**

VISIT PAGE 10 FOR MORE INFORMATION

NEW **FREE** download available from:

