

# seca 286

## Wireless ultrasonic measuring station for height and weight with voice guidance

### EMR-integrated

The seca 286 ultrasonic measuring station offers fully automatic measuring and was designed to be user-friendly for hands-free operation.

Voice guidance can be optionally selected, combined with the included poster and sticker on the measuring station illustrates the correct stance and posture. This makes it easy for the patient to use without any assistance from the nursing staff.



**\$3,295**

PN286

**2 Year  
Warranty**



seca 466




### Intelligently precise

- Capacity: 300 kg
- Graduation (g): 50 g < 150 kg > 100 g
- Measuring range in cm: 60-210 cm
- Graduation [Measuring Rod]: 1 mm
- Power supply: Power adapter
- Measure (W x H x D): 434 x 2270 x 466 mm
- Approval class: III (3)
- Wireless connection to EMR or Direct print to any printer via Seca SE466 dongle
- Optional printer (shown)



**seca**  
Precision for health



## seca 700

Mechanical Column Scales With Eye-Level Beam

- Capacity: 220kg
- Graduation: 50g
- Reset-to-zero function
- Height rod included
- 2 year warranty

**Measuring rod:**

- Measuring range: 60-200cm
- Graduation: 1mm
- 2 year warranty

### \$745

PN700



## seca 703S

When weight, height and BMI have to be obtained, the new seca 703s can do it all for toddlers to obese adults. The seca 703s with integrated stadiometer is a complete measuring station that can measure both weight and height in just one time-saving step.

- Capacity: 300kg
- Graduation: 50g < 150kg > 100 g
- Dimensions (WxHxD): 360 x 930 x 520 mm
- Platform (WxHxD): 335 x 80 x 345 mm
- Measuring range: 7-230 cm
- Weight: 17.6 kg
- Power supply: Batteries / power adapter (optional SE400)
- Data transmission: seca 360° wireless technology to PC or seca printer
- Functions: TARE, pre-TARE, mother/child-function,
- HOLD, auto-HOLD, BMI, auto-CLEAR, damping, automatic switch-off, kg/lbs switch-over, SEND/PRINT
- 2 year warranty

### \$1,595

PN703S



## seca 284

Wireless Measuring Station For Body Height And Weight



**Scale:**

- Capacity: 300kg
- Graduation: 50g
- Functions: TARE, pre-TARE, mother/child-function, HOLD, auto-HOLD, auto-BMI, CLEAR, auto-CLEAR, auto weighing range switch-over, kg/lbs switch-over, damping, SEND/PRINT
- EMR Compatible
- Optional: printer seca 466

**Measuring rod:**

- Measuring range: 30-220cm
- Graduation: 1mm
- Functions headpiece: HOLD, user-defined zero setting, automatic switch-off, SEND/PRINT
- 2 year warranty

### \$2,495

PN284



## seca 769

Electronic Column Scale With Bmi Function

- Capacity: 200kg
- Graduation: 100g
- Functions: TARE, HOLD, auto-HOLD, BMI, kg/lbs switch-over, automatic switch-off
- Optional: measuring rod seca 220 or seca 224, switch-mode power adapter seca 400
- 2 year warranty

### \$755

PN769





## seca 813

Electronic Flat Scale with Very High Capacity

**\$210**

PN813

- Capacity: 200kg
- Graduation: 100g
- Functions: auto-HOLD, step-off, kg/lb switch over, tip-on automatic switch-on, automatic switch-off.
- 2 year warranty



## seca sensa 804

With chrome electrodes for precise measuring

**\$139**

PN804

- Capacity: 150kg
- Graduation: 100g
- 24 electrodes. For analysis of body fat and body water (BF/BW)
- Automatic user recognition and automatic BMI calculation. For fast measuring results
- Functions: auto-HOLD, BMI, BF/BW-measurement, user memory, automatic user recognition, tip-on automatic switch-on, automatic switch-off
- 2 year warranty



## seca aura 807

Advanced technology for enhanced comfort.

**\$129**

SEC807

- Capacity: 150kg
- Graduation: 100g
- Functions: auto-HOLD, step-off, tip-on automatic switch-on, automatic switch-off
- Stable construction for users up to 150kg
- Roughened glass surface for maximum slip resistance
- Large LCD figures for easy reading
- 2 year warranty



## seca clara 803

Clarity in shape and measurements

|        |       |                |
|--------|-------|----------------|
| PN803W | White | <b>\$99.00</b> |
| PN803B | Black | <b>\$99.00</b> |

- Capacity: 150kg
- Graduation: 100g
- Functions: auto-HOLD, tip-on automatic switch-on, automatic switch-off
- Two-component rubber coating. For pleasant footing
- Four stable rubber feet. For safe standing
- 28-millimeter LCD figures. For easy reading
- 2 year warranty



## seca 876

Flat Scale For Mobile Use

**\$595**

SEC876

- Capacity: 250kg
- Graduation: 100g < 150kg > 200g
- Functions: mother/child function, kg/lbs switch-over, tap-on automatic switch-on, automatic switch-off.
- 2 year warranty
- Optional: carrying case seca 414

**seca**

## seca 874Dr

User-friendly and specially designed for use in medical practices.

- Capacity: 200kg
- Graduation: 50 g < 150 kg > 100 g
- Dimensions: 321(w) x 61(H) x 361(D) mm
- 6 batteries, type AA
- 2 year warranty
- Dual display for patient and doctor
- Mother/child-function, Auto-HOLD



**\$499**

SEC874

**seca**

## Soehnle Digital personal scales

Simple use for mobile weighing

- Max. Capacity 150 kg
- Less weight, 1.4kg
- Easy-to-read LCD display, 26 mm digit height
- Hold function: The weight value will be held for approx. 6 seconds
- Automatic zeroing at switch-on
- Units in kg, lb or st, switchable at the bottom of the scale
- Energy-saving auto-off function after approx. 8 seconds



**SOEHNLE**  
PROFESSIONAL

**\$99**

SO7810

**NEW**

# SOEHNLE

## PROFESSIONAL

# Superior Quality and Reliability



Soehnle is established globally as one of the leading brands of household and personal scales with a solid reputation for the highest quality, reliability, and inventive design.

Soehnle scales are synonymous with precise measurement. With a heritage extending back to 1868, Soehnle has a long history of producing accurate measuring devices in the private and commercial sector, and their strong, innovative spirit continues to drive the development of new concepts and designs.

Soehnle scales combine renowned German technological excellence with sophisticated design.

Distributed by





## Soehnle Digital Floor Scale

Simple to use

- Weighing range: 250 kg, Division: 100 g
- Simple use
- Easy to read LC-Display, 30 mm digit height
- Tap-on start-function
- Units in kg or lb
- Anti-slipping, structured weighing platform
- Energy-saving auto-off-function
- Batteries (4 x 1,5 V AA) included in scope of delivery

**\$205**

SO7815

**SOEHNLE**  
PROFESSIONAL



## seca culina 852

Digital portion and diet scales

- Capacity: 6.6 lbs / 3 kg
- Graduation: 0.05 oz / 1 g
- Dimensions (WxHxD): 210 x 140 x 210 mm
- Dimensions/scale platform (WxHxD): 140 x 53 x 140 mm
- Weight: 710 g
- Functions: TARE, automatic switch-off
- Power supply: Batteries
- 2 year warranty
- Includes Bowl with food safe surface

**\$125**

ECO852

**seca**



## seca 644

Digital multifunctional handrail scale with wireless transmission

- Capacity: 360kg
- Graduation: 50 g
- Dimensions: 640(w) x 1,150(H) x 660(D) mm
- Rechargeable batteries
- 2 year warranty

**\$2,795**

PN644

**seca**



## seca 719

With Four Load Cells For Increased Stability

- Capacity: 180kg
- Graduation: 100g
- Display on a high column for easy reading
- High stability for demanding users.
- Functions: auto-HOLD, step-off, tip-on automatic switch-on, automatic switch-off
- 2 year warranty

**\$315**

PN719

**seca**



## seca 954

Digital chair scale with wireless transmission

- Capacity: 300 kg
- Graduation: 50 g < 150 kg > 100 g
- Dimensions (WxHxD): 563 x 897 x 978 mm
- Weight: 24.5 kg
- Power supply: Power adapter/rechargeable batteries
- Data transmission: seca 360° wireless technology
- Functions: TARE, pre-TARE, HOLD, auto-HOLD, BMI, auto-CLEAR, lbs/kg switch-over, damping, automatic switch-off, SEND/PRINT
- 2 year warranty

**\$2,495**

SE954

**seca**



## seca 952

Chair scale for weighing while seated


- Capacity: 200kg
- Graduation: 100g
- Functions: TARE, auto-HOLD, kg/lbs switch over, automatic switch-off
- Battery operated
- 2 year warranty

**\$1,650**

PN952

**seca**





## Soehnle Mechanical Personal Scale

- Large and clear circle pointer display
- Easy to read in kg-steps
- Manual zero resetting possible
- Attractive chrome trim ring
- Sturdy steel housing
- Weighing platform with easy-to-clean plastic covering



**\$115**

SO6163

**SOEHNLE**  
PROFESSIONAL

## seca 786

Mechanical Column Scales With Large Round Dial

- Capacity: 150kg
- Graduation: 500g
- Functions: reset-to-zero function
- 2 year warranty
- Optional: measuring rod seca 224
- 2 year warranty



**\$625**

PN786

**seca**

## seca 750

Robust Mechanical Floor Scales

Pretty and stable with their pure, light design to make a good impression in any clinic. Powder coated steel casing, circular display and an easy to clean tread area of imitation leather.

- Capacity: 150kg
- Graduation: 1kg
- Dimensions (mm): 317W x 117H x 477D
- Functions: reset to zero function
- 2 year warranty



|         |       |                 |
|---------|-------|-----------------|
| PN750WW | White | <b>\$159.00</b> |
| PN750WB | Black | <b>\$159.00</b> |

**seca**

## seca 762

Robust Mechanical Personal Scales with Fine Graduation

- Capacity: 150kg
- Graduation: 500g
- Reset-to-zero function
- 2 year warranty



**\$199**

PN761

**seca**



## Opening an account with us is easy.

Call 1300 557 651 and speak to one of our Customer Service team members. Alternatively you can open an account online at [www.vitalmed.com.au/new-accounts](http://www.vitalmed.com.au/new-accounts)



## seca 310 Hanging Baby Scale

Mechanical circular dial scale for use far from any hospital

- Capacity: 25 kg
- Graduation: 50 g
- Dimensions (WxHxD): 220 x 225 x 52 mm
- Reset-to-zero function
- 2 year warranty



**\$135**

SEC310

**seca**

## seca 431 Backpack for Baby Scale

Digital baby scales can be transported safely in this compact, tailor-made backpack case made of water-resistant material, complete with shoulder strap cushioning. Ideal for mobile healthcare and hands-free transport of the baby scales seca 385 and seca 354.



**\$139**

SE431

**seca**

## Carts for Baby Scales

Sturdy mobile support carts with customised indentations to suit all seca baby scales

### seca 402:

- Dimensions (WxHxD): 715 x 800 x 435 mm
- Weight: 11.5 kg
- Max. load: 100 kg
- 4 bumpers
- 2 brakes
- Optional: paper roll holder code SE408

### seca 403:

- Dimensions (WxHxD): 715 x 800 x 435 mm
- Weight: 27.5 kg
- Max. load: 100 kg
- 4 bumpers
- 2 brakes
- Integrated drawer and shelf
- 2 year warranty

|       |                                    |                 |
|-------|------------------------------------|-----------------|
| SE402 | Cart Only                          | <b>\$725.00</b> |
| SE403 | Cart with shelf and drawer         | <b>\$995.00</b> |
| SE408 | Paper towel holder for scales cart | <b>\$169.00</b> |

seca 402



seca 403

**seca**

## seca 428 Carry Case

The baby scales seca 334 can be safely and efficiently transported in this customized carrying case of water-repellent material equipped with a shoulder strap and practical padding.

**\$175**

PN428



**seca**

## Soehnle Easy-Fold Baby Scale

Foldable baby scale – easy to transport

- Patented folding mechanism with magnetic lock
- Low space requirement, easy to transport
- Easy to read LCD display
- Tare function: Covers/paddings can be tared
- Hold function: The weight value will be held for a few seconds after unloading
- Switchable to kg/lb units
- Energy-saving auto-off function
- Battery operation with 2 x 1.5 V AAA alkaline batteries (included in scope of delivery)




**\$199**

SO8320

**SOEHNLE**  
PROFESSIONAL





## seca 374

Wireless Baby Scale With Extra Large Weighing Tray



- Capacity: 20kg
- Graduation: 5g < 10kg > 10g
- Functions: TARE, BMIF, auto-HOLD, auto-CLEAR, kg/lbs switch-over, damping, automatic switch-off, SEND/PRINT
- Optional: measuring rod seca 233, head and foot positioners seca 418, switch-mode power adapter seca 400
- 2 year warranty

### \$825

PN374



seca

## seca 354

Electronic baby scales with fine graduation, also useable as flat scales for children.

- Capacity: 20kg
- Graduation: 5g < 10kg > 10g
- Functions: TARE, BMIF, auto-HOLD, kg/lbs switch-over, automatic switch-off
- 2 year warranty
- Optional: measuring board seca 417, carrying case seca 414 (for combination with measuring board) and seca 413

### \$225

PN354



seca

## seca 334

Mobile Electronic Baby Scale

### Scale

- Capacity: 20kg
- Graduation: 5g < 10kg > 10g
- Functions: TARE, BMIF, auto-HOLD, kg/lbs switch-over, automatic switch-off
- 2 year warranty
- Optional: measuring rod seca 232, head positioner seca 419, switch-mode power adapter seca 400, carry case seca 428 for scale only. (Or larger carry case seca 414 for scale and measure rod)

### \$625

PN334



seca

## seca 232

Measuring Rod

- Measuring range: 35-80cm
- Graduation 1mm
- 2 year warranty

### \$399

SEC232



## seca 413

Carrying case for seca 354

### \$89

PN413



seca





## Soehnle Height Rod - Electric

Handy height rod for mobile use

- Ultrasound measuring technology, temperature-compensated
- Tilt sensor for correction if measuring instrument is not horizontal, resulting in an error message
- Low space requirement due to fold function
- Plastic housing
- LCD display, 9 mm digit height
- Acoustic acknowledgment at end of measurement and reading displayed
- Units in cm or inch
- Battery operation with 2 x 1.5 V AA batteries (included in scope of delivery instruments)



**\$195**

SO5003

**NEW**

**SOEHNLE**  
PROFESSIONAL

## Soehnle Measuring Tape - Mechanical

Small and handy

- Measuring tape for wall-mounting
- No wear parts and maintenance-free
- Units in mm



**\$22**

SO5002

**SOEHNLE**  
PROFESSIONAL

## seca 210

Mobile Measuring Mat For Babies And Toddlers

- Measuring range: 10-99 cm
- Graduation: 5mm
- 2 year warranty

**\$115**

PN210



## seca 218

Disposable Measuring Tape

- Made of high-quality paper to prevent cross infection
- Simple removal of measuring tape for immediate use
- Always at your fingertips thanks to the exclusive wall dispenser that can be simply installed anywhere
- Exact results read easily
- Scale from 0-100 cm
- Wall dispenser included
- 500 disposable tapes included
- Refills (Pkt/1000) also available

**\$115**

PN218



## seca 211

Disposable Head Circumference Tapes

- Made of sturdy material to prevent cross infection
- Always at your fingertips thanks to the exclusive wall dispenser that can be simply installed anywhere
- Exact results read easily
- Scale from 12-59 cm
- Wall dispenser included
- 500 disposable tapes included
- Refills (Pkt/1000) also available

**\$75**

PN211



## seca 212

Measuring Tape For Head Circumference of Babies and Toddlers

- Measuring range: 5-59cm
- Graduation: 1mm
- 2 year warranty

**\$8.50**

SEC212



**seca**

## Retractable Measuring Tape - 1.5 metres

**\$3.95**

SAS423



## seca 206

Mechanical Measuring Tape

- Measuring range: 0-220 cm
- Graduation: 1mm
- 2 year warranty

**\$35**

PN206



**seca**

## seca 201

Ergonomic Circumference Measuring Tape

- Measuring range: 0-205cm
- Graduation: 1mm
- 2 year warranty

**\$20**

SEC201



**seca**