



Nu-Tek Frequency Ultrasound Unit

- 1 & 3MHz ultrasound frequency
- Pulsed and continuous modes
- Variable duty cycles (16Hz, 48Hz and 100Hz)
- Ergonomically designed sound heads in 1cm², 5cm² sizes
- Contact control with both visual and audio feedback
- 20 preset programs & 80 custom programs provided



\$1,395

ACNTULT

ACNTRR	Nu-Tek Frequency Trolley	\$895.00
--------	--------------------------	-----------------

WHITELEY
ALLCARE

MCUBE BIOCON-900 Bladder Scanner

The BioCon-900 uses ultrasound signals to determine the bladder's outline and total volume of urine. The revolutionary real-time Pre-Scan function locates the bladder prior to scanning and improves the accuracy of results.

Immediate printout of results via the internal printer or use the included CubeScan PC software to upload saved data to a computer for later review. 24 Month Warranty.

\$9,850

EC-500



		Each
EC-900C	Roll Cart for Biocon-900 Bladder Scanner	\$1,075.00

Ultrasound Products

SCAN Ultrasound Gel

3.8Ltr

\$55

PKLUS1128S



Aquaflex Ultrasound Gel Pad

Ea

\$22

PKLUS0402



Transeptic Cleaning Solution

250ml

\$24

PKLUS0925



Acoustically correct for the broad range of ultrasound frequencies used today, Vue Ultrasound Gel is also water soluble and hypoallergenic. So it's a convenient, versatile way to provide optimal transmission in ultrasonic diagnosis, ultrasonic therapy, ECG, EEG and EMG.

As specialists in gels for clinical applications, we've formulated Vue Ultrasound Gel to offer the ideal viscosity for a wide variety of procedures. Working with our network of clinical contacts, and responding to their expert feedback, we've developed this product to be as effective and easy to use as possible

		Box/30	6+	12+
OM1156	20ml Sachets Vue ultrasound gel™	\$66.95	\$64.00	\$61.00
		Each	6+	12+
OM1157	250ml Bottle Vue ultrasound gel™	\$4.64	\$3.85	\$2.95
		Each		
OM1158	5L Bottle Vue ultrasound gel™	\$23.80		

Always read the label. Use only as directed.