



Experience the difference™



See the GentleCath™ Hydrophilic catheter and no-touch handling strip

GentleCath™ Hydrophilic Catheters

Help minimise the risk of infection with no-touch catheterisation

Hydrophilic catheters are coated with a hydrophilic polymer that reacts with water to create a smooth, slippery coating on the surface of the catheter. Compared to uncoated catheters, hydrophilic catheters have been shown to reduce the incidence of urinary tract infections.¹

REFERENCES: 1. Cardenas DD, Hoffman JM. Hydrophilic catheters versus noncoated catheters for reducing the incidence of urinary tract infections: a randomized controlled trial. *Arch Phys Med Rehabil.* 2009;90(10):1668-71.

ConvaTec (Australia) Pty Limited. ABN 70 131 232 570. PO Box 63, Mulgrave, VIC 3170. Phone: (03) 9239 2700
 Facsimile: (03) 9239 2742. www.convatec.com.au
 ConvaTec (New Zealand) Limited. PO Box 62663, Greenlane 1546 New Zealand. Phone: 0800 441 763. www.convatec.co.nz
 ®/™ Indicates a trademark of Unomedical a/s. ConvaTec Inc. February 2018. ©2018 ConvaTec Inc. AP-014361-MM CCC210

GentleCath™ Hydrophilic Urinary Catheters

The GentleCath™ Hydrophilic coated catheter includes

- a sterile water sachet
- wall attachment
- unique no-touch handling sleeve that helps minimise the risk of infection with no-touch catheterisation.



CO508988	Female FG8 21cm	\$4.60
CO508989	Female FG10 21cm	
CO508990	Female FG12 21cm	
CO508991	Female FG14 21cm	
CO508992	Female FG16 21cm	
CO508993	Male FG8 40cm	
CO508994	Male FG10 40cm	
CO508995	Male FG12 40cm	
CO508996	Male FG14 40cm	
CO508997	Male FG16 40cm	
CO508998	Male FG18 40cm	



Unomedical® 100% Silicone Foley Catheters

The improved range of Unomedical® silicone foley catheters offer a wide range of solutions for safe and comfortable, acute and long-term urinary drainage requirements.



- Recommended long-term use for up to 90 days after insertion

IN90742101	Catheters Foley SIL Paed FG06 1.5mL	\$28.00
IN74210308	Catheters Foley SIL Paed FG08 3mL	\$19.00
IP74210310	Catheters Foley SIL Paed FG10 5mL	\$19.00
IN90742112	Catheters Foley SIL FG12 10mL	\$9.50
UM41151410	Catheters Foley SIL FG14 10mL	
IN90742116	Catheters Foley SIL FG16 10mL	
IN90742118	Catheters Foley SIL FG18 10mL	
IN90742120	Catheters Foley SIL FG20 10mL	
IN90742122	Catheters Foley SIL FG22 10mL	
IN90742105	Catheters Foley SIL FG24 10mL	
IN90742130	Catheters Foley SIL FG16 30mL	
IN90742318	Catheters Foley SIL FG18 30mL	
IP74213020	Catheters Foley SIL FG20 30mL	
IN90742322	Catheters Foley SIL FG22 30mL	
IP74213024	Catheters Foley SIL FG24 30mL	





Unomedical® Intermittent Urinary Catheters

A broad range of intermittent urinary catheters are available to meet individual patient needs, manufactured from smooth PVC, featuring rounded eyes and a smooth profile for easy insertion.



IN00106040	Male Nelaton Catheter Firm 6FG 40cm		
IN00108040	Male Nelaton Catheter Firm 8FG 40cm		
IN00110040	Male Nelaton Catheter Firm 10FG 40cm		
IN00112040	Male Nelaton Catheter Firm 12FG 40cm	\$1.70	
IN00114040	Male Nelaton Catheter Firm 14FG 40cm		
IN00116040	Male Nelaton Catheter Firm 16FG 40cm		
IN00118040	Male Nelaton Catheter Firm 18FG 40cm		
IN00120040	Male Nelaton Catheter Firm 20FG 40cm		
IN00206040	Male Nelaton Catheter Firm 6FG 40cm DW		
IN00208040	Male Nelaton Catheter Firm 8FG 40cm DW		
IN00210040	Male Nelaton Catheter Firm 10FG 40cm DW		
IN00212040	Male Nelaton Catheter Firm 12FG 40cm DW	\$1.75	
IN00214040	Male Nelaton Catheter Firm 14FG 40cm DW		
IN00216040	Male Nelaton Catheter Firm 16FG 40cm DW		
IN00218040	Male Nelaton Catheter Firm 18FG 40cm DW		
IN00220040	Male Nelaton Catheter Firm 20FG 40cm DW		
IN00310040	Male Nelaton Catheter Soft 10FG 40cm		
IN00312040	Male Nelaton Catheter Soft 12FG 40cm		
IN00314040	Male Nelaton Catheter Soft 14FG 40cm	\$1.70	
IN00316040	Male Nelaton Catheter Soft 16FG 40cm		
IN00318040	Male Nelaton Catheter Soft 18FG 40cm		
IN00320040	Male Nelaton Catheter Soft 20FG 40cm		
IN00610040	Male Nelaton CatheterS Soft 10FG 40cm DW		
IN00612040	Male Nelaton CatheterS Soft 12FG 40cm DW		
IN00614040	Male Nelaton CatheterS Soft 14FG 40cm DW	\$1.75	
IN00616040	Male Nelaton CatheterS Soft 16FG 40cm DW		
IN00618040	Male Nelaton CatheterS Soft 18FG 40cm DW		
IN00620040	Male Nelaton CatheterS Soft 20FG 40cm DW		
IN01106023	Female Nelaton Catheter 6FG 22.5cm		
IN01108023	Female Nelaton Catheter 8FG 22.5cm		
IN01110023	Female Nelaton Catheter 10FG 22.5cm	\$1.85	
IN01112023	Female Nelaton Catheter 12FG 22.5cm		
IN01114023	Female Nelaton Catheter 14FG 22.5cm		
IN01116023	Female Nelaton Catheter 16FG 22.5cm		



Closed System Bag

Sterile 2 litre drainage bag with 120cm inlet tube

	Each	100+
SB340821	\$5.95	\$4.35

GST Free



NEW

- Soft silicone 'tack' layer provides immediate gentle adhesion to skin.
- Ideal for atraumatic removal without use of solutions to dissolve adhesive
- No sticky residue
- Effectively conforms to all skin types through multiple points of contact
- Re-adhesive tack enables easy repositioning and close skin inspection

Pack of 5

\$49.50

OM3004



NEW

- Atraumatic removal from the skin without use of solutions to dissolve adhesive
- No sticky residue left behind
- Can be repeatedly applied to same area of skin
- Re-adhesive tack enables easy repositioning
- No complex skin preparation required
- Reduced risk of allergic reactions due to silicone technology
- Showerproof, breathable and translucent



		Each
OM3031	Small 110 x 25mm	\$54.00
OM3030	Large 140 x 40mm	\$69.00



GREAT VALUE ON QUALITY MEDICAL CONSUMABLES

Urine Drainage Bags

Non Drainable 2 Litre Drainage Bag

	Each	200+	
SB331108	\$2.10	\$1.95	GST Free

2 Litre Drainage Bag with Outlet

	Each	250+	
SB320902	\$0.85	\$0.79	GST Free

