



Multi-purpose



3M™ Transpore™ Surgical Tape

# Transparent and tears easily.

- Clear, water-resistant plastic tape
- Good adhesion to skin and tubing
- Easy, straight, bi-directional tear

# 3M™ Transpore™ Surgical Tape

Application: **General Securement**

- Blood draws
- Dressings
- I.V. lines and tubing
- Non-critical tubes

**Ideal for intact skin** (for at-risk skin, use a barrier like 3M™ Cavilon™ No Sting Barrier Film)

**Tape material:** Plastic

Available in single-patient-use rolls ✓

 **Other multi-purpose 3M medical tapes:**

- 3M™ Transpore™ White Surgical Tape
- 3M™ Micropore™ Surgical Tape



## Ordering information

Catalogue Number	Size	Rolls/Box	Boxes/Case
1527-0	1.25 cm x 9.1 m	24	10
1527-1	2.5 cm x 9.1 m	12	10
1527-2	5 cm x 9.1 m	6	10
1527-3	7.6 cm x 9.1 m AUS only	4	10
Single-Patient-Use Rolls			
1527S-1	2.5 cm x 1.37 m AUS only	100	5



For more information, contact your 3M Sales Representative, or visit our website for more information.



### 3M Critical & Chronic Care Solutions Division

3M Australia Pty Limited  
 ABN 90 000 100 096  
 Building A, 1 Rivett Road  
 North Ryde, NSW 2113  
 1300 363 878  
[www.3M.com.au/healthcare](http://www.3M.com.au/healthcare)

3M New Zealand Limited  
 94 Apollo Drive  
 Rosedale Auckland 0632  
 New Zealand  
 0800 808 182  
[www.3M.co.nz/healthcare](http://www.3M.co.nz/healthcare)

3M, Transpore, Cavilon and Micropore are registered trademarks of 3M. Please recycle. © 3M 2017. All rights reserved.

AV011462395