

GAMMEX[®] Latex Textured

S U R G I C A L G L O V E S

Latex, powder-free surgical glove with a micro-textured finish for secure instrument handling and grip

This micro-textured surgical glove with special treatment on the external surface facilitates better grip on instruments under both dry and damp conditions while DERMASHIELD™ Technology allows for ease of donning especially during intra-operative changes.

KEY FEATURES & BENEFITS

- Micro-textured finish with special external surface treatment to ensure exceptional grip under dry and damp conditions
- Ideal thickness to optimize tactile sensitivity and protection
- DERMASHIELD™ Technology internal polymer coating facilitates dry and damp donning as well as intra-operative changes

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

RECOMMENDED FOR

- Sterile procedures in the O.R
- Enhanced grip in wet environments
- Dry and damp hand donning
- Extra comfort during long, sustained procedures
- Use as an outer glove when double gloving

FEATURED TECHNOLOGIES



GAMMEX[®] Innovation that Protects

Driven by innovation, GAMMEX[®] surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

GAMMEX[®] Latex Textured

S U R G I C A L G L O V E S

PRODUCT INFORMATION

Material	Natural rubber latex
Powder Content	Powder-free
Color	Natural
Shape	Anatomic with curved fingers
Internal Glove Surface	Polymer coated with DERMASHIELD™ Technology
External Glove Surface	Textured fingers and palm
Cuff Style	Straight with SUREFIT™ Technology
Thickness	Standard
Grip Level	High
Allergy Prevention	None
Double Gloving Recommendation	Outer glove
Freedom From Holes (AQL)	Manufacturing final release: 0.65 (Inspection level GI)
Tested For Use With Chemotherapy Drugs	No, not tested in accordance with ASTM D6978

ORDERING INFORMATION

Size	5 ½	6	6 ½	7	7 ½	8	8 ½	9
Ansell Product Code	331300655	331300660	331300665	331300670	331300675	331300680	331300685	331300690

PACKAGING AND STORAGE

Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

Product Standards	ASTM D3577, ASTM D7160, ASTM F739, EN 455 Part 1-4, EN ISO 374-1, EN 374-2 and -4, EN 16523-1, EN ISO 374-5, EN 420, ISO 10282
Quality / Environmental Standards	ISO 9001, ISO 13485, ISO 14001, ISO 11137-Part 1, EN 556

Contact your Ansell representative for ordering or more information.

North America
 US Tel: 800 952 9916
 CA Tel: 1-844-494-7854
 Email: insidesalesus@ansell.com

Central & South America
 Tel: +52(442) 296 20 50
 Email: cslac@ansell.com
Brazil
 Tel: +55-11-3356-3100
 Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
 Tel: +32 (0) 2 528 74 00
 Email: info@ansell.eu

Asia Pacific
 Tel: +603 8310 6688
 Email: apac.medical@ansell.com

Australia & New Zealand
 Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com