



BD Ultra-Fine™ Pen Needles

BD Ultra-Fine™ PRO 4mm Pen Needle

- BD Ultra-Fine™ 4mm Pen Needle is being enhanced with new ergonomic grips to make handling easier¹
- BD Ultra-Fine™ PRO 4mm Pen Needles will be identifiable through a refreshed packaging change.

Wider outer cover features a flat end and moulded ridges to support easier attachment and removal of the pen needle¹




Additional grips for easier handling.



4mm needle length helps ensure the medication is delivered into the fat layer and not the muscle

BD Ultra-Fine™ 5mm Pen Needle

- BD Ultra-Fine™ 5mm Pen Needle is being enhanced with EasyFlow™ Technology and 5-Bevel Comfort.
 - EasyFlow™ Technology is an extra thin wall cannula that improves the flow rate of insulin for a faster injection¹
 - 5-bevel needle designed to reduce penetration force and increase comfort.¹
- BD Ultra-Fine™ 5mm Pen Needles with EasyFlow™ Technology will be identifiable through a refreshed packaging change.



Vital Code:	Description	RRP	Price	SAVE
BD320566	ULTRA-FINE PEN NDLS 32G x 4mm 100s	\$29.70	\$22.00	\$7.70
BD320569	ULTRA-FINE PEN NDLS 31G x 5mm 100s	\$29.70	\$22.00	\$7.70

1. Compared to standard hub BD Pen Needles.
Offer valid from 1st - 31st July 2019



Vital Medical Supplies brings you great brands at a great price!

To place an order please contact Vital Medical Supplies on vitalconsultant@vitalmed.com.au